



**OsteoBiol**<sup>®</sup>  
by Tecnos

YEARS  
**20**  
ANNIVERSARY


## Bone Grafting Materials

REGENERATION SCIENCE


INSPIRED BY NATURE




SEM image of an OsteoBiol® Gen-Os® granule colonised by osteoblasts  
Author: Prof Ulf Nånmark, University of Göteborg, Sweden




4-months histology showing remnants of porcine dermal matrix surrounded by connective tissue  
Experimental study on dogs. Htx-eosine  
Author: Prof Ulf Nånmark, University of Göteborg, Sweden



LM image of an OsteoBiol® Lamina hydrated with blood: vascularisation enhanced by the presence of the original vascular canals  
Author: Prof Ulf Nånmark, University of Göteborg, Sweden



Part of a biopsy showing newly formed bone after treatment with OsteoBiol® Putty. Htx-eosine. Magnification x20  
Author: Prof Ulf Nånmark, University of Göteborg, Sweden






## **OUR MISSION**

*«To produce a xenogenic bone substitute as similar as possible to autogenous bone»*

**Giuseppe Oliva MD**  
R&D Director  
**Tecnoss S.r.l.**







# **THE OSTEObIOL® DUAL-PHASE HETEROLOGOUS BONE MATRIX**


*OsteoBiol® is the family of biomaterials produced by TecnoSS® for the dental and maxillo-facial surgeons.*

*In each OsteoBiol® granule, besides its mineral phase, the TecnoSS® process retains the xenogenic collagen phase with its precious biological properties, making it biocompatible and ideal for grafting and augmentation purposes.*

*Avoiding high process temperatures, the OsteoBiol® bone matrix prevents ceramization, maintaining a chemical composition extremely similar to autogenous bone<sup>(1)</sup>, and therefore gradually resorbable and replaceable by newly formed bone.*

(1) Figueiredo M, Henriques J, Martins G, Guerra F, Judas F, Figueiredo H  
**Physicochemical characterization of biomaterials commonly used in dentistry as bone substitutes - comparison with human bone**  
Journal of Biomedical Materials Research Part B: Applied Biomaterials, 2010 Feb; 92(2):409-19





SEM image of an OsteoBiol® Gen-Os® granule: osteoblastic colonisation. Magnification x3000  
Author: Prof Ulf Nånberg, University of Göteborg, Sweden




# **HIGH BIOCOMPATIBILITY**

*The chemical structure of each OsteoBiol® dual-phase granule, its ideal porosity and collagen content, make it a valid scaffold and substrate for osteoblasts anchorage, proliferation and new bone apposition<sup>(2)</sup>.*

(2) Nannmark U, Sennerby L

**The bone tissue responses to prehydrated and collagenated cortico-cancellous porcine bone grafts: a study in rabbit maxillary defects**  
Clinical Implant Dentistry and Related Research, 2008 Dec;10(4):264-70





**Biopsy from upper jaw region retrieved at four months showing gradual osteoclastic resorption of an OsteoBio!® GTO® granule**  
Author: Prof Ulf Nannmark, University of Göteborg, Sweden

# GRADUAL RESORPTION

*Autogenous bone is gradually replaced by newly formed bone: similarly, the OsteoBiol® bone matrix allows progressive osteoclastic resorption, with simultaneous new bone apposition.*

*Cells receive nutrients from newly formed vessels, that are able to colonize adequately the grafted site.*

*New bone grows in and around the OsteoBiol® granules<sup>(3)</sup>, which are partially but significantly replaced by vital bone at re-entry time.*

(3) Giuliani A, Iezzi G, Mazzoni S, Piattelli A, Perrotti V, Barone A

**Regenerative properties of collagenated porcine bone grafts in human maxilla: demonstrative study of the kinetics by synchrotron radiation microtomography and light microscopy**  
Clinical Oral Investigations, 2018 Jan;22(1):505-513






# **VASCULARIZATION IS THE KEY FOR CLINICAL SUCCESS**

*Dual-phase biomaterials are progressively resorbed by osteoclasts and replaced by new vital bone produced by osteoblasts, similarly to autogenous bone grafts. Both types of cells live thanks to blood supply, which is critical and essential for the success of any bone regeneration procedure.*

*The progressive resorption of OsteoBiol<sup>®</sup> granules allows an adequate colonization of the grafting site by new vessels, and is therefore a positive and significant factor within the regenerative process<sup>(4)</sup>.*

(4) Rombouts C, Jeanneau C, Camilleri J, Laurent P, About I  
Characterization and angiogenic potential of xenogeneic bone grafting materials: Role of periodontal ligament cells  
Dental Materials Journal, 2016 Dec 1;35(6):900-907







# THE ROLE OF COLLAGEN

*Collagen favours MSC differentiation and enhances osteoblasts proliferation<sup>(5,6)</sup>: it is considered as the ideal substrate for bone forming cells. OsteoBiol® dual-phase particulate bone substitutes contain approximately 22% collagen.*

*Furthermore, collagen gel mixed with dual-phase collagenated granules packed in syringes improves the handling and the stability of the graft, reducing also operatory time and risk of contamination.*

(5) Brunelli G, Sollazzo V, Carinci F, Palmieri A, Girardi A, Monguzzi R  
**OsteoBiol® influences osteogenic differentiation of adipose derived stem cells**  
European Journal of Inflammation, 2011, Vol. 9, no. 3 (S), 103-107

(6) Jeanneau C, Le Fournis C, About I  
**Xenogeneic bone filling materials modulate mesenchymal stem cell recruitment: role of the complement c5a**  
Clinical oral investigations; 2019 oct 23






Fig. 1




Fig. 2



Fig. 3




Fig. 4



Fig. 5




Fig. 6



Fig. 7



Fig. 8



Fig. 9

# **A SPECIFIC PRODUCT FOR EVERY CLINICAL INDICATION**

*OsteoBiol® is not only a marvellous collagenated bone matrix: it is a complete family of biomaterials specifically designed for bone and soft tissue augmentation in dentistry. For every clinical indication a dedicated product has been developed, with the goal of providing the best handling, the ideal granulometry and consistency, and finally optimal regenerative results in adequate re-entry time.*

*Enjoy one of the widest and most complete product ranges, with the security and support of 15 years of clinical research: you will experience that today it is finally possible to achieve predictable clinical success<sup>(7)</sup> without the availability limitations of autogenous bone.*

(7) Checchi V, Felice P, Zucchelli G, Barausse C, Piattelli M, Pistilli R, Grandi G, Esposito M

**Wide diameter immediate post-extractive implants vs delayed placement of normal-diameter implants in preserved sockets in the molar region: 1-year post-loading outcome of a randomised controlled trial**  
European Journal of Oral Implantology, 2017;10(3):263-278


Images authors pag 14: Fig.1, Dr Patrick Palacci, Marseille, France | Fig. 2-9, Dr Roberto Rossi, Genova, Italy



**OsteoBiol®** and natural human bone have the same density and very similar physico-chemical properties

Figueiredo et al. J Biomed Mater Res B: Appl Biomater, 2010 Feb; 92(2):409-19


### Mineral content



Figueiredo et al. J Biomed Mater Res B: Appl Biomater, 2010 Feb; 92(2):409-19


**Gen-Os®** has a higher angiogenic potential compared to anorganic xenografts

Rombouts et al. Dent mater J, 2016 Dec 1;35(6):900-907



In ridge preservation collagenated biomaterials show significant smaller volume reduction and basal area shrinkage compared to slowly resorbable xenografts

Barone et al. Clin Oral Implants Res, 2016 Nov;27(11):E105-E115



**OsteoBiol®** bone matrix promotes osteoblast differentiation and bone regeneration

Brunelli et al. Eur J Inflamm, 2011, Vol. 9, no. 3 (S), 103-107



**OsteoBiol®** bone scaffolds absorb growth factors secreted by MSCs and improve bone tissue repair

Mijiritsky et al. Materials, 2017 Sep 8;10(9)

## **KEY SCIENTIFIC DATA**

*Over 190 articles have been published on peer-reviewed journals during the last 15 years, proving with in-vitro, experimental and clinical studies the outstanding biological properties and clinical performance of the OsteoBiol® collagenated biomaterials.*





# ***PATIENTS FIRST***

*Combining the best skills and the best materials, within the limits and guidelines provided by scientific evidence, is the key for clinical success: however let us all remember that the patients are and will always be the center of all our attentions.*

*Meeting their expectations, helping them to recover function and aesthetics with long term success<sup>(8,9)</sup> is the greatest reward for any surgeon and fulfillment of our company mission.*

(8) Barone A, Orlando B, Tonelli P, Covani U

**Survival rate for implants placed in the posterior maxilla with and without sinus augmentation: a comparative cohort study**  
Journal of Periodontology, 2011 Feb; 82(2):219-26

(9) Scarano A, Piattelli A, Assenza B, Quaranta A, Perrotti V, Piattelli M, Iezzi G

**Porcine bone used in sinus augmentation procedures: a 5-year retrospective clinical evaluation**  
Journal of Oral and Maxillofacial Surgery, 2010 Aug;68(8):1869-73



# OsteoBiol® products clinical evidence\*

## Gen-Os®

Collagenated heterologous cortico-cancellous bone mix  
Granulometry 250-1000 µm  
For information on OsteoBiol® Gen-Os® see page 24

## mp3®

Pre-hydrated collagenated heterologous cortico-cancellous bone mix  
Granulometry 600-1000 µm  
For information on OsteoBiol® mp3® see page 32

## GTO®

Pre-hydrated collagenated heterologous cortico-cancellous bone mix  
Granulometry 600-1000 µm  
For information on OsteoBiol® GTO® see page 36

## Putty

Pre-hydrated collagenated heterologous cortico-cancellous bone paste  
Granulometry up to 300 µm  
For information on OsteoBiol® Putty see page 40


## Gel 40

Pre-hydrated collagenated heterologous cortico-cancellous bone gel  
Granulometry up to 300 µm  
For information on OsteoBiol® Gel 40 see page 44

### ALVEOLAR REGENERATION



### MAXILLARY SINUS LIFT




CRESTAL ACCESS ONLY



### PERI-IMPLANT DEFECTS



### HORIZONTAL AUGMENTATION



### VERTICAL AUGMENTATION

INLAY TECHNIQUE




### PERIODONTAL REGENERATION



### SOFT TISSUE AUGMENTATION

## Apatos

Cortico-cancellous and cortical bone  
Granulometry 600-1000  $\mu\text{m}$   
For information on OsteoBio<sup>®</sup> Apatos  
see page 48



WITH TSV GEL



IN ASSOCIATION



WITH LAMINA AND TSV GEL

## Sp-Block

Collagenated heterologous  
cancellous block  
For information on OsteoBio<sup>®</sup> Sp-Block  
see page 54

## Evolution

Heterologous collagen membrane  
For information on OsteoBio<sup>®</sup> Evolution  
see page 62




## Lamina

Collagenated heterologous cortical bone  
For information on OsteoBio<sup>®</sup> Lamina  
see page 70





## Derma

Collagen dermal matrix  
For information on OsteoBio<sup>®</sup> Derma  
see page 66



BONE LAYER TECHNIQUE





# BONE SUBSTITUTES



Fig. 1




Fig. 2



Fig. 3




Fig. 4




Fig. 5



Fig. 6




Fig. 7











Fig. 8

# OsteoBiol® Dual-Phase bone substitutes

## HETEROLOGOUS BONE MATRIX




Cortico-cancellous collagenated matrix					<b>Apatos Cortical</b> cortical bone
	TSV Gel	Collagen gel			
	Pre-hydrated	Pre-hydrated	Pre-hydrated	Pre-hydrated	
<b>Gen-Os®</b>	<b>GTO®</b>	<b>mp3®</b>	<b>Putty</b>	<b>Gel 40</b>	<b>Apatos Mix</b>
100% collagenated bone mix	~80% collagenated bone mix ~20% TSV Gel	~90% collagenated bone mix ~10% collagen gel	~80% collagenated bone mix ~20% collagen gel	~60% collagenated bone mix ~40% collagen gel	cortico-cancellous bone mix
					
Heterologous cortico-cancellous collagenated bone mix	Heterologous cortico-cancellous collagenated pre-hydrated bone mix	Heterologous cortico-cancellous collagenated pre-hydrated bone mix	Heterologous cortico-cancellous collagenated pre-hydrated bone paste	Heterologous cortico-cancellous collagenated pre-hydrated bone gel	Heterologous microcrystalline hydroxyapatite
For more information on OsteoBiol® Gen-Os® see page 24	For more information on OsteoBiol® GTO® see page 36	For more information on OsteoBiol® mp3® see page 32	For more information on OsteoBiol® Putty see page 40	For more information on OsteoBiol® Gel 40 see page 44	For more information on OsteoBiol® Apatos see page 48





# Gen-Os<sup>®</sup>



***The advantages of a dual-phase biomaterial***  
***Collagenated heterologous cortico-cancellous bone mix***



# Characteristics and handling



## Tissue of origin

Cortico-cancellous heterologous bone mix

## Tissue collagen

Preserved

## Physical form

Slightly radiopaque granules

## Composition

100% granulated mix

## Granulometry

250-1000  $\mu\text{m}$

1000-2000  $\mu\text{m}$

## Re-entry time

4/5 months, depending on grafting site characteristics

## Packaging

Vial: 0.25 g, 0.5 g, 1.0 g, 2.0 g

## Product codes

250-1000  $\mu\text{m}$

M1052FS | 1 Vial | 0.25 g | Porcine

M1052FE | 1 Vial | 0.25 g | Equine

M1005FS | 1 Vial | 0.5 g | Porcine

M1005FE | 1 Vial | 0.5 g | Equine

M1010FS | 1 Vial | 1.0 g | Porcine

M1010FE | 1 Vial | 1.0 g | Equine

M1020FS | 1 Vial | 2.0 g | Porcine

M1020FE | 1 Vial | 2.0 g | Equine

1000-2000  $\mu\text{m}$

M0210FS | 1 Vial | 1.0 g | Porcine

M0220FS | 1 Vial | 2.0 g | Porcine

## GMDN code

46425

## CND code

P900402

## CHARACTERISTICS

A natural replicate of autologous bone, Gen-Os<sup>®</sup> conserves the same intimate structures<sup>(1)</sup> (matrix and porous form) and presents highly osteoconductive properties<sup>(2,3)</sup>. It is biocompatible and bioavailable, as recognized by tests made according to the ISO 10993 method conducted at Eurofins Biolab. Gen-Os<sup>®</sup> is gradually resorbable and provides support in bone neoformation helping to preserve the original graft shape and volume<sup>(4)</sup>.

Moreover, thanks to its collagen content, the product facilitates blood clotting and the subsequent invasion of repairing and regenerative cells, favouring *restitutio ad integrum* of missing bone. Because of its marked hydrophilia<sup>(5)</sup>, it can function as a carrier for selected medications and drugs<sup>(6)</sup> and it is ideal to mix with GFs<sup>(7)</sup>.

## HANDLING

Gen-Os<sup>®</sup> must always be hydrated and thoroughly mixed with either a few drops of sterile physiological solution (or patient's blood) to activate its collagen matrix and to enhance its adhesivity or with TSV Gel to increase graft stability in not self-contained defects. If necessary, it can as well be mixed with the drug selected for surgery.



SEM image of OsteoBioL<sup>®</sup> Gen-Os<sup>®</sup> granules. Magnif. x50

Author: Prof Ulf Nannmark, University of Göteborg, Sweden



Source: TecnoSS<sup>®</sup> Dental Media Library

Gen-Os<sup>®</sup>, a cortico-cancellous bone mix, has been the first product developed with the TecnoSS<sup>®</sup> innovative biotechnology and, due to its universal use, still is today the most demanded from the market. Gen-Os<sup>®</sup> has been successfully used and documented for alveolar ridge preservation<sup>(8)</sup> in combination with *Evolution* membranes: the application of this biomaterial limits the alveolar ridge width reduction that would naturally occur with spontaneous healing, preserving thus the alveolar ridge volume and allowing a correct second stage implant placement<sup>(9)</sup>. Gen-Os<sup>®</sup> has been used for lateral access maxillary sinus lift<sup>(3)</sup> and dehiscence regeneration<sup>(10)</sup>, always in association with *Evolution* membranes. Gen-Os<sup>®</sup> has been documented in periodontal regeneration of deep infrabony defects<sup>(11)</sup>. Due to its collagen content, once hydrated Gen-Os<sup>®</sup> becomes very sticky and hydrophilic<sup>(5)</sup>: it combines therefore extremely well with blood and is very stable once applied into the grafting site.

Its cortico-cancellous composition allows a progressive resorption of osteoclastic type, with in parallel a similar rate of new bone formation<sup>(2)</sup>: these unique properties allow a very good graft volume preservation, a healthy and well vascularized new bony tissue and, ultimately, a successful implant rehabilitation. Gen-Os<sup>®</sup> is in fact able to boost vascularization: in vitro<sup>(12)</sup> assays proved an increase in the secretion of VEGF by periodontal ligament cells (PDL) in the presence of Gen-Os<sup>®</sup>, as well as an enhanced proliferation of endothelial cells.



**LATERAL ACCESS SINUS LIFT**  
maxillary sinus floor augmentation




**PERIODONTAL REGENERATION**  
intrabony defects



**HORIZONTAL AUGMENTATION**  
two-wall defects



**DEHISCENCES AND FENESTRATIONS**  
peri-implant lesions



**CRESTAL ACCESS SINUS LIFT**  
osteotome technique



**ALVEOLAR REGENERATION**  
socket preservation

free animated videos  
on OsteoBiol<sup>®</sup> APP




## BIBLIOGRAPHY


- (1) FIGUEIREDO M, HENRIQUES J, MARTINS G, GUERRA F, JUDAS F, FIGUEIREDO H  
**PHYSICO-CHEMICAL CHARACTERIZATION OF BIOMATERIALS COMMONLY USED IN DENTISTRY AS BONE SUBSTITUTES - COMPARISON WITH HUMAN BONE**  
J BIOMED MATER RES B APPL BIOMATER, 2010 FEB; 92(2):409-19
- (2) NANNMARK U, SENNERBY L  
**THE BONE TISSUE RESPONSES TO PREHYDRATED AND COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS: A STUDY IN RABBIT MAXILLARY DEFECTS**  
CLIN IMPLANT DENT RELAT RES, 2008 DEC;10(4):264-70
- (3) CASSETTA M, PERROTTI V, CALASSO S, PIATTELLI A, SINJARI B, IEZZI G  
**BONE FORMATION IN SINUS AUGMENTATION PROCEDURES USING AUTOLOGOUS BONE, PORCINE BONE, AND A 50 : 50 MIXTURE: A HUMAN CLINICAL AND HISTOLOGICAL EVALUATION AT 2 MONTHS**  
CLIN ORAL IMPLANTS RES, 2015 OCT;26(10):1180-4
- (4) CARDAROPOLI D, CARDAROPOLI G  
**PRESERVATION OF THE POSTEXTRACTION ALVEOLAR RIDGE: A CLINICAL AND HISTOLOGIC STUDY**  
INT J PERIODONTICS RESTORATIVE DENT, 2008 OCT; 28(5):469-77
- (5) FIGUEIREDO A, COIMBRA P, CABRITA A, GUERRA F, FIGUEIREDO M  
**COMPARISON OF A XENOGENIC AND AN ALLOPLASTIC MATERIAL USED IN DENTAL IMPLANTS IN TERMS OF PHYSICO-CHEMICAL CHARACTERISTICS AND IN VIVO INFLAMMATORY RESPONSE**  
MATER SCI ENG C MATER BIOL APPL, 2013 AUG 1;33(6):3506-13
- (6) FISCHER KR, STAVROPOULOS A, CALVO GUIRADO JL, SCHNEIDER D, FICKL S  
**INFLUENCE OF LOCAL ADMINISTRATION OF PAMIDRONATE ON EXTRACTION SOCKET HEALING - A HISTOMORPHOMETRIC PROOF-OF-PRINCIPLE PRE-CLINICAL IN VIVO EVALUATION**  
CLIN ORAL IMPLANTS RES, 2015 OCT;26(10):1135-42
- (7) MIJIRITSKY E, FERRONI L, GARDIN C, BRESSAN E, ZANETTE G, PIATTELLI A, ZAVAN B  
**PORCINE BONE SCAFFOLDS ADSORB GROWTH FACTORS SECRETED BY MSCS AND IMPROVE BONE TISSUE REPAIR**  
MATERIALS, 2017 SEP 8;10(9)
- (8) CHECCHI V, FELICE P, ZUCCHELLI G, BARAUSSE C, PIATTELLI M, PISTILLI R, GRANDI G, ESPOSITO M  
**WIDE DIAMETER IMMEDIATE POST-EXTRACTIVE IMPLANTS VS DELAYED PLACEMENT OF NORMAL-DIAMETER IMPLANTS IN PRESERVED SOCKETS IN THE MOLAR REGION: 1-YEAR POST-LOADING OUTCOME OF A RANDOMISED CONTROLLED TRIAL**  
EUR J ORAL IMPLANTOL, 2017;10(3):263-278
- (9) FESTA VM, ADDABBO F, LAINO L, FEMIANO F, RULLO R  
**PORCINE-DERIVED XENOGRAFT COMBINED WITH A SOFT CORTICAL MEMBRANE VERSUS EXTRACTION ALONE FOR IMPLANT SITE DEVELOPMENT: A CLINICAL STUDY IN HUMANS**  
CLIN IMPLANT DENT RELAT RES, 2013 OCT;15(5):707-13
- (10) CASSETTA M, RICCI L, IEZZI G, DELL'AQUILA D, PIATTELLI A, PERROTTI V  
**RESONANCE FREQUENCY ANALYSIS OF IMPLANTS INSERTED WITH A SIMULTANEOUS GRAFTING PROCEDURE: A 5-YEAR FOLLOW-UP STUDY IN MAN**  
INT J PERIODONTICS RESTORATIVE DENT, 2012 OCT;32(5):581-9
- (11) ESPOSITO M, GRUSOVIN MG, LAMBERT F, MATOS S, PIETRUSKA M, ROSSI R, SALHI L, BUTI J  
**THE EFFECTIVENESS OF A RESORBABLE BONE SUBSTITUTE WITH A RESORBABLE MEMBRANE IN THE TREATMENT OF PERIODONTAL INFRA-BONY DEFECT - A MULTICENTER RANDOMISED CONTROLLED TRIAL**  
EUR J ORAL IMPLANTOL, 2015;8(3):233-244
- (12) ROMBOUITS C, JEANNEAU C, CAMILLERI J, LAURENT P, ABOUT I  
**CHARACTERIZATION AND ANGIOGENIC POTENTIAL OF XENOGENIC BONE GRAFTING MATERIALS: ROLE OF PERIODONTAL LIGAMENT CELLS**  
DENT MATER J, 2016 DEC 1;35(6):900-907

For further information see the complete literature on p. 92





# TSV Gel



***The resorbable solution for ideal graft stability***  
***Thermosensitive resorbable gel for graft stabilization***



# Characteristics and handling



### Composition

Heterologous type I and III collagen gel  
Thermogelling synthetic biocompatible copolymer

### Physical form

LV phase at <math>+4^{\circ}\text{C}</math>  
Gel viscosity at >math>+13^{\circ}\text{C}</math>

### Packaging

Syringe: 0.5 cc, 1.0 cc

Available only in combination with OsteoBiol® Gen-Os® and Apatos  
0.5 g, 1.0 g

### Product codes

TSV005S | 1 Syringe | 0.5 cc | Porcine  
TSV005E | 1 Syringe | 0.5 cc | Equine  
TSV010S | 1 Syringe | 1.0 cc | Porcine  
TSV010E | 1 Syringe | 1.0 cc | Equine

### GMDN code

46425

### CND code

P900402



## CHARACTERISTICS

The purpose of *TSV Gel* is to provide mechanical stability to bone substitutes and barrier membranes.


*TSV Gel* is sterilized by Gamma irradiation and is radio-transparent. It contains heterologous type I and III collagen gel with polyunsaturated fatty acids diluted in aqueous solution containing a biocompatible synthetic copolymer that gives *TSV Gel* thermo-reversible and thermo-gelling properties. At low temperature (+4°C) the gel is relatively flowable and easy to mix and manipulate with the graft but it becomes more viscous when *in situ* and exposed to body temperature.

## HANDLING

*TSV Gel* must be refrigerated for at least 20 minutes at +4°C before use, in order to reach the low viscosity (LV) phase, which makes it easier to mix with *Gen-Os®* or *Apatos*.

At room temperature, the product remains at LV phase for few minutes, whereas once *in situ* its viscosity quickly increases with body temperature. *TSV Gel* in LV phase can be used instead of saline for hydrating and mixing with *Gen-Os®* or *Apatos*. The result will be a sticky mixture easy to place and extremely stable once *in situ*.

*TSV Gel* can also be applied to the rough side of the *Evolution* membrane to stabilize it during graft covering and whilst suturing.



Part of a biopsy showing newly formed bone around a particle of OsteoBiol® Gen-Os® mixed with OsteoBiol® TSV Gel two weeks after grafting in rabbit. Htx-eosine.

Author: Prof Ulf Nannmark, University of Göteborg, Sweden



Source: Tecnos® Dental Media Library



Source: Tecnos® Dental Media Library

TSV Gel can be used in GBR procedures together with OsteoBiol® bone substitutes and membranes to enhance graft stability. The viscosity reached by TSV Gel at body temperature improves the stability of Gen-Os® or Apatos granules and has proven particularly beneficial in cases where there is little bony support around the defect i.e. lateral augmentation, sockets with a compromised buccal wall, dehiscences and periodontal two and one wall defects.

Additionally, the viscosity of OsteoBiol® TSV Gel improves the stability and handling of Evolution membranes, particularly during the delicate phase of flap closure.

The above clinical information is based on the experience of expert surgeons



**DEHISCENCES AND FENESTRATIONS**  
peri-implant lesions




**ALVEOLAR REGENERATION**  
socket preservation



**PERIODONTAL REGENERATION**  
intrabony defects




**HORIZONTAL AUGMENTATION**  
two-wall defects



**Peri-implant defect treated with OsteoBiol® Gen-Os® mixed with TSV Gel**

Author: Dr Roberto Rossi, Genova, Italy

## OsteoBiol® TSV Gel GELIFICATION KINETICS



Source: Politecnico di Torino, Italy


The graph shows the effect of temperature change on 3 TSV Gel samples.

As temperature increases from 0°C (1°C/min), the viscosity of the gel reaches its minimum at 4°C.


It then increases rapidly until it plateaus at 13°C. At room and body temperature TSV Gel is gel-like. It does not harden but keeps a soft consistency that allows the mixture with Gen-Os® or Apatos granules. Thanks to the hydrophilic properties of OsteoBiol® bone substitutes, the mixture becomes a sticky, stable conglomerate that can easily be placed in the defect site.

TSV Gel is biocompatible and rapidly resorbed.





# mp3<sup>®</sup>



**Ultimate performance and handling**

*Pre-hydrated collagenated heterologous cortico-cancellous bone mix*





#### Tissue of origin

Cortico-cancellous heterologous bone mix

#### Tissue collagen

Preserved plus an additional 10% collagen gel

#### Physical form

Pre-hydrated granules and collagen gel

#### Composition

90% granulated mix, 10% collagen gel

#### Granulometry

600-1000  $\mu\text{m}$

1000-2000  $\mu\text{m}$

#### Re-entry time

About 5 months

#### Packaging

Syringe: 0.5 cc, 1.0 cc, 3x0.25 cc, 3x0.5 cc, 3x1.0 cc

Wide tip syringe: 2.0 cc

#### Product codes

600-1000  $\mu\text{m}$

A3095FS | 1 Syringe | 0.5 cc | Porcine

A3095FE | 1 Syringe | 0.5 cc | Equine

A3005FS | 1 Syringe | 1.0 cc | Porcine

A3005FE | 1 Syringe | 1.0 cc | Equine

A3075FS | 3 Syringes | 3x0.25 cc | Porcine

A3015FS | 3 Syringes | 3x0.5 cc | Porcine

A3015FE | 3 Syringes | 3x0.5 cc | Equine

A3030FS | 3 Syringes | 3x1.0 cc | Porcine

A3030FE | 3 Syringes | 3x1.0 cc | Equine

A3010FS | 1 Wide tip syringe | 2.0 cc | Porcine

A3010FE | 1 Wide tip syringe | 2.0 cc | Equine

1000-2000  $\mu\text{m}$

A3210FS | 1 Wide tip syringe | 2.0 cc | Porcine

A3210FE | 1 Wide tip syringe | 2.0 cc | Equine

#### GMDN code

46425

#### CND code

P900402

## Characteristics and handling

### CHARACTERISTICS

Heterologous origin biomaterial made of 600-1000  $\mu\text{m}$  or 1000-2000  $\mu\text{m}$  pre-hydrated collagenated cortico-cancellous granules, properly mixed with collagen gel. Thus, it is possible both skipping the hydration phase and decreasing the risk of accidental exposure of the material to pathogens during manipulation and grafting phases; furthermore, the syringe is flexible and ideal to simplify grafting in the receiving site.

The granules are endowed with characteristics very similar to human mineral bone, and can be used as an alternative to autologous bone.

Their natural micro-porous consistency facilitates new bone tissue formation<sup>(1)</sup> in defect sites and accelerates the regeneration process.


Gradually resorbable<sup>(2,3)</sup>, it preserves the original graft shape and volume (osteoconductive property)<sup>(4,5)</sup>.

Moreover, thanks to its collagen content, the product facilitates blood clotting and the subsequent invasion of repairing and regenerative cells.

### HANDLING

*mp3*<sup>®</sup> is available in ready-to-use syringes and can be easily grafted avoiding the hydration and manipulation phases.

After adapting the material to the defect shape, it is necessary to remove non-stable residues before proceeding to soft tissue suture. It is recommended to always compact *mp3*<sup>®</sup> after grafting to achieve optimal stabilization.



**Histology on maxillary sinus biopsy taken at 24 months. 48% new bone formation, 13% residual granules**

Author: Biopsy by Dr Roberto Rossi, Genova, Italy. Histology by Prof Ulf Nannmark, University of Göteborg, Sweden



Source: Tecnos<sup>®</sup> Dental Media Library


mp3® is a pre-hydrated cortico-cancellous bone mix with 10% collagen gel. It has been developed with this innovative biotechnology and is a "ready-to-use" product.

mp3® is commonly used for lateral access maxillary sinus lift<sup>(1,6)</sup>, always in association with *Evolution* membranes, to cover the antrotomy: the mp3® syringe can be directly applied into the bony window without having to mix the mp3® granules with saline.


Due to its collagen gel content, mp3® allows an excellent graft stability while its hydrophilia guarantees quick blood absorption and therefore the necessary graft vascularization. mp3® has also been successfully used in combination with *Evolution* membranes for alveolar ridge preservation<sup>(3,7,8)</sup>: the application of this biomaterial limits the alveolar ridge width and height reduction that would naturally occur with spontaneous healing, preserving thus the alveolar ridge volume and allowing a correct second stage implant placement.

mp3® has been documented for horizontal augmentation (two wall defects) in combination with autogenous bone blocks or with OsteoBioL® Lamina<sup>(9,10)</sup>: its cortico-cancellous composition allows a progressive resorption of osteoclastic type, and in parallel a similar rate of new bone formation<sup>(2)</sup>.

These unique properties allow a very good graft volume preservation<sup>(11)</sup>, a healthy new bony tissue and ultimately, a successful implant rehabilitation.



Socket grafted with OsteoBioL® mp3®



Periodontal defect grafted with OsteoBioL® mp3®

Author: Dr Roberto Rossi, Genova, Italy

free animated videos  
on OsteoBioL® APP



**LATERAL ACCESS SINUS LIFT**  
maxillary sinus floor augmentation



**ALVEOLAR REGENERATION**  
post-extractive sockets



**HORIZONTAL AUGMENTATION**  
two-wall defects

## BIBLIOGRAPHY

(1) RAMIREZ FERNANDEZ MP, CALVO GUIRADO JL, MATÉ SANCHEZ DE VAL JE, DELGADO RUIZ RA, NEGRI B, BARONA DORADO C  
**ULTRASTRUCTURAL STUDY BY BACKSCATTERED ELECTRON IMAGING AND ELEMENTAL MICROANALYSIS OF BONE-TO-BIOMATERIAL INTERFACE AND MINERAL DEGRADATION OF PORCINE XENOGRAPTS USED IN MAXILLARY SINUS FLOOR ELEVATION**  
CLIN ORAL IMPLANTS RES, 2013 MAY;24(5):523-30

(2) NANNMARK U, SENNERBY L  
**THE BONE TISSUE RESPONSES TO PREHYDRATED AND COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS: A STUDY IN RABBIT MAXILLARY DEFECTS**  
CLIN IMPLANT DENT RELAT RES, 2008 DEC;10(4):264-70

(3) GIULIANI A, IEZZI G, MAZZONI S, PIATTELLI A, PERROTTI V, BARONE A  
**REGENERATIVE PROPERTIES OF COLLAGENATED PORCINE BONE GRAFTS IN HUMAN MAXILLA: DEMONSTRATIVE STUDY OF THE KINETICS BY SYNCHROTRON RADIATION MICROTOMOGRAPHY AND LIGHT MICROSCOPY**  
CLIN ORAL INVESTIG, 2017 2018 JAN;22(1):505-513

(4) SCARANO A, LORUSSO F, RAVERA L, MORTELLARO C, PIATTELLI A  
**BONE REGENERATION IN ILIAC CRESTAL DEFECTS: AN EXPERIMENTAL STUDY ON SHEEP**  
BIOMED RES INT, 2016;2016:4086870

(5) IEZZI G, PIATTELLI A, GIULIANI A, MANGANO C, BARONE A, MANZONI L, DEGIDI M, SCARANO A, FILIPPONE A, PERROTTI V  
**MOLECULAR, CELLULAR AND PHARMACEUTICAL ASPECTS OF FILLING BIOMATERIALS DURING MAXILLARY SINUS-LIFT PROCEDURES. PART 2: DETAILED CHARACTERISTICS OF THE MATERIALS**  
CURR PHARM BIOTECHNOL, 2017, 18, 33-44

(6) SILVESTRI M, MARTEGANI P, D'AVENIA F, FARNETI M, CAPRI D, PAOLANTONI G, LANDI L  
**SIMULTANEOUS SINUS AUGMENTATION WITH IMPLANT PLACEMENT: HISTOMORPHOMETRIC COMPARISON OF TWO DIFFERENT GRAFTING MATERIALS. A MULTICENTER DOUBLE-BLIND PROSPECTIVE RANDOMIZED CONTROLLED CLINICAL TRIAL**  
INT J ORAL MAXILLOFAC IMPLANTS, 2013 MAR-APR;28(2):543-9

(7) BARONE A, BORGIA V, COVANI U, RICCI M, PIATTELLI A, IEZZI G  
**FLAP VERSUS FLAPLESS PROCEDURE FOR RIDGE PRESERVATION IN ALVEOLAR EXTRACTION SOCKETS: A HISTOLOGICAL EVALUATION IN A RANDOMIZED CLINICAL TRIAL**  
CLIN ORAL IMPLANTS RES, 2015 JUL;26(7):806-13

(8) BARONE A, RICCI M, TONELLI P, SANTINI S, COVANI U  
**TISSUE CHANGES OF EXTRACTION SOCKETS IN HUMANS: A COMPARISON OF SPONTANEOUS HEALING VS. RIDGE PRESERVATION WITH SECONDARY SOFT TISSUE HEALING**  
CLIN ORAL IMPLANTS RES, 2013 NOV;24(11):1231-7



(9) WACHTEL H, FICKL S, HINZE M, BOLZ W, THALMAIR T  
**THE BONE LAMINA TECHNIQUE: A NOVEL APPROACH FOR LATERAL RIDGE AUGMENTATION - A CASE SERIES**  
INT J PERIODONTICS RESTORATIVE DENT, 2013 JUL-AUG;33(4):491-7

(10) ROSSI R, RANCITELLI D, POLI PP, RASIA DAL POLO M, NANNMARK U, MAIORANA C  
**THE USE OF A COLLAGENATED PORCINE CORTICAL LAMINA IN THE RECONSTRUCTION OF ALVEOLAR RIDGE DEFECTS. A CLINICAL AND HISTOLOGICAL STUDY**  
MINERVA STOMATOL, 2016 OCT;65(5):257-68

(11) BARONE A, TOTI P, MENCHINI-FABRIS GB, DERCHI G, MARCONCINI S, COVANI U  
**EXTRA ORAL DIGITAL SCANNING AND IMAGING SUPERIMPOSITION FOR VOLUME ANALYSIS OF BONE REMODELING AFTER TOOTH EXTRACTION WITH AND WITHOUT 2 TYPES OF PARTICULATE PORCINE MINERAL INSERTION: A RANDOMIZED CONTROLLED TRIAL**  
CLIN IMPLANT DENT RELAT RES, 2017 AUG;19(4):750-759

For further information see the complete literature on p. 92





# GTO®



***The new standard of excellence in biomaterials***

***Collagenated heterologous cortico-cancellous bone mix + TSV Gel***



# Characteristics and handling



### Tissue of origin

Cortico-cancellous heterologous bone mix

### Tissue collagen

Preserved

### Physical form

Pre-hydrated granules and TSV Gel

### Composition

80% granulated mix, 20% TSV Gel

### Granulometry

600-1000  $\mu\text{m}$

### Re-entry time

About 5 months

### Packaging

Syringe: 0.5 cc

Wide tip syringe: 2.0 cc

### Product codes

MU0005S | 1 Syringe | 0.5 cc | Porcine

MU0005E | 1 Syringe | 0.5 cc | Equine

MU0020S | 1 Wide tip syringe | 2.0 cc | Porcine

MU0020E | 1 Wide tip syringe | 2.0 cc | Equine

### GMDN code

46425

### CND code


P900402

## CHARACTERISTICS

Heterologous bone grafting material made of a mix of collagenated cortico-cancellous granules of size ranging from 600 to 1000  $\mu\text{m}$ , properly blended with TSV Gel, which is a mixture of heterologous type I and III collagen gel with polyunsaturated fatty acids and a biocompatible synthetic copolymer diluted in aqueous solution. GTO<sup>®</sup> is gradually resorbed and it is extremely osteoconductive. Moreover, the preserved collagen matrix characterizing the granules facilitates blood clotting and the subsequent invasion of repairing and regenerative cells. These unique properties guarantee an excellent rate of new bone formation, delivering adequate graft volume preservation, a healthy new bony tissue and, ultimately, a successful implant rehabilitation. The presence of the same kind of granules of its progenitor, mp3<sup>®</sup>, which are very similar to human mineral bone, assures a similar biological response of the host tissue. GTO<sup>®</sup> can be used as alternative to autologous bone.

## HANDLING

GTO<sup>®</sup> is available in two sizes (0.5 and 2.0 cc) as ready-to-use pre-hydrated biomaterial and can be easily grafted to the defect site. Thus, clinicians can skip the hydration step with saline or blood, saving time and decreasing the risk of accidental exposure to pathogens. The presence of TSV Gel ensures optimal stickiness of the material, which is also easily adaptable to the recipient site and extremely stable.



Part of a biopsy showing newly formed bone 4 months after treatment with OsteoBio<sup>®</sup> GTO<sup>®</sup>

Author: Prof Ulf Nannmark, University of Göteborg, Sweden





Source: Tecno<sup>®</sup> Dental Media Library

Author: Dr Patrick Palacci, Marseille, France




GTO® has been conceived as a universal biomaterial, easily adaptable to any bone defect, in association with *Evolution* membranes or *Lamina* to protect the graft. Nonetheless, thanks to its stickiness, has proved particularly effective for horizontal augmentation procedures (e.g.: two-walls defects, when the crest is resorbed) and for socket preservation cases with compromised buccal plate. GTO® can also be successfully used to treat peri-implant lesions. In case of open defects, GTO® should be grafted in consecutive layers compacting each layer with a sterile gauze.

The above clinical information is based on the experience of expert surgeons



**Peri-implant defect treated with OsteoBiol® GTO®**  
Author: Dr Patrick Palacci, Marseille, France



**DEHISCENCES AND FENESTRATIONS**  
peri-implant grafting




**ALVEOLAR REGENERATION**  
post-extractive sockets



**HORIZONTAL AUGMENTATION**  
two-wall defects

## OsteoBiol® TSV Gel GELIFICATION KINETICS



Source: Politecnico di Torino, Italy


The graph shows the effect of temperature change on 3 TSV Gel samples.

As temperature increases from 0°C (1°C/min), the viscosity of the gel reaches its minimum at 4°C.

It then increases rapidly until it plateaus at 13°C. At room and body temperature TSV Gel is gel-like. It does not harden but keeps a soft consistency that allows the mixture with Gen-Os® or Apatos granules. Thanks to the hydrophilic properties of OsteoBiol® bone substitutes, the mixture becomes a sticky, stable conglomerate that can easily be placed in the defect site.

TSV Gel is biocompatible and rapidly resorbed.





# Putty



***Engineered for peri-implant defects***  
*Pre-hydrated collagenated heterologous cortico-cancellous bone paste*



# Characteristics and handling



## Tissue of origin

Cortico-cancellous heterologous bone mix

## Tissue collagen

Preserved plus an additional 20% collagen gel

## Physical form

Plastic consistency composed of collagen gel loaded with 80% micronized bone mix

## Composition

80% granulated mix, 20% collagen gel

## Granulometry

Up to 300  $\mu\text{m}$

## Re-entry time

About 4 months

## Packaging

Syringe: 0.25 cc, 0.5 cc, 3x0.5 cc, 3x0.25 cc  
Wide tip syringe: 1.0 cc

## Product codes

HPT52S	1 Syringe	0.25 cc	Porcine
HPT09S	1 Syringe	0.5 cc	Porcine
HPT09E	1 Syringe	0.5 cc	Equine
HPT32S	3 Syringes	3x0.25 cc	Porcine
HPT32E	3 Syringes	3x0.25 cc	Equine
HPT35S	3 Syringes	3x0.5 cc	Porcine
HPT35E	3 Syringes	3x0.5 cc	Equine

HPT61S	1 Wide tip syringe	1.0 cc	Porcine
HPT61E	1 Wide tip syringe	1.0 cc	Equine

## GMDN code

46425

## CND code


P900402

## CHARACTERISTICS

*Putty* is a bone paste with at least 80% micronized heterologous bone (granulometry up to 300  $\mu\text{m}$ ) and collagen gel. It is made with an exclusive process that provides the product with exceptional malleability and plasticity, making it easy to apply into peri-implant defects with walls. Thanks to its collagen component, the product facilitates blood clotting and the subsequent invasion of repairing and regenerative cells, showing an osteoconductive behaviour<sup>(1)</sup>. Successful grafting needs complete stability of the biomaterial: for this reason *Putty* must be used only in cavities able to firmly contain it. Therefore, *Putty* must not be grafted in two wall defects or in lateral access sinus lift procedures.

## HANDLING

Inject the product and adapt it to defect morphology without compression; any non-stable residue must be removed before soft tissue suture. An *Evolution* membrane is recommended to protect *Putty* grafted in peri-implant defects.



Part of a biopsy showing newly formed bone after treatment with OsteoBiol® Putty


Author: Prof Ulf Nannmark, University of Göteborg, Sweden



Source: Tecnos® Dental Media Library



The extraordinary handling properties of *Putty* syringe make this product the ideal choice for self-contained peri-implant defects<sup>(2)</sup> and all small defects that present a self-contained cavity. Furthermore, the TecnoSS® manufacturing process avoids granules ceramization, allowing a progressive resorption of the biomaterial and, at the same time, an adequate new-bone formation rate<sup>(3)</sup>. *Putty*'s "soft" consistency also guarantees an easy and healthy soft-tissues healing. Thanks to these unique characteristics, *Putty* has been effectively used for peri-implant defects regeneration: following immediate post-extractive implants placement, *Putty* can be injected between the defect walls and the implant, guaranteeing a perfect filling of the entire defect volume<sup>(4)</sup>.



**SEM image of OsteoBio® Putty**  
Author: Prof Ulf Nannmark, University of Göteborg, Sweden

The product versatility also makes *Putty* the ideal solution when bone tissue has been lost due to peri-implantitis as long as the containing walls are present. In fact, the primary condition for gaining a successful regeneration is to achieve the biomaterial initial stability. Therefore, *Putty* must be used only in defects where the surrounding walls guarantee such condition: for example inside the bone crest when ridge-split technique is adopted<sup>(5)</sup>, or with horizontally resorbed crests, in association with OsteoBio® Lamina (Bone Layer technique)<sup>(6)</sup>.



**Peri-implant lesion grafted with OsteoBio® Putty**  
Author: Dr Roberto Rossi, Genova, Italy



**DEHISCENCES AND FENESTRATIONS**  
peri-implant defects




**HORIZONTAL AUGMENTATION**  
ridge split

## BIBLIOGRAPHY

- (1) ARCURI C, CECCHETTI F, GERMANO F, MOTTA A, SANTACROCE C  
**CLINICAL AND HISTOLOGICAL STUDY OF A XENOGENIC BONE SUBSTITUTE USED AS A FILLER IN POSTEXTRACTIVE ALVEOLUS**  
MINERVA STOMATOL, 2005 JUN;54(6):351-62
- (2) BARONE A, AMERI S, COVANI U  
**IMMEDIATE POSTEXTRACTION IMPLANTS: TREATMENT OF RESIDUAL PERI-IMPLANT DEFECTS. A RETROSPECTIVE ANALYSIS**  
EUR J IMPLANT PROSTHODONTICS, 2006,2: 99-106
- (3) NANNMARK U, AZARMEHR I  
**SHORT COMMUNICATION: COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS. A STUDY IN RABBIT MAXILLARY DEFECTS**  
CLIN IMPLANT DENT RELAT RES, 2010 JUN 1; 12(2):161-3
- (4) CASSETTA M, RICCI L, IEZZI G, DELL'AQUILA D, PIATTELLI A, PERROTTI V  
**RESONANCE FREQUENCY ANALYSIS OF IMPLANTS INSERTED WITH A SIMULTANEOUS GRAFTING PROCEDURE: A 5-YEAR FOLLOW-UP STUDY IN MAN**  
INT J PERIODONTICS RESTORATIVE DENT, 2012 OCT;32(5):581-9
- (5) SANTAGATA M, GUARINIELLO L, TARTARO G  
**A MODIFIED EDENTULOUS EXPANSION (MERE) TECHNIQUE FOR IMMEDIATE PLACEMENT OF IMPLANTS. A CASE REPORT**  
J ORAL IMPLANTOL, 2011 MAR;37 SPEC N.:114-9
- (6) LOPEZ MA, ANDREASI BASSI M, CONFALONE L, CARINCI F, ORMIANER Z, LAURITANO D  
**THE USE OF RESORBABLE CORTICAL LAMINA AND MICRONIZED COLLAGENATED BONE IN THE REGENERATION OF ATROPHIC CRESTAL RIDGES: A SURGICAL TECHNIQUE. CASE SERIES**  
J BIOL REGUL HOMEOST AGENTS, 2016 APR-JUN;30(2 SUPPL 1):81-85





# Gel 40



***A unique heterologous bone gel***  
***Collagenated heterologous cortico-cancellous bone mix***





## Characteristics and handling


### CHARACTERISTICS

Gel 40 is made of a collagen matrix (type I and III) obtained using an exclusive Tecnos<sup>®</sup> process, loaded for 60% of its volume with micronized heterologous bone (granulometry up to 300  $\mu\text{m}$ ). Thanks to its collagen component, Gel 40 facilitates the formation of primary blood clot and the subsequent invasion of repairing and regenerative cells; moreover, the cortico-cancellous component provides the necessary scaffold function.

The collagen gel component contained in Gel 40 is rapidly and totally resorbed; it is also endowed with exceptional anti-inflammatory, eutrophic and cicatrizing properties. This lipophilia is due mainly to a percentage of polyunsaturated fatty acids of the oleic-linoleic series (to which Omega 3 also belongs) directly derived from the raw material. Such components possess a valuable antioxidant action on the free radicals and therefore aid tissue regeneration.

### HANDLING

The distinctive characteristics of viscosity and density of Gel 40 facilitate the handling of the product by the operator, providing a glue-like support. If viscosity is excessive, add a few drops of sterile lukewarm saline and then re-mix thoroughly to obtain the desired density.



Part of a biopsy showing newly formed bone after treatment with OsteoBiol<sup>®</sup> Gel 40. Biopsies were taken 5 weeks after implantation in rabbit maxillae. Htx-eosine. Original magnification x20  
Author: Prof Ulf Nannmark, University of Göteborg, Sweden



Source: Tecnos<sup>®</sup> Dental Media Library

#### Tissue of origin

Cortico-cancellous heterologous bone mix

#### Tissue collagen

Preserved plus an additional 40% collagen gel

#### Physical form

Collagen gel type I and III loaded with 60% bone mix

#### Composition

60% granulated mix, 40% collagen gel

#### Granulometry

Up to 300  $\mu\text{m}$

#### Re-entry time

About 4 months

#### Packaging

Syringe: 0.5 cc, 3x0.5 cc

#### Product codes

05GEL40S | 1 Syringe | 0.5 cc | Porcine  
05GEL40E | 1 Syringe | 0.5 cc | Equine  
15GEL40S | 3 Syringes | 3x0.5 cc | Porcine  
15GEL40E | 3 Syringes | 3x0.5 cc | Equine

#### GMDN code

46425

#### CND code

P900402



The exclusive Tecnos® manufacturing process guarantees an exceptional malleability and plasticity: furthermore, the syringe packaging provides Gel 40 extraordinary handling properties making this product the ideal choice for crestal access sinus lift<sup>(1,2)</sup>, deep and narrow peri-implant defects<sup>(3)</sup>, three-wall intrabony defects and, in combination with Evolution membranes, for treating gingival recessions<sup>(4)</sup>.

Furthermore, the Tecnos® manufacturing process avoids granules ceramization, allowing a progressive resorption of the biomaterial and, at the same time, an adequate new-bone formation rate<sup>(5,6)</sup>.


Gel 40 "soft" consistency also guarantees an easy and healthy soft-tissues healing.



Crestal access sinus lift with OsteoBiol® Gel 40  
Source: Tecnos® Dental Media Library



**PERIODONTAL REGENERATION**  
intrabony defects and gingival recessions




**CRESTAL ACCESS SINUS LIFT**  
crestal sinus floor augmentation

## BIBLIOGRAPHY

- (1) BARONE A, CORNELINI R, CIAGLIA R, COVANI U  
**IMPLANT PLACEMENT IN FRESH EXTRACTION SOCKETS AND SIMULTANEOUS OSTEOTOME SINUS FLOOR ELEVATION: A CASE SERIES**  
INT J PERIODONTICS RESTORATIVE DENT, 2008 JUN; 28(3):283-9
- (2) SANTAGATA M, GUARINIELLO L, RAUSO R, TARTARO G  
**IMMEDIATE LOADING OF DENTAL IMPLANT AFTER SINUS FLOOR ELEVATION WITH OSTEOTOME TECHNIQUE: A CLINICAL REPORT AND PRELIMINARY RADIOGRAPHIC RESULTS**  
J ORAL IMPLANTOL, 2010; 36(6):485-489
- (3) COVANI U, CORNELINI R, BARONE A  
**BUCCAL BONE AUGMENTATION AROUND IMMEDIATE IMPLANTS WITH AND WITHOUT FLAP ELEVATION: A MODIFIED APPROACH**  
INT J ORAL MAXILLOFAC IMPLANTS, 2008 SEP-OCT; 23(5):841-6
- (4) CARDAROPOLI D, CARDAROPOLI G  
**HEALING OF GINGIVAL RECESSIONS USING A COLLAGEN MEMBRANE WITH A DEMINERALIZED XENOGRIFT: A RANDOMIZED CONTROLLED CLINICAL TRIAL**  
INT J PERIODONTICS RESTORATIVE DENT, 2009 FEB; 29(1):59-67
- (5) NANNMARK U, AZARMEHR I  
**SHORT COMMUNICATION: COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS. A STUDY IN RABBIT MAXILLARY DEFECTS**  
CLIN IMPLANT DENT RELAT RES, 2010 JUN 1; 12(2):161-3
- (6) LORENZON G, BUTTARELLO GM, CHESSA G  
**CASE REPORT: IMPLANT PLACEMENT AND IMMEDIATE LOADING WITH SIMULTANEOUS BONE REGENERATION FOLLOWING JAW ODONTOGENIC CYST ENUCLEATION**  
DENTISTRY, 2015, 5:2

For further information see the complete literature on p. 92





# Apatos



***Microcrystalline hydroxyapatite***  
*Heterologous cortico-cancellous and cortical bone*



# Characteristics and handling



## Tissue of origin

Apatos Mix: cortico-cancellous heterologous bone mix  
Apatos Cortical: heterologous cortical bone

## Tissue collagen

Degraded

## Physical form

Radiopaque granules of mineral hydroxyapatite

## Composition

Apatos Mix: 100% cortico-cancellous mix  
Apatos Cortical: 100% cortical bone

## Granulometry

600-1000  $\mu\text{m}$   
1000-2000  $\mu\text{m}$

## Re-entry time

About 5 months

## Packaging

Mix | Vial: 0.5 g, 1.0 g, 2.0 g  
Cortical | Vial: 0.5 g, 1.0 g

## Product codes

600-1000  $\mu\text{m}$   
Mix | A1005FS | 1 Vial | 0.5 g | Porcine  
Mix | A1005FE | 1 Vial | 0.5 g | Equine  
Mix | A1010FS | 1 Vial | 1.0 g | Porcine  
Mix | A1010FE | 1 Vial | 1.0 g | Equine  
Mix | A1020FS | 1 Vial | 2.0 g | Porcine  
Mix | A1020FE | 1 Vial | 2.0 g | Equine  
Cortical | AC1005FS | 1 Vial | 0.5 g | Porcine  
Cortical | AC1010FS | 1 Vial | 1.0 g | Porcine  
1000-2000  $\mu\text{m}$   
Mix | A0210FS | 1 Vial | 1.0 g | Porcine  
Mix | A0210FE | 1 Vial | 1.0 g | Equine

## GMDN code

46425

## CND code


P900402

## CHARACTERISTICS

Apatos is a biocompatible<sup>(1,2)</sup>, osteoconductive<sup>(3,4)</sup> biomaterial of heterologous origin with characteristics similar to mineralized human bone<sup>(5)</sup>; it can therefore be used as an alternative to autologous bone. The natural microporous consistency of Apatos facilitates the formation of new bone tissue in bone defect area<sup>(6)</sup>, accelerating the process. Apatos microcrystalline hydroxyapatite is available in cortical and mixed granules.

## HANDLING

Apatos must always be hydrated and thoroughly mixed with a few drops of sterile saline or with *TSV Gel* to increase graft stability in not self-contained defects; it can also be mixed with patient's blood. Finally it can be mixed if necessary with the drug selected for surgery; the mixture thus obtained should be positioned with a sterile spatula or syringe for biomaterials.



SEM image of OsteoBiol® Apatos, cancellous granules

Source: Nobil Bio Ricerche, Villafranca d'Asti, Italy



Source: Tecnos® Dental Media Library




Apatos is a universal filler, that can be used to treat peri-implant defects and two-wall defects<sup>(7,8)</sup>. Because of its granulometry, Apatos cannot be used in narrow defects, but it fits well in big sockets, e.g. after molar extractions<sup>(9)</sup>. Both sinus lift procedures (with crestal or lateral access)<sup>(2,10)</sup> can be performed with Apatos as bone substitute, as well as surgeries for horizontal regenerations.

Apatos Cortical is characterized by a very long resorption time<sup>(11)</sup>, guaranteeing adequate preservation of the grafted volume.

When needed, Apatos grafts can be protected with OsteoBiol® Evolution membrane<sup>(12)</sup> or stabilized with Cortical Lamina.

The above clinical information is based on the experience of expert surgeons




**LATERAL ACCESS SINUS LIFT**  
maxillary sinus floor augmentation




**ALVEOLAR REGENERATION**  
socket preservation



**HORIZONTAL AUGMENTATION**  
two-wall defects



**DEHISCENCES AND FENESTRATIONS**  
peri-implant grafting



**CRESTAL ACCESS SINUS LIFT**  
osteotome sinus floor augmentation






## BIBLIOGRAPHY

- (1) TRUBIANI O, SCARANO A, ORSINI G, DI IORIO D, D'ARCANGELO C, PICCIRILLI M, SIGISMONDO M, CAPUTI S  
**THE PERFORMANCE OF HUMAN PERIODONTAL LIGAMENT MESENCHYMAL STEM CELLS ON XENOGENIC BIOMATERIALS**  
INT J IMMUNOPATHOL PHARMACOL, 2007 JAN-MAR; 20(1 SUPPL 1):87-91
- (2)ORSINI G, SCARANO A, PIATTELLI M, PICCIRILLI M, CAPUTI S, PIATTELLI A  
**HISTOLOGIC AND ULTRASTRUCTURAL ANALYSIS OF REGENERATED BONE IN MAXILLARY SINUS AUGMENTATION USING A PORCINE BONE-DERIVED BIOMATERIAL**  
J PERIODONTOL, 2006 DEC;77(12):1984-90
- (3) BRUNELLI G, SOLLAZZO V, CARINCI F, PALMIERI A, GIRARDI A, MONGUZZI R  
**OSTEOBIOL® INFLUENCES OSTEOGENIC DIFFERENTIATION OF ADIPOSE DERIVED STEM CELLS**  
EUR J INFLAMM, 2011, VOL. 9, NO. 3(S), 103-107
- (4) CAKIR M, KARACA IR, AYSEGÜL F, KAYMAZ F, BOZKAYA S  
**EXPERIMENTAL EVALUATION OF THE EFFECTS OF ANKAFERD BLOOD STOPPER AND COLLAGENATED HETEROLOGOUS BONE GRAFT ON BONE HEALING IN SINUS FLOOR AUGMENTATION**  
CLIN ORAL IMPLANTS RES, 2015 MAR-APR;30(2):279-85
- (5) KOLMAS J, SZWAJA M, KOLODZIEJSKI W  
**SOLID-STATE NMR AND IR CHARACTERIZATION OF COMMERCIAL XENOGENIC BIOMATERIALS USED AS BONE SUBSTITUTES**  
J PHARM BIOMED ANAL, 2012 MAR 5;61:136-41
- (6) BARONE A, TOTI P, QUARANTA A, ALFONSI F, CUCCHI A, NEGRI B, DI FELICE R, MARCHIONNI S, CALVO GUIRADO JL, COVANI U, NANNMARK U  
**CLINICAL AND HISTOLOGICAL CHANGES AFTER RIDGE PRESERVATION WITH TWO XENOGRAFTS: PRELIMINARY RESULTS FROM A MULTICENTER RANDOMIZED CONTROLLED CLINICAL TRIAL**  
J CLIN PERIODONTOL, 2017 FEB;44(2):204-214
- (7) BARONE A, AMERI S, COVANI U  
**IMMEDIATE POSTEXTRACTION IMPLANTS: TREATMENT OF RESIDUAL PERI-IMPLANT DEFECTS. A RETROSPECTIVE ANALYSIS**  
EUR J IMPLANT PROSTHODONTICS, 2006,2: 99-106
- (8) BARONE A, TOTI P, QUARANTA A, DERCHI G, COVANI U  
**THE CLINICAL OUTCOMES OF IMMEDIATE VERSUS DELAYED RESTORATION PROCEDURES ON IMMEDIATE IMPLANTS: A COMPARATIVE COHORT STUDY FOR SINGLE-TOOTH REPLACEMENT**  
CLIN IMPLANT DENT RELAT RES, 2015 DEC;17(6):1114-26
- (9) BARONE A, TOTI P, QUARANTA A, ALFONSI F, CUCCHI A, CALVO GUIRADO JL, NEGRI B, DI FELICE R, COVANI U  
**VOLUMETRIC ANALYSIS OF REMODELLING PATTERN AFTER RIDGE PRESERVATION COMPARING USE OF TWO TYPES OF XENOGRAFTS. A MULTICENTRE RANDOMIZED CLINICAL TRIAL**  
CLIN IMPLANT DENT RELAT RES, 2015 DEC;17(6):1114-26
- (10) IEZZI G, DEGIDI M, PIATTELLI A, MANGANO C, SCARANO A, SHIBLI JA, PERROTTI V  
**COMPARATIVE HISTOLOGICAL RESULTS OF DIFFERENT BIOMATERIALS USED IN SINUS AUGMENTATION PROCEDURES: A HUMAN STUDY AT 6 MONTHS**  
CLIN ORAL IMPLANTS RES, 2012 DEC;23(12):1369-76
- (11) SCARANO A, PIATTELLI A, PERROTTI V, MANZON L, IEZZI G  
**MAXILLARY SINUS AUGMENTATION IN HUMANS USING CORTICAL PORCINE BONE: A HISTOLOGICAL AND HISTOMORPHOMETRICAL EVALUATION AFTER 4 AND 6 MONTHS**  
CLIN IMPLANT DENT RELAT RES, 2011 MAR; 13(1):13-18
- 12 | MARCONCINI S, GIAMMARINARO E, DERCHI G, ALFONSI F, COVANI U, BARONE A  
**CLINICAL OUTCOMES OF IMPLANTS PLACED IN RIDGE-PRESERVED VERSUS NONPRESERVED SITES: A 4-YEAR RANDOMIZED CLINICAL TRIAL**  
CLIN IMPL DENT RELAT RES, 2018 Dec;20(6):906-914

For further information see the complete literature on p. 92




# BLOCKS



# OsteoBiol® bone blocks

## Sp-Block


collagenated cancellous bone



SEM image of OsteoBiol® Sp-Block. Magnification 200x.  
Source: Politecnico di Torino, Italy  
For more information on OsteoBiol® Sp-Block see page 54


## Dual-Block

collagenated cortico-cancellous bone



SEM image of OsteoBiol® Dual-Block. Magnification 20x.  
Source: Politecnico di Torino, Italy  
For more information on OsteoBiol® Dual-Block see page 54






# Sp-Block

**Cancellous block for the inlay technique in the mandible**



**Highly osteoconductive properties**



# Dual-Block

**Cortico-cancellous scaffold for horizontal augmentation in the maxilla**





# Characteristics, handling and clinical information

free animated videos  
on OsteoBiol® APP



## CHARACTERISTICS


*Sp-Block* is a cancellous block of xenogenic bone produced with an exclusive Tecnos® process which avoids ceramization of the hydroxyapatite crystals, thus accelerating physiological resorption. *Sp-Block* supports new bone formation<sup>(1)</sup>: thanks to its rigid consistency it is able to maintain the original graft volume, which is particularly important in case of large regenerations. Moreover, its collagen content facilitates blood clotting and the subsequent invasion of regenerative and repairing cells, favoring the *restitutio ad integrum* of missing bone.

## HANDLING

*Sp-Block* must be hydrated before use for 5/10 minutes with sterile lukewarm physiological solution or with antibiotics. Afterwards, it can be adapted to the receiving site; the block must always be fixed with osteosynthesis microscrews and should be protected with a resorbable membrane (*Evolution*).

## CLINICAL INFORMATION

*Sp-Block* has been documented in cases where a vertical gain in the posterior mandible is required, to achieve an augmentation of maximum 5 mm, by means of the inlay technique<sup>(2-4)</sup>. The gaps around the block can be filled with a biomaterial in granules; the augmented area is stabilized with mini-plates and screws and covered with an *Evolution* membrane.



SEM image of OsteoBiol® cancellous block  
Author: Prof Ulf Nannmark, University of Göteborg, Sweden

### Tissue of origin

Cancellous bone

### Tissue collagen

Preserved

### Physical form

Rigid dried block

### Composition

Collagenated cancellous bone

### Re-entry time

About 8 months, variable depending on characteristics and irradiation grade of grafting site and on clinical conditions of the patient

### Packaging

Sterile blister

### Product codes

BNOE | 10x10x10 mm | Equine  
BN1E | 10x10x20 mm | Equine  
BN2E | 10x20x20 mm | Equine  
BN8E | 35x10x5 mm | Equine

### GMDN code

46425

### CND code

P900402

## BIBLIOGRAPHY

(1) SCARANO A, LORUSSO F, RAVERA L, MORTELLARO C, PIATTELLI A  
BONE REGENERATION IN ILIAC CRESTAL DEFECTS: AN EXPERIMENTAL STUDY ON SHEEP  
BIOMED RES INT, 2016;2016:4086870

(2) ESPOSITO M, BARAUSSE C, PISTILLI R, PIATTELLI M, DI SIMONE S, IPPOLITO DR, FELICE P

POSTERIOR ATROPHIC JAWS REHABILITATED WITH PROSTHESES SUPPORTED BY 5 X 5 MM IMPLANTS WITH A NOVEL NANOSTRUCTURED CALCIUM-INCORPORATED TITANIUM SURFACE OR BY LONGER IMPLANTS IN AUGMENTED BONE. FIVE-YEAR RESULTS FROM A RANDOMISED CONTROLLED TRIAL  
INT J OF ORAL IMPLANTOL, 2019;12(1):39-54

(3) FELICE P, BARAUSSE C, BARONE A, ZUCHELLI G, PIATTELLI M, PISTILLI R, IPPOLITO DR, SIMION M  
INTERPOSITIONAL AUGMENTATION TECHNIQUE IN THE TREATMENT OF POSTERIOR MANDIBULAR ATROPHIES: A RETROSPECTIVE STUDY COMPARING 129 AUTOGENOUS AND HETEROLOGOUS BONE BLOCKS WITH 2 TO 7 YEARS FOLLOW-UP  
INT J PERIODONTICS RESTORATIVE DENT, 2017 JUL/AUG;37(4):469-480

(4) BARONE A, TOTI P, MENCHINI FABRIS GB, MARCHIONNI S, COVANI U  
EARLY VOLUMETRIC CHANGES AFTER VERTICAL AUGMENTATION OF THE ATROPHIC POSTERIOR MANDIBLE WITH INTERPOSITIONAL BLOCK GRAFT VERSUS ONLAY BONE GRAFT: A RETROSPECTIVE RADIOLOGICAL STUDY  
J CRANIO-MAXILLOFAC, 2017 SEP;45(9):1438-1447




VERTICAL AUGMENTATION  
inlay technique

## CHARACTERISTICS

*Dual-Block* is a cortico-cancellous block of xenogenic bone with osteoconductive characteristics. It can be used when the regeneration of big volumes is needed: thanks to the collagen content that promotes blood clotting and migration of regenerative and repairing cells<sup>(1)</sup>, the graft offers an adequate support for tissue reconstruction and is gradually resorbed, while new bone is produced by osteoblasts.

## HANDLING

*Dual-Block* must be hydrated before use with sterile lukewarm physiological solution or with antibiotics (5/10 minutes for Soft version; up to 40 minutes for Norm version). Afterwards, the block can be adapted to the receiving site which must be accurately decorticated in order to guarantee maximum contact; the block should always be fixed with osteosynthesis microscrews and protected with *Evolution* membrane.



SEM image of OsteoBiol® *Dual-Block*  
Source: Politecnico di Torino, Italy

## CLINICAL INFORMATION

*Dual-Block* can be grafted with the onlay technique only to augment horizontally heavily resorbed maxilla. The gaps around the block can be filled with a biomaterial in granules to achieve the desired volume and contour of the augmented recipient site.


*The above clinical information is based on the experience of expert surgeons*

## BIBLIOGRAPHY


(1) MANESCU A, GIULIANI A, MOHAMMADI S, TROMBA G, MAZZONI S, DIOMEDE F, ZINI N, PIATTELLI A, TRUBIANI O  
OSTEOGENIC POTENTIAL OF DUAL-BLOCKS CULTURED WITH HUMAN PERIODONTAL LIGAMENT STEM CELLS: IN VITRO AND SYNCHROTRON  
J PERIODONTAL RES, 2016 Feb;51(1):112-24



OsteoBiol® *Dual-Block*  
Source: TecnoSS® Dental Media Library



**HORIZONTAL AUGMENTATION**  
onlay technique



### Tissue of origin

Cortico-cancellous bone

### Tissue collagen

Preserved

### Physical form

Rigid dried block

### Composition

Collagenated cortico-cancellous bone

### Re-entry time

About 8 months, variable depending on characteristics and irradiation grade of grafting site and on clinical conditions of the patient

### Packaging

Sterile blister

### Product codes

STS7S | 20x15x5 mm | Soft | Porcine curved  
STN5S | 20x10x5 mm | Norm | Porcine curved

### GMDN code


46425

### CND code

P900402




# Histological results in alveolar regeneration



- **A)** Barone A et al. - **Flap versus flapless procedure for ridge preservation in alveolar extraction sockets: a histological evaluation in a randomized clinical trial**  
Clinical Oral Implants Research, 2015 Jul;26(7):806-13
- **B)** Giuliani A et al. - **Regenerative properties of collagenated porcine bone grafts in human maxilla: demonstrative study of the kinetics by synchrotron radiation microtomography and light microscopy**  
Clinical Oral Investigations, 2018 Jan;22(1):505-513
- **C)** Barone A et al. - **Clinical and histological changes after ridge preservation with two xenografts: preliminary results from a multicenter randomized controlled clinical trial**  
Journal of Clinical Periodontology, 2017 Feb;44(2):204-214
- **D)** Crespi R et al. - **Corticancellous porcine bone in the healing of human extraction sockets: combining histomorphometry with osteoblast gene expression profiles in vivo**  
Int Journal of Oral and Maxillofacial Implants, 2011 Jul - Aug; 26(4):866-72
- **E)** Crespi R et al. - **Comparison of magnesium-enriched hydroxyapatite and porcine bone in human extraction socket healing: a histologic and histomorphometric evaluation**  
Int Journal of Oral and Maxillofacial Implants, 2011 Sep-Oct;26(5):1057-62
- **F)** Barone A et al. - **Xenograft versus extraction alone for ridge preservation after tooth removal: a clinical and histomorphometric study**  
Journal of Periodontology, 2008 Aug; 79(8):1370-7

# Histological results in sinus lift



- **A)** Cassetta M et al. - **Bone formation in sinus augmentation procedures using autologous bone, porcine bone, and a 50 : 50 mixture: a human clinical and histological evaluation at 2 months**  
Clinical Oral Implants Research, 2015 Oct; 26(10):1180-4
- **B)** Scarano A et al. - **Maxillary sinus augmentation in humans using cortical porcine bone: a histological and histomorphometrical evaluation after 4 and 6 months**  
Clinical Implant Dentistry and Related Research, 2011 Mar; 13(1):13-18
- **C)** Orsini G et al. - **Histologic and ultrastructural analysis of regenerated bone in maxillary sinus augmentation using a porcine bone-derived biomaterial**  
Journal of Periodontology, 2006 Dec; 77(12):1984-90
- **D)** Iezzi G et al. - **Comparative histological results of different biomaterials used in sinus augmentation procedures: a human study at 6 months**  
Clinical Oral Implants Research, 2012 Dec;23(12):1369-76
- **E)** Tanaka K et al. - **Sinus floor elevation and antrostomy healing: a histomorphometric clinical study in humans**  
Implant dentistry, 2019 Dec; 28(6):537-542



# MEMBRANES AND BARRIERS



Fig. 1




Fig. 2




Fig. 3



Fig. 4




Fig. 5




Fig. 6




Fig. 7




Fig. 8

# OsteoBiol® membranes and barriers

## MEMBRANES

## BARRIERS

### Evolution

*Heterologous  
mesenchymal tissue*



**Dried membrane with one  
smooth side and one  
micro-rough side**

For more information on OsteoBiol® Evolution  
see page 62

### Derma

*Porcine derma*



**Dried  
membrane**

For more information on OsteoBiol® Derma  
see page 66

### Special

*Heterologous  
pericardium*



**Translucent  
dried membrane**

For more information on OsteoBiol® Special  
see page 74

### Lamina

*Cortical bone*




**Semi-rigid and rigid  
dried lamina**


For more information on OsteoBiol® Lamina  
see page 70

SEM image showing collagenic matrix of OsteoBiol® membranes  
Source: Courtesy of Nobil Bio Ricerche, Villafranca D'Asti, Italy





# *Evolution*



***The natural Evolution of collagen membranes***  
*Heterologous mesenchymal tissue*



# Characteristics and handling



### Tissue of origin

Heterologous mesenchymal tissue

### Tissue collagen

Preserved

### Physical form

Dried membrane with one smooth side and one micro-rough side

### Thickness

X-Fine: 0.2 mm  
Fine: 0.3 mm  
Standard: 0.4 mm

### Estimated resorption time

X-Fine: about 2 months  
Fine: about 3 months  
Standard: about 4 months

### Size

20x20 mm, 30x30 mm, 25x35 mm (oval), 40x40 mm, 80x60 mm

### Product codes

EM33XS | 3 pcs | 30x30 mm | X-Fine | Porcine  
EV02LLE | 20x20 mm | Fine | Equine  
EV03LLE | 30x30 mm | Fine | Equine  
EVOLLE | 25x35 mm (oval) | Fine | Equine  
EV04LLE | 40x40 mm | Fine | Equine  
EV06LLE | 80x60 mm | Fine | Equine  
EM02HS | 20x20 mm | Standard | Porcine  
EV02HHE | 20x20 mm | Standard | Equine  
EM03HS | 30x30 mm | Standard | Porcine  
EV03HHE | 30x30 mm | Standard | Equine  
EM00HS | 25x35 mm (oval) | Standard | Porcine

### GMDN code

47184

### CND code

P900402

## CHARACTERISTICS


Obtained from heterologous mesenchymal tissue, the *Evolution* membrane is gradually resorbable<sup>(1)</sup>. Its structure is made of dense collagen fibers of high consistency and of extraordinary resistance that offer the specialist surgeon:

- maximum adaptability to bone tissue and soft tissues
- easy and secure suturability to nearby tissues
- best membrane-bone and membrane-periosteum interface
- stability and prolonged protection of the underlying graft
- clot stabilization and isolation<sup>(2)</sup>

## HANDLING


The membrane can be shaped with sterile scissors until the desired size is reached; unless the grafting site is already bleeding, the membrane should be rehydrated with lukewarm physiological solution. Once it acquires the desired plasticity, it must be adapted to the grafting site.

NB: in case of accidental exposure, the dense collagenic matrix of *Evolution* protects the graft from infection; the membrane itself will also not be infected, allowing second intention healing<sup>(3-5)</sup>.



SEM image of an OsteoBiol® Evolution standard membrane

Source: Politecnico di Torino, Italy



Source: Tecnos® Dental Media Library




*Evolution* is obtained from heterologous mesenchymal tissue and is completely resorbable. Experimental studies have shown histological evidence of the prolonged barrier effect of this membrane, which lasts at least eight weeks<sup>(1)</sup>, protecting the graft from external agents.

This property is particularly important in case of flapless regeneration<sup>(3)</sup> of large posterior sockets<sup>(5)</sup>: in these cases, the standard model has proved to be the most effective.

In lateral access sinus lift, *Evolution* membranes have been documented for antrostomy coverage (standard model)<sup>(6,7)</sup> and for protection of the sinus membrane from cutting risk due to graft pressure (fine model)<sup>(8)</sup>.


*Evolution* can be used to protect peri-implant regenerations<sup>(9)</sup> and periodontal grafts<sup>(10)</sup>. Furthermore, *Evolution* fine has been successfully used to protect *Sp-Block* in vertical augmentation with the inlay technique<sup>(11)</sup>.




**LATERAL ACCESS SINUS LIFT**  
maxillary sinus floor augmentation



**PERIODONTAL REGENERATION**  
intrabony defects



**HORIZONTAL AUGMENTATION**  
two-wall defects



**DEHISCENCES AND FENESTRATIONS**  
peri-implant lesions



**ALVEOLAR REGENERATION**  
graft protection




**VERTICAL AUGMENTATION**  
inlay technique

## BIBLIOGRAPHY


- (1) NANNMARK U, SENNERBY L  
THE BONE TISSUE RESPONSES TO PREHYDRATED AND COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS: A STUDY IN RABBIT MAXILLARY DEFECTS  
CLIN IMPLANT DENT RELAT RES, 2008 DEC;10(4):264-70
- (2) KILINC A, ATAOL M  
HOW EFFECTIVE IS COLLAGEN RESORBABLE MEMBRANE PLACEMENT AFTER PARTIALLY IMPACTED MANDIBULAR THIRD MOLAR SURGERY ON POSTOPERATIVE MORBIDITY? A PROSPECTIVE RANDOMIZED COMPARATIVE STUDY  
BMC ORAL HEALTH, 2017 OCT 5;17(1):126
- (3) BARONE A, BORGIA V, COVANI U, RICCI M, PIATTELLI A, IEZZI G  
FLAP VERSUS FLAPLESS PROCEDURE FOR RIDGE PRESERVATION IN ALVEOLAR EXTRACTION SOCKETS: A HISTOLOGICAL EVALUATION IN A RANDOMIZED CLINICAL TRIAL  
CLIN ORAL IMPLANTS RES, 2015 JUL;26(7):806-13
- (4) BARONE A, RICCI M, TONELLI P, SANTINI S, COVANI U  
TISSUE CHANGES OF EXTRACTION SOCKETS IN HUMANS: A COMPARISON OF SPONTANEOUS HEALING VS. RIDGE PRESERVATION WITH SECONDARY SOFT TISSUE HEALING  
CLIN ORAL IMPLANTS RES, 2013 NOV;24(11):1231-7
- (5) GIULIANI A, IEZZI G, MAZZONI S, PIATTELLI A, PERROTTI V, BARONE A  
REGENERATIVE PROPERTIES OF COLLAGENATED PORCINE BONE GRAFTS IN HUMAN MAXILLA: DEMONSTRATIVE STUDY OF THE KINETICS BY SYNCHROTRON RADIATION MICROTOMOGRAPHY AND LIGHT MICROSCOPY  
CLIN ORAL INVEST, 2017 2018 JAN;22(1):505-513
- (6) BARONE A, RICCI M, GRASSI RF, NANNMARK U, QUARANTA A, COVANI U  
A 6-MONTH HISTOLOGICAL ANALYSIS ON MAXILLARY SINUS AUGMENTATION WITH AND WITHOUT USE OF COLLAGEN MEMBRANES OVER THE OSTEOTOMY WINDOW: RANDOMIZED CLINICAL TRIAL  
CLIN ORAL IMPLANTS RES, 2013 JAN;24(1):1-6
- (7) SCARANO A, PIATTELLI A, PERROTTI V, MANZON L, IEZZI G  
MAXILLARY SINUS AUGMENTATION IN HUMANS USING CORTICAL PORCINE BONE: A HISTOLOGICAL AND HISTOMORPHOMETRICAL EVALUATION AFTER 4 AND 6 MONTHS  
CLIN IMPLANT DENT RELAT RES, 2011 MAR; 13(1):13-18
- (8) CASSETTA M, RICCI L, IEZZI G, CALASSO S, PIATTELLI A, PERROTTI V  
USE OF PIEZOSURGERY DURING MAXILLARY SINUS ELEVATION: CLINICAL RESULTS OF 40 CONSECUTIVE CASES  
INT J PERIODONTICS RESTORATIVE DENT, 2012 DEC;32(6):E182-8
- (9) BARONE A, MARCONCINI S, GIAMMARINARO E, MIJIRITSKY E, GELPI F, COVANI U  
CLINICAL OUTCOMES OF IMPLANTS PLACED IN EXTRACTION SOCKETS AND IMMEDIATELY RESTORED: A 7-YEAR SINGLE-COHORT PROSPECTIVE STUDY  
CLIN IMPLANT DENT RELAT RES, 2016 DEC;18(6):1103-1112
- (10) ESPOSITO M, GRUSOVIN MG, LAMBERT F, MATOS S, PIETRUSKA M, ROSSI R, SALHI L, BUTI J  
THE EFFECTIVENESS OF A RESORBABLE BONE SUBSTITUTE WITH A RESORBABLE MEMBRANE IN THE TREATMENT OF PERIODONTAL INFRABONY DEFECT - A MULTICENTER RANDOMISED CONTROLLED TRIAL  
EUR J ORAL IMPLANTOL, 2015;8(3):233-244
- (11) FELICE P, PIANA L, CHECCHI L, CORVINO V, NANNMARK U, PIATTELLI M  
VERTICAL RIDGE AUGMENTATION OF ATROPHIC POSTERIOR MANDIBLE WITH AN INLAY TECHNIQUE AND CANCELLOUS EQUINE BONE BLOCK: A CASE REPORT  
INT J PERIODONTICS RESTORATIVE DENT, 2013 MAR;33(2):159-66

For further information see the complete literature on p. 92





# Derma



**A xenogenic matrix for soft tissue augmentation**  
**Collagen dermal membrane**



# Characteristics and handling



## Tissue of origin

Porcine derma

## Tissue collagen

Preserved

## Physical form

Dried membrane

## Composition

100% derma

## Thickness

0.5 mm ( $\pm 0.1$  mm)

0.9 mm ( $\pm 0.1$  mm)

2.0 mm ( $\pm 0.2$  mm)

## Estimated resorption time

0.5 mm: about 1 month

0.9 mm: about 3 months

2.0 mm: about 5 months

## Size

0.5 mm: 20x20 mm

0.9 mm: 25x25 mm, 12x8 mm, 50x50 mm

2.0 mm: 7x5 mm, 15x5 mm, 30x30 mm, 50x50 mm

## Product codes

ED02LS | 20x20 mm | 0.5 mm | Porcine

ED21FS | 12x8 mm | 0.9 mm | Porcine

ED25FS | 25x25 mm | 0.9 mm | Porcine

ED05FS | 50x50 mm | 0.9 mm | Porcine

ED75SS | 7x5 mm | 2.0 mm | Porcine

ED15SS | 15x5 mm | 2.0 mm | Porcine

ED03SS | 30x30 mm | 2.0 mm | Porcine

ED05SS | 50x50 mm | 2.0 mm | Porcine

## GMDN code

47184

## CND code


P900402

## CHARACTERISTICS

Obtained from derma of porcine origin, using an exclusive Tecnos<sup>®</sup> process that preserves the natural collagen fibers<sup>(1)</sup>, *Derma* membranes are gradually integrated<sup>(2)</sup> with the autologous soft tissues. Their strong consistency and resistance allow a perfect stabilization and a prolonged protection of the underlying graft<sup>(3)</sup> in socket regeneration procedures, together with a strong barrier action to guide the growth of epithelium and preventing its invagination.


## HANDLING

*Derma* membrane can be shaped with scissors until the desired size is reached; then it must be thoroughly hydrated in sterile lukewarm physiological solution until the desired consistency is obtained. Once it acquires the desired plasticity, it must be adapted to the grafting site. It is always recommendable to prepare a pocket with an elevator in order to stabilize the membrane in the site after stitching the flaps.




SEM image of OsteoBiol<sup>®</sup> Derma

Source: Politecnico di Torino, Italy



4-months histology showing remnants of porcine dermal matrix surrounded by connective tissue. Experimental study on dogs. Htx-eosine

Author: Prof Ulf Nannmark, University of Göteborg, Sweden




Source: Tecnos<sup>®</sup> Dental Media Library

*Derma* membrane is a collagen resorbable matrix useful to augment soft tissues and to protect and stabilize bone grafting materials; only in this specific indication it can be used also in open healing<sup>(3)</sup> situations due to its perfect tissue integration characteristics.

If a residual band of keratinized tissue is still present around teeth or implants, *Derma* membrane can be used as an alternative to connective tissue graft<sup>(2)</sup> to improve the quality of keratinized tissues<sup>(4)</sup>. *Derma* has been also documented for horizontal soft tissue augmentation around implants<sup>(5,8)</sup>.

Mild gingival recessions<sup>(6,7)</sup> can be treated with *Derma* to avoid patient morbidity and discomfort due to connective tissue graft harvesting. To avoid membrane exposure, usually *Derma* is completely covered by the coronally advanced flap. A properly shaped *Derma* membrane with rounded edges has been also documented for the tunnel technique<sup>(6)</sup>.



Graft protection using OsteoBioL® *Derma*



Positioning of OsteoBioL® *Derma* with the tunneling technique



OsteoBioL® *Derma* shaped for a gingival recession treatment




Treatment of a gingival recession using OsteoBioL® *Derma*  
Author: Dr Roberto Rossi, Genova, Italy



**SOFT TISSUE AUGMENTATION**  
soft tissue improvement



**PERIODONTAL REGENERATION**  
gingival recessions



**ALVEOLAR REGENERATION**  
graft protection

## BIBLIOGRAPHY

(1) DE MARCO P, ZARA S, DE COLLI M, RADUNOVIC M, LAZOVIC V, ETTORRE V, DI CRESCENZO A, PIATELLI A, CATALDI A, FONTANA A  
**GRAPHENE OXIDE IMPROVES THE BIOCOMPATIBILITY OF COLLAGEN MEMBRANES IN AN IN VITRO MODEL OF HUMAN PRIMARY GINGIVAL FIBROBLASTS**  
BIOMED MATER, 2017 SEP 13;12(5):055005

(2) FICKL S, NANNMARK U, SCHLAGENHAUF U, HÜRZELER M, KEBSCHULL M  
**PORCINE DERMAL MATRIX IN THE TREATMENT OF DEHISCENCE-TYPE DEFECTS – AN EXPERIMENTAL SPLIT-MOUTH ANIMAL TRIAL**  
CLIN ORAL IMPLANTS RES, 2015 JUL;26(7):799-805

(3) TALLARICO M, XHANARI E, PISANO M, DE RIU G, TULLIO A, MELONI SM  
**SINGLE POST-EXTRACTIVE ULTRA-WIDE 7 MM-DIAMETER IMPLANTS VERSUS IMPLANTS PLACED IN MOLAR HEALED SITES AFTER SOCKET PRESERVATION FOR MOLAR REPLACEMENT: 6-MONTH POST-LOADING RESULTS FROM A RANDOMISED CONTROLLED TRIAL**  
EUR J ORAL IMPLANTOL, 2016;9(3):263-275

(4) FISCHER KR, FICKL S, MARDAS N, BOZEC L, DONOS N  
**STAGE-TWO SURGERY USING COLLAGEN SOFT TISSUE GRAFTS: CLINICAL CASES AND ULTRASTRUCTURAL ANALYSIS**  
QUINTESSENCE INT, 2014 NOV-DEC;45(10):853-60

(5) FISCHER K R, TESTORI T, WACHTEL H, MÜHLEMANN S, HAPPE A, DEL FABBRO M, DEL FABBRO M  
**SOFT TISSUE AUGMENTATION APPLYING A COLLAGENATED PORCINE DERMAL MATRIX DURING SECOND STAGE SURGERY: A PROSPECTIVE MULTICENTER CASE SERIES**  
CLIN IMPLANT DENT RELAT RES, 2019;1-8


(6) FICKL S, JOCKEL-SCHNEIDER Y, LINCKE T, BECHTOLD M, FISCHER KR, SCHLAGENHAUF U  
**PORCINE DERMAL MATRIX FOR COVERING OF RECESSION TYPE DEFECTS: A CASE SERIES**  
QUINTESSENCE INT, 2013;44(3):243-6

(7) MATOH U, PETELIN M, GASPERIC R  
**SPLIT-MOUTH COMPARISON OF CORONALLY ADVANCED FLAP WITH CONNECTIVE TISSUE GRAFT OR COLLAGEN MATRIX FOR TREATMENT OF ISOLATED GINGIVAL RECESSIONS**  
INT J PERIODONTICS RESTORATIVE DENT, 2019;39(3):439-446

(8) VERARDI S, ORSINI M, LOMBARDI T, AUSENDA F, TESTORI T, PULICI A, OREGLIA F, VALENTE NA, STACCHI C  
**COMPARISON BETWEEN TWO DIFFERENT TECHNIQUES FOR PERI-IMPLANT SOFT TISSUE AUGMENTATION: PORCINE DERMAL MATRIX GRAFT VS. TENTING SCREW**  
J PERIODONTOL. 2020; ACCEPTED, IN PUBLICATION

For further information see the complete literature on p. 92





# Lamina



**A unique cortical bone barrier**  
*Heterologous collagenated cortical bone*






# Characteristics and handling

## CHARACTERISTICS

*Lamina* barriers are made of cortical bone of heterologous origin produced with an exclusive TecnoSS® process that avoids the ceramization of hydroxyapatite crystals, thus allowing gradual resorption. After a process of superficial decalcification, *Lamina soft* acquires an elastic consistency, nevertheless maintaining the typical compactness of the bone tissue from which it originates; the margins are soft in order not to cause micro-traumas to the surrounding tissues. *Curved soft Lamina* has a semi-rigid consistency and should be grafted without hydration, provided that it is previously shaped to fit the defect morphology. *Rigid Lamina* undergoes a process of superficial semi-decalcification (50% vs *Lamina soft*) therefore increasing its consistency, typical of the cortical bone tissue<sup>(1-2)</sup>.



SEM image of OsteoBiol® Lamina  
Source: Politecnico di Torino, Italy

## HANDLING

*Lamina soft* can be shaped with sterile scissors until the desired size is reached, then it must be hydrated for 5/10 minutes in sterile physiological solution. Once it acquires the desired plasticity, it must be adapted to the grafting site; it should always be immobilized either with titanium microscrews or sutured (fine model) directly to the surrounding tissues with a triangular section non-traumatic needle. **Curved soft Lamina should not be hydrated** in order to maintain its tenting effect but can also be shaped with sterile scissors, and must be fixated with osteosynthesis screws. In case of exposure, *Lamina* should only be removed if there is a clear suprainfection, because its consistency allows to achieve a complete second intention healing of the wound.



Source: TecnoSS® Dental Media Library

### Tissue of origin

Cortical bone

### Tissue collagen

Preserved

### Physical form

Lamina soft: semi-rigid flexible dried lamina

Lamina: rigid dried lamina, flexible after re-hydration

### Composition

100% cortical bone

### Thickness

0.5 mm (±0.1 mm)

0.7 mm (±0.1 mm)

1.0 mm (±0.1 mm)

3.0 mm (±1 mm)

### Estimated re-entry time

0.5 mm: about 5 months

0.7 mm: about 7 months

1.0 mm: about 6 months

3.0 mm: about 8 months

### Size

0.5 mm: 25x25 mm, 25x35 mm (oval)

0.7 mm: 35x15 mm

1.0 mm: 35x35 mm (Curved), 20x40 mm

3.0 mm: 30x30 mm

### Product codes

LS25FS | 25x25 mm | 0.5 mm | soft | Porcine

LS25FE | 25x25 mm | 0.5 mm | soft | Equine

LS23FS | 25x35 mm (Oval) | 0.5 mm | soft | Porcine

LS23FE | 25x35 mm (Oval) | 0.5 mm | soft | Equine

LS24LS | 20x40 mm | 1.0 mm | soft | Porcine

LS10HS | 35x35 mm (Curved) | 1.0 mm | soft | Porcine

LS10HE | 35x35 mm (Curved) | 1.0 mm | soft | Equine

LS03SS | 30x30 mm | 3.0 mm | soft | Porcine

LS03SE | 30x30 mm | 3.0 mm | soft | Equine

LS15LS | 35x15 mm | 0.7 mm | rigid | Porcine

LS35LS | 35x35 mm | 1.0 mm | rigid | Porcine

### GMDN code

46425

### CND code

P900402



*Lamina soft* becomes flexible after hydration and can be shaped<sup>(3)</sup> and adapted to the defect morphology creating, once fixated with osteosynthesis screws, a semi-rigid covering to the underlying graft<sup>(4-6)</sup>. This property is particularly useful when it is necessary to maintain the graft volume in aesthetic areas, as well as in horizontal augmentation<sup>(6-8)</sup> of two wall defects and in lateral access sinus lift procedures<sup>(5,9,10)</sup>. *Lamina* can also be used in regenerations with risks of exposure. *Curved soft Lamina* has a ~1.0 mm thickness and must be directly grafted without hydration<sup>(11)</sup>; it can be particularly effective in association with GTO® for regeneration of ridges with compromised buccal plate.


*Rigid Lamina* (ref. LS35LS) has been documented for orbital floor and wall reconstruction<sup>(12)</sup>.

The new 0.7 mm thickness rigid *Lamina* (ref. LS15LS) represents a viable alternative to autogenous cortical bone plates in the reconstruction of three-dimensional crestal defects with the shell technique.

The above clinical information is based on the experience of expert surgeons




OsteoBiol® Lamina positioning  
Source: Tecness® Dental Media Library



**HORIZONTAL AUGMENTATION**  
two-wall defects



**HORIZONTAL AUGMENTATION**  
bone-layer technique




**ORBITAL FLOOR RESTORATION**

## BIBLIOGRAPHY


- (1) DI CARLO R, ZARA S, VENTRELLA A, SIANI G, DA ROS T, IEZZI G, CATALDI A, FONTANA A  
**COVALENT DECORATION OF CORTICAL MEMBRANES WITH GRAPHENE OXIDE AS A SUBSTRATE FOR DENTAL PULP STEM CELLS**  
NANOMATERIALS, 2019;9:604
- (2) CABALLÉ-SERRANO J, MUNAR-FRAUJA A, DELGADO L, PÉREZ R, HERNÁNDEZ-ALFARO F  
**PHYSICO-CHEMICAL CHARACTERIZATION OF BARRIER MEMBRANES FOR BONE REGENERATION**  
J MECH BEHAV BIOMED, 2019;97:13-20
- (3) ROSSI R, RANCITELLI D, POLI PP, RASIA DAL POLO M, NANNMARK U, MAIORANA C  
**THE USE OF A COLLAGENATED PORCINE CORTICAL LAMINA IN THE RECONSTRUCTION OF ALVEOLAR RIDGE DEFECTS. A CLINICAL AND HISTOLOGICAL STUDY**  
MINERVA STOMATOL, 2016 OCT;65(5):257-68
- (4) PAGLIANI L, ANDERSSON P, LANZA M, NAPPO A, VERROCCHI D, VOLPE S, SENNERBY L  
**A COLLAGENATED PORCINE BONE SUBSTITUTE FOR AUGMENTATION AT NEOSSE IMPLANT SITES: A PROSPECTIVE 1-YEAR MULTICENTER CASE SERIES STUDY WITH HISTOLOGY**  
CLIN IMPLANT DENT RELAT RES, 2012 OCT;14(5):746-58
- (5) FESTA VM, ADDABBO F, LAINO L, FEMIANO F, RULLO R  
**PORCINE-DERIVED XENOGRAFT COMBINED WITH A SOFT CORTICAL MEMBRANE VERSUS EXTRACTION ALONE FOR IMPLANT SITE DEVELOPMENT: A CLINICAL STUDY IN HUMANS**  
CLIN IMPLANT DENT AND RELAT RES, 2013 OCT;15(5):707-13
- (6) WACHTEL H, FICKL S, HINZE M, BOLZ W, THALMAIR T  
**THE BONE LAMINA TECHNIQUE: A NOVEL APPROACH FOR LATERAL RIDGE AUGMENTATION - A CASE SERIES**  
INT J PERIODONTICS RESTORATIVE DENT, 2013 JUL-AUG;33(4):491-7
- (7) LOPEZ MA, ANDREASI BASSI M, CONFALONE L, CARINCI F, ORMIANER Z, LAURITANO D  
**THE USE OF RESORBABLE CORTICAL LAMINA AND MICRONIZED COLLAGENATED BONE IN THE REGENERATION OF ATROPHIC CRESTAL RIDGES: A SURGICAL TECHNIQUE. CASE SERIES**  
J BIOL REGUL HOMEOST AGENTS, 2016 APR-JUN;30(2 SUPPL 1):81-85
- (8) POLIS-YANES C, CADENAS-SEBASTIÁN C, GUAL-VAQUÉS P, AYUSO-MONTERO R, MARÍ-ROIG A, LÓPEZ-LÓPEZ J  
**GUIDED BONE REGENERATION OF AN ATROPHIC MAXILLA USING HETEROLOGOUS CORTICAL LAMINA**  
CASE REP DENT, 2019; 5216362
- (9) HINZE M, VRIELINCK L, THALMAIR T, WACHTEL H, BOLZ W  
**ZYGOMATIC IMPLANT PLACEMENT IN CONJUNCTION WITH SINUS BONE GRAFTING: THE "EXTENDED SINUS ELEVATION TECHNIQUE". A CASE-COHORT STUDY**  
ORAL CRANIOFAC TISSUE ENG, 2011;1:188-197
- (10) SCARANO A, MURMURA G, MASTRANGELO F, LORUSSO F, GRECO LUCCHINA A, CARINCI F  
**A NOVEL TECHNIQUE TO PREVENT SINUS MEMBRANE COLLAPSE DURING MAXILLARY SINUS FLOOR AUGMENTATION WITHOUT BONE GRAFT: TECHNICAL NOTE**  
J BIOL REGUL HOMEOST AGENTS, 2018 NOV-DEC;32(6):1589-1592
- (11) ROSSI R, FOCE E, SCOLAVINO S  
**THE CORTICAL LAMINA TECHNIQUE: A NEW OPTION FOR ALVEOLAR RIDGE AUGMENTATION. PROCEDURE, PROTOCOL, AND CASE REPORT**  
J LEBANESE DENTAL ASS, 2017 JAN-JUN; 52(1):35-41
- (12) RINNA C, REALE G, FORESTA E, MUSTAZZA MC  
**MEDIAL ORBITAL WALL RECONSTRUCTION WITH SWINE BONE CORTEX**  
J CRANIOFAC SURG, 2009 MAY; 20(3):881-4

For further information see the complete literature on p. 92





# Special




***A translucent membrane to  
separate bone and soft tissues***

***Engineered to protect hard and soft tissue grafts***



# Characteristics, handling and clinical information



## Tissue of origin

Heterologous pericardium

## Tissue collagen

Preserved

## Physical form

Translucent dried membrane

## Composition

100% pericardium

## Thickness

0.2 mm

## Resorption time

About 40 days

## Size

20x20 mm, 30x30 mm

## Product codes

EM02LE | 20x20 mm | Equine

EM03LE | 30x30 mm | Equine

## GMDN code

47184

## CND code

P900402

## CHARACTERISTICS

Obtained from pericardium of heterologous origin, using an exclusive TecnoSS® process, the dried *Special* membranes are completely resorbable. Once hydrated, they become translucent and flexible, guiding the growth of epithelium and preventing its invagination: their action favors therefore an optimal regeneration of the underlying bone tissue.

## HANDLING


The membrane can be shaped with sterile scissors until the desired size is reached; it must then be rehydrated with lukewarm physiological solution. Once it acquires the desired plasticity, it must be adapted to the grafting site. It is recommended to prepare a pocket with an elevator in order to stabilize the membrane in the site after stitching the flaps. If this is not possible, the membrane can be stabilized with envelope sutures which bridle it with the gingival flaps.

## CLINICAL INFORMATION

In periodontology, the *Special* membrane can be used to protect and stabilize the biomaterial in the treatment of intrabony defects.

*Special* can be used to protect the sinus membrane before the insertion of the grafting material, to close sinus membrane perforations. Grafts placed in post-extractive sockets with closed healing procedure can also be protected with this membrane.

*The above clinical information is based on the experience of expert surgeons*



SEM images of OsteoBio® Special  
Source: Nobil Bio Ricerche, Villafranca d'Asti, Italy



**PERIODONTAL REGENERATION**  
intrabony defects



**LATERAL ACCESS SINUS LIFT**  
Schneider membrane protection

# Bone, Biomaterials & Beyond

Prof Antonio Barone, Prof Ulf Nannmark

## CONTENTS

The introduction of osseointegrated dental implants soon 50 years ago has indeed revolutionized dentistry.


The scientific evaluation of their use has shown good and increasingly successful treatment outcomes.

A prerequisite though is the availability of sufficient bone volumes to ensure integration and acceptable aesthetic results.

In this book, various surgical techniques using different augmentation materials are described and explained.

The aim is to highlight minimally invasive surgical techniques, which lead to less risk of morbidity and reduce treatment time.

Readers will enjoy a comprehensive atlas providing some practical advice for every day surgical practice based on solid scientific evidence.



## EDRA Editions

Publication: March 2014

Pages: 200

Images: 786 color

Binding: hardcover

Size: 21x29,7 cm

Euro: 59,00

ISBN 978.88.214.3758.8

For information and orders

please contact: [ordiniedra@lswr.it](mailto:ordiniedra@lswr.it)

Available in the following languages:

English

Italian

German

French

Spanish

Russian

Korean



# Success through innovation:

Products launch



Gen-Os®



Apatos



Gel 40



Special



Lamina



Sp-Block



Dual-Block



Derma

Year

2000

2001

2002

2003

2004

2005

2006

2007

2008

2009

2010

Worldwide distribution countries

1

2

8

12

19

28

41

44

Publications on international journals

1

4

9

12

19

26

32

# history of the **OsteoBiol<sup>®</sup>** brand by Tecoss



**Putty**



**Evolution**



**mp3<sup>®</sup>**



**TSV Gel**



**GTO<sup>®</sup>**



**Lamina**  
for shell technique





**INNOVATION**

A close-up photograph of a scientist wearing a white lab coat, a white surgical cap, and a white face mask. The scientist is looking through the eyepieces of a white and black microscope. The background is a blurred laboratory setting with various pieces of equipment. The word "INNOVATION" is overlaid in large, bold, blue capital letters on the left side of the image.

# Tecnoss® bone vs human bone

Studies and researches have demonstrated that gold standard in bone regeneration is autologous bone<sup>(1,2)</sup>.

It is also well known, though, what disadvantages are related to the harvesting and grafting of autogenous bone<sup>(2-4)</sup>.

The goal of bone regeneration is to heal bone deficits with newly-formed quality tissue, in order to achieve a functional recovery and esthetics. To obtain these results, hundreds of studies have been conducted about the clinical performance of biomaterials. The examination of clinical results and the commercial diffusion of various kinds of products developed by the biomedical industry show a

clear superiority of products of natural origin over those of synthetic derivation.

The structure of animal bone is morphologically more similar to human bone than any synthesized product, the latter presenting a morphological pattern and properties artificially created, which differ in various ways from the structure of natural bone<sup>(5)</sup>.

Over the last thirty years several processes have been developed to allow the grafting of heterologous

origin products in the human body without adverse reaction<sup>(6,7)</sup>.

The first products developed through these technologies have shown encouraging clinical results, even if made of bone mineral matrix only.

The OsteoBiol® new generation of biomaterials, thanks to a revolutionary technology, goes beyond the simple role of aiding natural bone regrowth by stimulating and accelerating contact osteogenesis, with a behaviour similar to that of autogenous bone<sup>(8-10)</sup>.



## BIBLIOGRAPHY

- (1) ORSINI G, SCARANO A, PIATTELLI M, PICCIRILLI M, CAPUTI S, PIATTELLI A  
**HISTOLOGIC AND ULTRASTRUCTURAL ANALYSIS OF REGENERATED BONE IN MAXILLARY SINUS AUGMENTATION USING A PORCINE BONE-DERIVED BIOMATERIAL**  
J PERIODONTOL, 2006 DEC;77(12):1984-90
- (2) BARONE A, ALFONSI F, BORGIA V, IEZZI G, PIATTELLI A, COVANI U, TONELLI P  
**MOLECULAR, CELLULAR AND PHARMACEUTICAL ASPECTS OF FILLING BIOMATERIALS DURING THE MANAGEMENT OF EXTRACTION SOCKETS**  
CURR PHARM BIOTECHNOL, 2017;18(1):64-75
- (3) IEZZI G, PIATTELLI A, GIULIANI A, MANGANO C, BARONE A, MANZONI L, DEGIDI M, SCARANO A, FILIPPONE A, PERROTTI V  
**MOLECULAR, CELLULAR AND PHARMACEUTICAL ASPECTS OF FILLING BIOMATERIALS DURING MAXILLARY SINUS-LIFT PROCEDURES. PART 2: DETAILED CHARACTERISTICS OF THE MATERIALS**  
CURR PHARM BIOTECHNOL, 2017, 18, 33-44
- (4) BARONE A, CRESPI R, ALDINI NN, FINI M, GIARDINO R, COVANI U  
**MAXILLARY SINUS AUGMENTATION: HISTOLOGIC AND HISTOMORPHOMETRIC ANALYSIS**  
INT J ORAL MAXILLOFAC IMPLANTS, 2005 JUL-AUG; 20(4):519-25
- (5) FIGUEIREDO A, COIMBRA P, CABRITA A, GUERRA F, FIGUEIREDO M  
**COMPARISON OF A XENOGENIC AND AN ALLOPLASTIC MATERIAL USED IN DENTAL IMPLANTS IN TERMS OF PHYSICO-CHEMICAL CHARACTERISTICS AND IN VIVO INFLAMMATORY RESPONSE**  
MATER SCI ENG C MATER BIOL APPL, 2013 AUG 1;33(6):3506-13
- (6) SCARANO A, PIATTELLI A, ASSENZA B, QUARANTA A, PERROTTI V, PIATTELLI M, IEZZI G  
**PORCINE BONE USED IN SINUS AUGMENTATION PROCEDURES: A 5-YEAR RETROSPECTIVE CLINICAL EVALUATION**  
J ORAL MAXILLOFAC SURG, 2010 AUG; 68(8):1869-73
- (7) RAMIREZ FERNANDEZ MP, CALVO GUIRADO JL, MATÉ SANCHEZ DE VAL JE, DELGADO RUIZ RA, NEGRI B, BARONA DORADO C  
**ULTRASTRUCTURAL STUDY BY BACKSCATTERED ELECTRON IMAGING AND ELEMENTAL MICROANALYSIS OF BONE-TO-BIOMATERIAL INTERFACE AND MINERAL DEGRADATION OF PORCINE XENOGRAPTS USED IN MAXILLARY SINUS FLOOR ELEVATION**  
CLIN ORAL IMPLANTS RES, 2013 MAY;24(5):523-30
- (8) CASSETTA M, RICCI L, IEZZI G, DELLAQUILA D, PIATTELLI A, PERROTTI V  
**RESONANCE FREQUENCY ANALYSIS OF IMPLANTS INSERTED WITH A SIMULTANEOUS GRAFTING PROCEDURE: A 5-YEAR FOLLOW-UP STUDY IN MAN**  
INT J PERIODONTICS RESTORATIVE DENT, 2012 OCT;32(5):581-9
- (9) CASSETTA M, PERROTTI V, CALASSO S, PIATTELLI A, SINJARI B, IEZZI G  
**BONE FORMATION IN SINUS AUGMENTATION PROCEDURES USING AUTOLOGOUS BONE, PORCINE BONE, AND A 50 : 50 MIXTURE: A HUMAN CLINICAL AND HISTOLOGICAL EVALUATION AT 2 MONTHS**  
CLIN ORAL IMPLANTS RES, 2015 OCT;26(10):1180-4
- (10) JEANNEAU C, LE FOURNIS C, ABOUT I  
**XENOGENIC BONE FILLING MATERIALS MODULATE MESENCHYMAL STEM CELL RECRUITMENT: ROLE OF THE COMPLEMENT C5A**  
CLINICAL ORAL INVESTIGATIONS, 2019 OCT 23



## Why xenografts?



Xenografts are the most used biomaterials worldwide.

This is because:

- tissues of origin are extremely safe and available in unlimited quantities
- xenogenic bone surface and porosity are extremely similar to autogenous bone
- there is no need to harvest autogenous bone in extraoral sites, with the related risk of morbidity and post-operative complications
- sterile xenografts are completely biocompatible and safe
- no adverse reactions after grafting deriving from biomaterial degradation
- easy to handle, quick learning curve
- collagenated xenografts enhance osteoblasts and osteoclasts activity
- wide scientific documentation
- excellent clinical performance
- storage can be done at room temperature
- long shelf life (5 years from production date)
- excellent price/quality ratio

# Characteristics of Tecnos® process

Tecnos® has developed manufacturing processes for the treatment of tissues from various animal species, allowing to obtain the biocompatibility of these tissues, preserving at the same time their collagen matrix<sup>(1)</sup>.

The protein components of animal tissues are determinant to make every individual unique. They activate the cells of the immune system of the receiving organism by interacting with receptors of the Major Histocompatibility Complex (MHC).

Their neutralization/denaturation allows collagen mineral bone matrix to be transferred from animal to man without any dangerous adverse reaction outbreak.

Successful Guided Bone Regeneration (GBR) depends both on stimulation of tissues involved in new bone formation and on the characteristics of grafted biomaterials, which can determine the quality of bone/graft interface. The development of OsteoBio® product line has thus been driven by the ideal

biomaterial concept: a material with the highest affinity to the new endogenous bone.

To pursue this aim, Tecnos® developed a biotechnology able to preserve the structure of natural hydroxyapatite, avoiding the high temperature ceramization phase, therefore allowing a bone turnover time of the grafted site similar to the one of the physiologic natural process<sup>(2)</sup>.

Thanks to this innovative technology, the OsteoBio® line has the following important characteristics:

1. Cell growth support and differentiation<sup>(3)</sup>
2. Absence of a foreign body response<sup>(4,5)</sup>
3. Gradual resorption over time<sup>(2,6)</sup>
4. Stimulation of the physiological tissue regeneration process<sup>(7,8)</sup>
5. Protection of the grafting site from infection (membranes)<sup>(5,9)</sup>
6. Capability of carrying medication to the surgical site<sup>(10)</sup>
7. Absorption and release over time of growth factors<sup>(11)</sup>
8. Enhancement of endothelial cells proliferation<sup>(7)</sup>



## BIBLIOGRAPHY

(1) FIGUEIREDO M, HENRIQUES J, MARTINS G, GUERRA F, JUDAS F, FIGUEIREDO H  
**PHYSICO-CHEMICAL CHARACTERIZATION OF BIOMATERIALS COMMONLY USED IN DENTISTRY AS BONE SUBSTITUTES - COMPARISON WITH HUMAN BONE**  
J BIOMED MATER RES B APPL BIOMATER, 2010 FEB; 92(2):409-19

(2) NANNMARK U, SENNERBY L  
**THE BONE TISSUE RESPONSES TO PREHYDRATED AND COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS: A STUDY IN RABBIT MAXILLARY DEFECTS**  
CLIN IMPLANT DENT RELAT RES, 2008 DEC;10(4):264-70

(3) TRUBIANI O, SCARANO A, ORSINI G, DI IORIO D, D'ARCANGELO C, PICCIRILLI M, SIGISMONDO M, CAPUTI S  
**THE PERFORMANCE OF HUMAN PERIODONTAL LIGAMENT MESENCHYMAL STEM CELLS ON XENOGENIC BIOMATERIALS**  
INT J IMMUNOPATHOL PHARMACOL, 2007 JAN-MAR; 20(1 SUPPL 1):87-91

(4) BARONE A, RICCI M, GRASSI RF, NANNMARK U, QUARANTA A, COVANI U  
**A 6-MONTH HISTOLOGICAL ANALYSIS ON MAXILLARY SINUS AUGMENTATION WITH AND WITHOUT USE OF COLLAGEN MEMBRANES OVER THE OSTEOTOMY WINDOW: RANDOMIZED CLINICAL TRIAL**  
CLIN ORAL IMPLANTS RES, 2013 JAN;24(1):1-6

(5) BARONE A, BORGIA V, COVANI U, RICCI M, PIATTELLI A, IEZZI G  
**FLAP VERSUS FLAPLESS PROCEDURE FOR RIDGE PRESERVATION IN ALVEOLAR EXTRACTION SOCKETS: A HISTOLOGICAL EVALUATION IN A RANDOMIZED CLINICAL TRIAL**  
CLIN ORAL IMPLANTS RES, 2015 JUL;26(7):806-13

(6) BARONE A, RICCI M, COVANI U, NANNMARK U, AZARMEHR I, CALVO GUIRADO JL  
**MAXILLARY SINUS AUGMENTATION USING PREHYDRATED CORTICOCANCELLOUS PORCINE BONE: HISTOMORPHOMETRIC EVALUATION AFTER 6 MONTHS**  
CLIN IMPLANT DENT RELAT RES, 2012 JUN;14(3):373-9

(7) ROMBOUITS C, JEANNEAU C, CAMILLERI J, LAURENT P, ABOUT I  
**CHARACTERIZATION AND ANGIOGENIC POTENTIAL OF XENOGENIC BONE GRAFTING MATERIALS: ROLE OF PERIODONTAL LIGAMENT CELLS**  
DENT MATER J, 2016 DEC 1;35(6):900-907

(8) JEANNEAU C, LE FOURNIS C, ABOUT I  
**XENOGENIC BONE FILLING MATERIALS MODULATE MESENCHYMAL STEM CELL RECRUITMENT: ROLE OF THE COMPLEMENT C5A**  
CLINICAL ORAL INVESTIGATIONS, 2019 OCT 23

(9) BARONE A, RICCI M, TONELLI P, SANTINI S, COVANI U  
**TISSUE CHANGES OF EXTRACTION SOCKETS IN HUMANS: A COMPARISON OF SPONTANEOUS HEALING VS. RIDGE PRESERVATION WITH SECONDARY SOFT TISSUE HEALING**  
CLIN ORAL IMPLANTS RES, 2013 NOV;24(11):1231-7

(10) FISCHER KR, STAVROPOULOS A, CALVO GUIRADO JL, SCHNEIDER D, FICKL S  
**INFLUENCE OF LOCAL ADMINISTRATION OF PAMIDRONATE ON EXTRACTION SOCKET HEALING – A HISTOMORPHOMETRIC PROOF-OF-PRINCIPLE PRE-CLINICAL IN VIVO EVALUATION**  
CLIN ORAL IMPLANTS RES, 2015 OCT;26(10):1135-42

(11) MJIRITSKY E, FERRONI L, GARDIN C, BRESSAN E, ZANETTE G, PIATTELLI A, ZAVAN B  
**PORCINE BONE SCAFFOLDS ADSORB GROWTH FACTORS SECRETED BY MSCS AND IMPROVE BONE TISSUE REPAIR MATERIALS**, 2017 SEP 8;10(9)



# Collagen: a key factor for clinical success

## BIBLIOGRAPHY

(1) FIGUEIREDO M, HENRIQUES J, MARTINS G, GUERRA F, JUDAS F, FIGUEIREDO H

**PHYSICO-CHEMICAL CHARACTERIZATION OF BIOMATERIALS COMMONLY USED IN DENTISTRY AS BONE SUBSTITUTES - COMPARISON WITH HUMAN BONE**

J BIOMED MATER RES B APPL BIOMATER, 2010 FEB; 92(2):409-19

(2) ORSINI G, SCARANO A, PIATTELLI M, PICCIRILLI M, CAPUTI S, PIATTELLI A

**HISTOLOGIC AND ULTRASTRUCTURAL ANALYSIS OF REGENERATED BONE IN MAXILLARY SINUS AUGMENTATION USING A PORCINE BONE-DERIVED BIOMATERIAL**

J PERIODONTOL, 2006 DEC;77(12):1984-90

(3) RAMIREZ FERNANDEZ MP, CALVO GUIRADO JL, MATÉ SANCHEZ DE VAL JE, DELGADO RUIZ RA, NEGRI B, BARONA DORADO C

**ULTRASTRUCTURAL STUDY BY BACKSCATTERED ELECTRON IMAGING AND ELEMENTAL MICROANALYSIS OF BONE-TO-BIOMATERIAL INTERFACE AND MINERAL DEGRADATION OF PORCINE XENOGRAFTS USED IN MAXILLARY SINUS FLOOR ELEVATION**

CLIN ORAL IMPLANTS RES, 2013 MAY;24(5):523-30

(4) FELICE P, PIANA L, CHECCHI L, CORVINO V, NANNMARK U, PIATTELLI M

**VERTICAL RIDGE AUGMENTATION OF ATROPHIC POSTERIOR MANDIBLE WITH AN INLAY TECHNIQUE AND CANCELLOUS EQUINE BONE BLOCK: A CASE REPORT**

INT J PERIODONTICS RESTORATIVE DENT, 2013 MAR-APR;33(2):159-66

(5) IEZZI G, PIATTELLI A, GIULIANI A, MANGANO C, BARONE A, MANZON L, DEGIDI M, SCARANO A, FILIPPONE A, PERROTTI V

**MOLECULAR, CELLULAR AND PHARMACEUTICAL ASPECTS OF FILLING BIOMATERIALS DURING MAXILLARY SINUS-LIFT PROCEDURES. PART 2: DETAILED CHARACTERISTICS OF THE MATERIALS**

CURR PHARM BIOTECHNOL, 2017,18,33-44

(6) MIZUNO M, FUJISAWA R, KUBOKI Y

**TYPE I COLLAGEN-INDUCED OSTEOBLASTIC DIFFERENTIATION OF BONE-MARROW CELLS MEDIATED BY COLLAGEN-A2B1 INTEGRIN INTERACTION**

J CELL PHYSIOL. 2000 AUG;184(2):207-13

(7) HSU FY, CHUEH SC, WANG YJ

**MICROSPHERES OF HYDROXYAPATITE/RECONSTITUTED COLLAGEN AS SUPPORTS FOR OSTEOBLAST CELL GROWTH**

BIOMATERIALS 1999, 20:1931-1936

(8) ABDELGAWAD ME, SØE K, ANDERSEN TL, MERRILL DM, CHRISTIANSEN P, KJÆRGAARD-ANDERSEN P, DELAISSE JM

**DOES COLLAGEN TRIGGER THE RECRUITMENT OF OSTEOBLASTS INTO VACATED BONE RESORPTION LACUNAE DURING BONE REMODELING?**

BONE, 2014 OCT;67:181-8

(9) ROMBOUITS C, JEANNEAU C, CAMILLERI J, LAURENT P, ABOUT I

**CHARACTERIZATION AND ANGIOGENIC POTENTIAL OF XENOGENEIC BONE GRAFTING MATERIALS: ROLE OF PERIODONTAL LIGAMENT CELLS**

DENT MATER J, 2016 DEC 1;35(6):900-907

(10) JEANNEAU C, LE FOURNIS C, ABOUT I

**XENOGENEIC BONE FILLING MATERIALS MODULATE MESENCHYMAL STEM CELL RECRUITMENT: ROLE OF THE COMPLEMENT C5A**

CLINICAL ORAL INVESTIGATIONS, 2019 OCT 23

Tecnoss<sup>®</sup> exclusive manufacturing process is able to neutralize the antigenic components present in heterologous bone achieving biocompatibility and preserving the collagen matrix inside the biomaterial. Moreover, the molecular structure of natural hydroxyapatite is not significantly altered thanks to the mild process temperature<sup>(1)</sup>.

These characteristics of OsteoBiol<sup>®</sup> products allow a consistent bone neo-formation and a close contact between mature neo-formed bone and biomaterial particles<sup>(2-5)</sup>.

Collagen has a key role in bone regeneration process in that:

- it acts as a valid substrate for platelet activation and aggregation
- it serves to attract and differentiate the mesenchymal stem cells present in the bone marrow<sup>(6)</sup>
- it increases the proliferation rate of the

osteoblasts up to 2/3 times<sup>(7)</sup>


- it stimulates the activation of the platelets, osteoblasts and osteoclasts in the bone healing process<sup>(8)</sup>.

The presence of collagen inside each granule makes OsteoBiol<sup>®</sup> Gen-Os<sup>®</sup> hydrophilic and facilitates further mixing with collagen gel and TSV Gel.


This technology has permitted the development of several versatile and innovative products: OsteoBiol<sup>®</sup> GTO<sup>®</sup>, OsteoBiol<sup>®</sup> mp3<sup>®</sup>, OsteoBiol<sup>®</sup> Putty and OsteoBiol<sup>®</sup> Gel 40. Their consistency allows an ideal filling of bone defects and guarantees simple handling and fast application.

The OsteoBiol<sup>®</sup> new generation of biomaterials, thanks to a revolutionary technology, goes beyond the simple role of aiding natural bone regrowth by stimulating and accelerating this vital physiological process<sup>(9,10)</sup>.

Composition of **OsteoBiol<sup>®</sup> Gen-Os<sup>®</sup>**



Source: University of Duisburg-Essen, Germany



Guided bone regeneration (GBR) is necessary to treat bone deficits due to lesions or bacterial infections.

Bone defect recovery occurs through the general mechanisms of tissue healing: complex dynamic mechanisms directed towards the repair of tissue function and anatomic integrity.

The discovery of the events pathway leading to tissue healing has helped to clearly identify the main actors in bone healing process; the concomitant presence of the following three components is necessary for the formation of “*de novo*” bone tissue:

- the platelets represent the principal actors during the first phase of the healing process, when, subsequent to a lesion, an initial deposition of fibrin and the formation of blood clot take place. This phase is characterized by significant activation of the chemical signals mediated by cytokines and growth factors.

In fact, the primary post-haemorrhagic clot formation process through platelet aggregation and lysis causes the release of both the coagulation cascade factors and growth factors, such as PDGF, IGF 1, IGF 2 and VEGF which are known for their activating effect on osteoblasts and osteoclasts, and TGF- $\beta$  (Bone Morphogenetic Proteins belong to this superfamily) which starts bony callus formation.

- the osteoblastic precursors deriving from bone marrow mesenchymal stem cells are responsible, after cell differentiation in osteoblasts, for the second phase of the healing process (enchondral and/or intramembranous ossification) thanks to the synthesis of collagen and other components of the

extracellular matrix.

- an insoluble substrate, suitable carrier for osteoinductive signal and able to support and guide new bone tissue formation. Sampath and Reddi (1980) demonstrated crosslinked type I collagen to be the most appropriate carrier for promoting osteoinductive signal activity. The continuous progresses in comprehension of biological mechanisms regulating bone tissue morphogenesis can be exploited also for elaboration of natural or artificial products able to restore or maintain the function of damaged tissues and organs (tissue engineering)<sup>(1-3)</sup>.

In vitro studies demonstrated that heterologous collagen is able to induce differentiation of mesenchymal osteoprogenitor stem cells into osteoblasts<sup>(4)</sup>, and that association of collagen type I with a scaffold of hydroxyapatite significantly enhances osteoblasts proliferation rate.

This important scientific evidence provides the rationale behind OsteoBioI<sup>®</sup> product line: a complete series of biomaterials with collagen base.

Collagen, in addition to its well-known structural action carried on connective tissues, is endowed with the following important properties, useful in tissue reparation processes:

## 1. Haemostasis

Collagen is able to activate the receptors present on cellular membranes of platelets, responsible for their aggregation and lysis process; moreover, during the first week, it reinforces the action of fibrin in the formation of the primary clot, and then, in the second week, it replaces the

function of fibrin.

## 2. Debridement

Collagen has a chemotactic action on monocyte/macrophage cell lines, from which osteoclasts derive; these cells, through their action on mineral component resorption of both bone tissue and OsteoBioI<sup>®</sup> biomaterials, can draw, activate and collaborate with osteoblasts in bone rearranging and remodeling.

## 3. Angiogenesis

The drawn monocytes/macrophages, in their turn, stimulate both osteoblastic activity and angiogenesis process in grafting site.

## 4. Osteoblastic activity

Collagen, binding to fibronectin, promotes the anchorage of mesenchymal stem progenitors, on which it exerts its chemotactic action, and induces differentiation into osteoblasts<sup>(4,5)</sup>.

## 5. Receiving site remodeling

Exogenous collagen grafting can contribute in decreasing remodeling times of immature bone tissue.

## 6. Osteoconduction and guided regeneration

Naturally integrated with mineral component, collagen is able to increase osteoblasts proliferation rate while as a resorbable membrane it is able to guide connective tissue regeneration.

## BIBLIOGRAPHY


(1) GRIFFITH LG, NAUGHTON G  
**TISSUE ENGINEERING-CURRENT CHALLENGES AND EXPANDING OPPORTUNITIES**  
SCIENCE 2002, 295:1009-14

(2) REDDI AH  
**MORPHOGENESIS AND TISSUE ENGINEERING OF BONE AND CARTILAGE: INDUCTIVE SIGNALS, STEM CELLS, AND BIOMIMETIC BIOMATERIALS**  
TISSUE ENG 2000, 6(4):351-59

(3) NAKASHIMA N, REDDI AH  
**THE APPLICATION OF BONE MORPHOGENETIC PROTEINS TO DENTAL TISSUE ENGINEERING**  
NAT BIOTECHNOL 2003, 9:1025-32

(4) SALASZNYK RM, WILLIAMS WA, BOSKEY A, BATORSKY A, PLOPPER GE  
**ADHESION TO VITRONECTIN AND COLLAGEN I PROMOTES OSTEOGENIC DIFFERENTIATION OF HUMAN MESENCHYMAL STEM CELLS**  
J BIOMED BIOTECHNOL 2004, 1:24-34

(5) BRUNELLI G, SOLLAZZO V, CARINCI F, PALMIERI A, GIRARDI A, MONGUZZI R  
**OSTEOBIOL<sup>®</sup> INFLUENCES OSTEOGENIC DIFFERENTIATION OF ADIPOSE DERIVED STEM CELLS**  
EUR J INFLAMMAT, 2011, VOL. 9, NO. 3(S), 103-107





# From heterologous bone to biomaterial

## RESULTS OF **INORGANIC** CHEMICAL ANALYSES PERFORMED ON OSTEObIOL® GEN-OS®

Chemical element	OsteoBiol® Gen-OS® (% in weight)	
Ca	25.7%	
PO <sub>4</sub> <sup>3-</sup>	35.2%	
C	13.6%	
H	2.2%	
N	2.9%	
O (not in PO <sub>4</sub> <sup>3-</sup> )	20.4%	
<b>TOTAL</b>	<b>100.0%</b>	<b>Mineral component</b> 73.6%
Ca/P (n:n)	1.73	<b>Organic matrix</b> 22.4%
		<b>Water</b> 4.0%

### Inorganic chemical analyses results

Source: University of Duisburg-Essen, Germany

## RESULTS OF **ORGANIC** CHEMICAL ANALYSES PERFORMED ON OSTEObIOL® GEN-OS®



“The separated proteins (one lane) were fractionated in ten portions and analysed with nano-LC-ESI MS/MS. In the fractions 1-5 in the range from 20-200kDa we found ONLY COLLAGEN. In the fractions 6-10 we identify NO PROTEIN”

### Organic chemical analyses results

Source: Proteome Factory, Germany

A biomaterial for the reconstruction of bone defects must be biocompatible and have good handling and modeling properties; in specific clinical situations, it must also provide sufficient mechanical resistance. TecnoSS® laboratories are specialized in processing heterologous bony and collagenic tissues. OsteoBiol® bone process, in particular, has been developed to modify while preserving the original collagen matrix of heterologous tissue, in order to maintain its positive biological functions, and complete biocompatibility<sup>(1)</sup>. Most biomaterials are inert products that do not interfere, or rather, do not take

part in the physiology of bone remodeling: since they have been developed according to the sole concept of biocompatibility, their function is limited to the preservation of the graft volume (scaffold). The biocompatibility concept has an essential purpose in the implant of permanent prosthetic elements inside the human body, but it is extremely restrictive in case of materials used for bone reconstruction. OsteoBiol® biomaterials, being gradually resorbed and replaced by abundant newly formed bone over time, create the ideal conditions for the osseointegration of dental implants at re-entry<sup>(2)</sup>.

(1) NANNMARK U, SENNERBY L  
THE BONE TISSUE RESPONSES TO PREHYDRATED AND COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS: A STUDY IN RABBIT MAXILLARY DEFECTS  
CLIN IMPLANT DENT RELAT RES, 2008 DEC;10(4):264-70.

(2) CHECCHI V, FELICE P, ZUCHELLI G, BARAUSSE C, PIATTELLI M, PISTILLI R, GRANDI G, ESPOSITO M  
WIDE DIAMETER IMMEDIATE POST-EXTRACTIVE IMPLANTS VS DELAYED PLACEMENT OF NORMAL-DIAMETER IMPLANTS IN PRESERVED SOCKETS IN THE MOLAR REGION: 1-YEAR POST-LOADING OUTCOME OF A RANDOMISED CONTROLLED TRIAL  
EUR J ORAL IMPLANTOL, 2017;10(3):263-278



# OsteoBiol®: the most complete products range



The extensive OsteoBiol® range of products are engineered to help surgeons making the right decision when it comes to choose the perfect product for a specific clinical indication, both in dental and maxillofacial surgery.

Tecnoss® development of new products or improvement of existing ones, is focussed on supporting the technical capabilities of the practitioner to improve both intraoperative techniques and clinical results.

Specialists and researchers share their experience, blending clinical background and hands-on experience with the most advanced bio-technologies: the main goal is to obtain a specific solution to satisfy each clinical need.

OsteoBiol® collagenated grafting materials contribute to mineral deposition, vascular ingrowth and growth factor binding, thus providing a favourable environment for bone regeneration. The scientific literature has demonstrated that OsteoBiol® bone matrix is similar to human bone, and it has been reported in humans to be osteoconductive, well integrated in the host site and partially resorbed after 3-6 months, with no signs of adverse reaction<sup>(1)</sup>.

(1) GIULIANI A, IEZZI G, MAZZONI S, PIATTELLI A, PERROTTI V, BARONE A  
**REGENERATIVE PROPERTIES OF COLLAGENATED PORCINE BONE GRAFTS IN HUMAN MAXILLA: DEMONSTRATIVE STUDY OF THE KINETICS BY SYNCHROTRON RADIATION MICROTOMOGRAPHY AND LIGHT MICROSCOPY**  
CLIN. ORAL INVESTIG., 2017 2018 JAN;22(1):505-513



# CERTIFICATIONS





**APPROVAL CERTIFICATE Bone Matrix**  
Source: TecnoSS® s.r.l.



**APPROVAL CERTIFICATE Membranes**  
Source: TecnoSS® s.r.l.



**UNI EN ISO 13485 KIWA CERMET** quality certificate  
Source: TecnoSS® s.r.l.



**EC DESIGN-EXAMINATION CERTIFICATE Bone Matrix**  
Source: TecnoSS® s.r.l.



**EC DESIGN-EXAMINATION CERTIFICATE Membranes**  
Source: TecnoSS® s.r.l.





**FREE APP**  
for smartphone, iPhone, tablet  
and iPad including:

8 animations to show your patients  
the main GBR techniques

Information about the full range of  
OsteoBiol® biomaterials

Direct access to the database of videos  
and cases on osteobiol.com



This App may be too large to download over a mobile connection,  
or may exceed data usage limits. Wi-Fi connection recommended.





• international distributors

## Collagenated biomaterials

Distributed in 81 countries  
200 international scientific publications  
20 years of clinical success  
Over 1.000.000 surgeries performed

**OsteoBiol**<sup>®</sup>  
by Tecnos



REGENERATION SCIENCE

INSPIRED BY NATURE



# LITERATURE



**OsteoBiol®**  
by TecnoSS

SCIENTIFIC ABSTRACTS

REGENERATION SCIENCE

INSPIRED BY NATURE

Evolution:

Michela Casetta  
Laura Ricci  
Giovanna Izzzi  
Sabrina Calosso  
Adriano Pignelli  
Vittoria Perrotti

LATERAL ACCESS SINUS LIFT

ORIGINAL ARTICLE

International Journal of Periodontics &  
Oral Implants  
2012; 2(4): 120-126

filled with  
SUBSTITUTE  
BioBiol® Gen-Ox

LASL

filled with  
SUBSTITUTE  
BioBiol® Gen-Ox

LAE

filled with  
SUBSTITUTE  
BioBiol® mp2®

MEMBRANE  
OsteoBiol® Evolution

ALR

COVANI U, AMERI S, CRESPI R, BARONE A  
**PRESERVAZIONE DEL PROCESSO ALVEOLARE CON OSSO ETEROLOGO. CONSIDERAZIONI ISTOLOGICHE**  
ITALIAN ORAL SURGERY, 2004, VOL 3, 1: 17-23

CASSETTA M, CALASSO S, VOZZA I, DELL'AQUILA D  
**REHABILITATION OF ATROPHIC ALVEOLAR CRESTS WITH CYLINDRICAL SANDBLASTED AND ACID ETCHED IMPLANTS: A PILOT STUDY**  
EUR J IMPLANT PROSTHODONTICS, 2005;(3)1:133-144

ARCURI C, CECCHETTI F, GERMANO F, MOTTA A, SANTACROCE C  
**CLINICAL AND HISTOLOGICAL STUDY OF A XENOGENIC BONE SUBSTITUTE USED AS A FILLER IN POSTEXTRACTIVE ALVEOLUS**  
MINERVA STOMATOL, 2005 JUN;54(6):351-62

BARONE A, CRESPI R, ALDINI NN, FINI M, GIARDINO R, COVANI U  
**MAXILLARY SINUS AUGMENTATION: HISTOLOGIC AND HISTOMORPHOMETRIC ANALYSIS**  
INT J ORAL MAXILLOFAC IMPLANTS, 2005 JUL-AUG; 20(4):519-25

RINNA C, UNGARI C, SALTARELLA A, CASSONI A, REALE G  
**ORBITAL FLOOR RESTORATION**  
J CRANIOFAC SURG, 2005 NOV; 16(6):968-72

BARONE A, AMERI S, COVANI U  
**IMMEDIATE POSTEXTRACTION IMPLANTS: TREATMENT OF RESIDUAL PERI-IMPLANT DEFECTS. A RETROSPECTIVE ANALYSIS**  
EUR J IMPLANT PROSTHODONTICS, 2006,2: 99-106

BARONE A, SANTINI S, SBORDONE L, CRESPI R, COVANI U  
**A CLINICAL STUDY OF THE OUTCOMES AND COMPLICATIONS ASSOCIATED WITH MAXILLARY SINUS AUGMENTATION**  
INT J ORAL MAXILLOFAC IMPLANTS, 2006 JAN-FEB; 21(1):81-5

COVANI U, BARONE A, CORNELINI R, CRESPI R  
**CLINICAL OUTCOME OF IMPLANTS PLACED IMMEDIATELY AFTER IMPLANT REMOVAL**  
J PERIODONTOL, 2006 APR;77(4):722-7

ORSINI G, SCARANO A, PIATTELLI M, PICCIRILLI M, CAPUTI S, PIATTELLI A  
**HISTOLOGIC AND ULTRASTRUCTURAL ANALYSIS OF REGENERATED BONE IN MAXILLARY SINUS AUGMENTATION USING A PORCINE BONE-DERIVED BIOMATERIAL**  
J PERIODONTOL, 2006 DEC;77(12):1984-90

TRUBIANI O, SCARANO A, ORSINI G, DI IORIO D, D'ARCANGELO C, PICCIRILLI M, SIGISMONDO M, CAPUTI S  
**THE PERFORMANCE OF HUMAN PERIODONTAL LIGAMENT MESENCHYMAL STEM CELLS ON XENOGENIC BIOMATERIALS**  
INT J IMMUNOPATHOL PHARMACOL, 2007 JAN-MAR; 20 (1 SUPPL 1):87-91

BARONE A, COVANI U  
**MAXILLARY ALVEOLAR RIDGE RECONSTRUCTION WITH NON- VASCULARIZED AUTOGENOUS BLOCK BONE: CLINICAL RESULTS**  
J ORAL MAXILLOFAC SURG, 2007 OCT;65(10):2039-46

DEL CORSO M  
**SOFT TISSUE RESPONSE TO PLATELET RICH FIBRIN: CLINICAL EVIDENCES**  
COSMETIC DENTISTRY, 2008, 3:16-20

BARONE A, SANTINI S, MARCONCINI S, GIACOMELLI L, GHERLONE E, COVANI U  
**OSTEOTOMY AND MEMBRANE ELEVATION DURING THE MAXILLARY SINUS AUGMENTATION PROCEDURE. A COMPARATIVE STUDY: PIEZOELECTRIC DEVICE VS. CONVENTIONAL ROTATIVE INSTRUMENTS**  
CLIN ORAL IMPLANTS RES, 2008 MAY;19(5):511-5

BARONE A, CORNELINI R, CIAGLIA R, COVANI U  
**IMPLANT PLACEMENT IN FRESH EXTRACTION SOCKETS AND SIMULTANEOUS OSTEOTOME SINUS FLOOR ELEVATION: A CASE SERIES**  
INT J PERIODONTICS RESTORATIVE DENT, 2008 JUN; 28(3):283-9

BARONE A, ALDINI NN, FINI M, GIARDINO R, CALVO GUIRADO JL, COVANI U  
**XENOGRAFT VERSUS EXTRACTION ALONE FOR RIDGE PRESERVATION AFTER TOOTH REMOVAL: A CLINICAL AND HISTOMORPHOMETRIC STUDY**  
J PERIODONTOL, 2008 AUG;79(8):1370-7

COVANI U, CORNELINI R, BARONE A  
**BUCCAL BONE AUGMENTATION AROUND IMMEDIATE IMPLANTS WITH AND WITHOUT FLAP ELEVATION: A MODIFIED APPROACH**  
INT J ORAL MAXILLOFAC IMPLANTS, 2008 SEP-OCT; 23(5):841-6

CARDAROPOLI D, CARDAROPOLI G  
**PRESERVATION OF THE POSTEXTRACTION ALVEOLAR RIDGE: A CLINICAL AND HISTOLOGIC STUDY**  
INT J PERIODONTICS RESTORATIVE DENT, 2008 OCT; 28(5):469-77

NANNMARK U, SENNERBY L  
**THE BONE TISSUE RESPONSE TO PREHYDRATED AND COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS: A STUDY IN RABBIT MAXILLARY DEFECTS**  
CLIN IMPLANT DENT RELAT RES, 2008 DEC;10(4):264-70

SCARANO A, PIATTELLI M, CARINCI F, PERROTTI V  
**REMOVAL, AFTER 7 YEARS, OF AN IMPLANT DISPLACED INTO THE MAXILLARY SINUS. A CLINICAL AND HISTOLOGIC CASE REPORT**  
J OSSEOINTEGR, 2009;1(1):35-40

COVANI U, MARCONCINI S, CRESPI R, BARONE A  
**IMMEDIATE IMPLANT PLACEMENT AFTER REMOVAL OF A FAILED IMPLANT: A CLINICAL AND HISTOLOGICAL CASE REPORT**  
J ORAL IMPLANTOL, 2009; 35(4):189-95

FIGUEIREDO M, HENRIQUES J, MARTINS G, GUERRA F, JUDAS F, FIGUEIREDO H  
**PHYSICO-CHEMICAL CHARACTERIZATION OF BIOMATERIALS COMMONLY USED IN DENTISTRY AS BONE SUBSTITUTES - COMPARISON WITH HUMAN BONE**  
J BIOMED MATER RES B APPL BIOMATER, 2010FEB; 92(2):409-19

GRENGA PL, REALE G, COFONE C, MEDURI A, CERUTI P, GRENGA R  
**HESS AREA RATIO AND DIPLOPIA: EVALUATION OF 30 PATIENTS UNDERGOING SURGICAL REPAIR FOR ORBITAL BLOW-OUT FRACTURE**  
OPHTHAL PLAST RECONSTR SURG, 2009 MAR-APR; 25(2):123-5

CRESPI R, CAPPARÈ P, GHERLONE E  
**DENTAL IMPLANTS PLACED IN EXTRACTION SITES GRAFTED WITH DIFFERENT BONE SUBSTITUTES: RADIOGRAPHIC EVALUATION AT 24 MONTHS**  
J PERIODONTOL, 2009 OCT; 80(10):1616-1621

RINNA C, REALE G, FORESTA E, MUSTAZZA MC  
**MEDIAL ORBITAL WALL RECONSTRUCTION WITH SWINE BONE CORTEX**  
J CRANIOFAC SURG, 2009 MAY; 20(3): 881-4

CARDAROPOLI D, CARDAROPOLI G  
**HEALING OF GINGIVAL RECESSIONS USING A COLLAGEN MEMBRANE WITH A THE MINERALIZED XENOGRAFT: A RANDOMIZED CONTROLLED CLINICAL TRIAL**  
INT J PERIODONTICS RESTORATIVE DENT, 2009 FEB; 29(1):59-67

NANNMARK U, AZARMEHR I  
**SHORT COMMUNICATION: COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS. A STUDY IN RABBIT MAXILLARY DEFECTS**  
CLIN IMPLANT DENT RELAT RES, 2010 JUN 1; 12(2):161-3

SCARANO A, PIATTELLI A, ASSENZA B, QUARANTA A, PERROTTI V, PIATTELLI M, IEZZI G  
**PORCINE BONE USED IN SINUS AUGMENTATION PROCEDURES: A 5-YEAR RETROSPECTIVE CLINICAL EVALUATION**  
J ORAL MAXILLOFAC SURG, 2010 AUG; 68(8):1869-73

ROSSI R, MORALES RS, FRASCARIA M, BENZI R, SQUADRITO N  
**PLANNING IMPLANTS IN THE ESTHETIC ZONE USING A NEW IMPLANT 3D NAVIGATION SYSTEM**  
EUR J ESTHETIC DENT, 2010 SUMMER; 5(2):172-88

SCARANO A, CARINCI F, ASSENZA B, PIATTELLI M, MURMURA G, PIATTELLI A  
**VERTICAL RIDGE AUGMENTATION OF ATROPHIC POSTERIOR MANDIBLE USING AN INLAY TECHNIQUE WITH A XENOGRAFT WITHOUT MINISCREWS AND MINIPLATES: CASE SERIES**  
CLIN ORAL IMPLANTS RES, 2011 OCT;22(10):1125-30

PAGLIANI L, ANDERSSON P, LANZA M, NAPPO A, VERROCCHI D, VOLPE S, SENNERBY L  
**A COLLAGENATED PORCINE BONE SUBSTITUTE FOR AUGMENTATION AT NEOS IMPLANT SITES: A PROSPECTIVE 1-YEAR MULTICENTER CASE SERIES STUDY WITH HISTOLOGY**  
CLIN IMPLANT DENT RELAT RES, 2012 OCT;14(5):746-58

SANTAGATA M, GUARINIELLO L, TARTARO G  
**A MODIFIED EDENTULOUS RIDGE EXPANSION (MERE) TECHNIQUE FOR IMMEDIATE PLACEMENT OF IMPLANTS. A CASE REPORT**  
J ORAL IMPLANTOL, 2011 MAR;37 SPEC N°:114-9

SCARANO A, PIATTELLI A, PERROTTI V, MANZON I, IEZZI G  
**MAXILLARY SINUS AUGMENTATION IN HUMANS USING CORTICAL PORCINE BONE: A HISTOLOGICAL AND HISTOMORPHOMETRIC EVALUATION AFTER 4 AND 6 MONTHS**  
CLIN IMPLANT DENT RELAT RES, 2011 MAR;13(1):13-18

CRESPI R, CAPPARÈ P, ROMANOS GE, MARIANI E, BENASCIIUTTI E, GHERLONE E  
**CORTICOCANCELLOUS PORCINE BONE IN THE HEALING OF HUMAN EXTRACTION SOCKETS: COMBINING HISTOMORPHOMETRY WITH OSTEOBLAST GENE EXPRESSION PROFILES IN VIVO**  
INT J ORAL MAXILLOFAC IMPLANTS, 2011 JUL-AUG; 26(4):866-72

HINZE M, VRIELINCK L, THALMAIR T, WACHTEL H, BOLZ W  
**ZYGOMATIC IMPLANT PLACEMENT IN CONJUNCTION WITH SINUS BONE GRAFTING: THE "EXTENDED SINUS ELEVATION TECHNIQUE". A CASE-COHORT STUDY**  
ORAL CRANIOFAC TISSUE ENG, 2011; 1:188-197

IEZZI G, DEGIDI M, PIATTELLI A, MANGANO C, SCARANO A, SHIBLI JA, PERROTTI V  
**COMPARATIVE HISTOLOGICAL RESULTS OF DIFFERENT BIOMATERIALS USED IN SINUS AUGMENTATION PROCEDURES: A HUMAN STUDY AT 6 MONTHS**  
CLIN ORAL IMPLANTS RES, 2012 DEC;23(12):1369-76

SLOTTE C, LINDFORS N, NANNMARK U  
**SURGICAL RECONSTRUCTION OF PERI-IMPLANT BONE DEFECTS WITH PREHYDRATED AND COLLAGENATED PORCINE BONE AND COLLAGEN BARRIERS: CASE PRESENTATIONS**  
CLIN IMPLANT DENT RELAT RES, 2013 OCT;15(5):714-23

BARONE A, RICCI M, GRASSI RF, NANNMARK U, QUARANTA A, COVANI U  
**A 6-MONTH HISTOLOGICAL ANALYSIS ON MAXILLARY SINUS AUGMENTATION WITH AND WITHOUT USE OF COLLAGEN MEMBRANES OVER THE OSTEOTOMY WINDOW: RANDOMIZED CLINICAL TRIAL**  
CLIN ORAL IMPLANTS RES, 2013 JAN; 24(1):1-6

SANTAGATA M, GUARINIELLO L, RAUSO R, TARTARO G  
**IMMEDIATE LOADING OF DENTAL IMPLANT AFTER SINUS FLOOR ELEVATION WITH OSTEOTOME TECHNIQUE: A CLINICAL REPORT AND PRELIMINARY RADIOGRAPHIC RESULTS**  
J ORAL IMPLANTOL, 2010 DEC; 36(6):485-489

FESTA VM, ADDABBO F, LAINO L, FEMIANO F, RULLO R  
**PORCINE-DERIVED XENOGRAFT COMBINED WITH A SOFT CORTICAL MEMBRANE VERSUS EXTRACTION ALONE FOR IMPLANT SITE DEVELOPMENT: A CLINICAL STUDY IN HUMANS**  
CLIN IMPLANT DENT RELAT RES, 2013 OCT;15(5):707-13

RAMIREZ FERNANDEZ MP, CALVO GUIRADO JL, MATÉ SANCHEZ DE VAL JE, DELGADO RUIZ RA, NEGRI B, BARONA DORADO C  
**ULTRASTRUCTURAL STUDY BY BACKSCATTERED ELECTRON IMAGING AND ELEMENTAL MICROANALYSIS OF BONE-TO-BIOMATERIAL INTERFACE AND MINERAL DEGRADATION OF PORCINE XENOGRAFTS USED IN MAXILLARY SINUS FLOOR ELEVATION**  
CLIN ORAL IMPLANTS RES, 2013 MAY;24(5):523-30

CASSETTA M, RICCI L, IEZZI G, DELL'AQUILA D, PIATTELLI A, PERROTTI V  
**RESONANCE FREQUENCY ANALYSIS OF IMPLANTS INSERTED WITH A SIMULTANEOUS GRAFTING PROCEDURE: A 5-YEAR FOLLOW-UP STUDY IN MAN**  
INT J PERIODONTICS RESTORATIVE DENT, 2012 OCT;32(5):581-9

BARONE A, ORLANDO B, CINGANO L, MARCONCINI S, DERCHI G, COVANI U  
**A RANDOMIZED CLINICAL TRIAL TO EVALUATE AND COMPARE IMPLANTS PLACED IN AUGMENTED VS. NON-AUGMENTED EXTRACTION SOCKETS. A 3-YEAR EVALUATION**  
J PERIODONTOL, 2012 JUL;83(7):836-46



ESPOSITO M, CANNIZZARO G, SOARDI E, PISTILLI R, PIATTELLI M, CORVINO V, FELICE P  
**POSTERIOR ATROPHIC JAWS REHABILITATED WITH PROSTHESES SUPPORTED BY 6 MM-LONG, 4 MM-WIDE IMPLANTS OR BY LONGER IMPLANTS IN AUGMENTED BONE. PRELIMINARY RESULTS FROM A PILOT RANDOMISED CONTROLLED TRIAL**  
 EUR J ORAL IMPLANTOL, 2012 SPRING;5(1):19-33

FELICE P, PIANA L, CHECCHI L, PISTILLI R, PELLEGRINO G  
**VERTICAL RIDGE AUGMENTATION OF THE ATROPHIC POSTERIOR MANDIBLE WITH A 2-STAGE INLAY TECHNIQUE: A CASE REPORT**  
 IMPLANT DENT, 2012 JUN;21(3):190-5

BARONE A, RICCI M, TONELLI P, SANTINI S, COVANI U  
**TISSUE CHANGES OF EXTRACTION SOCKETS IN HUMANS: A COMPARISON OF SPONTANEOUS HEALING VS. RIDGE PRESERVATION WITH SECONDARY SOFT TISSUE HEALING**  
 CLIN ORAL IMPLANTS RES, 2013 NOV;24(11):1231-7

CASSETTA M, RICCI L, IEZZI G, CALASSO S, PIATTELLI A, PERROTTI V  
**USE OF PIEZOSURGERY DURING MAXILLARY SINUS ELEVATION: CLINICAL RESULTS OF 40 CONSECUTIVE CASES**  
 INT J PERIODONTICS RESTORATIVE DENT, 2012 DEC;32(6):E182-8

BRUNELLI G, SOLLAZZO V, CARINCI F, PALMIERI A, GIRARDI A, MONGUZZI R  
**OSTEOBIOL® INFLUENCES OSTEOGENIC DIFFERENTIATION OF ADIPOSE DERIVED STEM CELLS**  
 EUR J INFLAMMAT, 2011, VOL. 9, NO. 3 (S), 103-107

FELICE P, PIANA L, CHECCHI L, CORVINO V, NANNMARK U, PIATTELLI M  
**VERTICAL RIDGE AUGMENTATION OF ATROPHIC POSTERIOR MANDIBLE WITH AN INLAY TECHNIQUE AND CANCELLOUS EQUINE BONE BLOCK: A CASE REPORT**  
 INT J PERIODONTICS RESTORATIVE DENT, 2013 MAR;33(2):159-66

FICKL S, JOCKEL-SCHNEIDER Y, LINCKE T, BECHTOLD M, FISCHER KR, SCHLAGENHAUF U  
**PORCINE DERMAL MATRIX FOR COVERING OF RESSION TYPE DEFECTS: A CASE SERIES**  
 QUINTESSENCE INT, 2013;44(3):243-6

SILVESTRI M, MARTEGANI P, D'AVENIA F, FARNETI M, CAPRI D, PAOLANTONI G, LANDI L  
**SIMULTANEOUS SINUS AUGMENTATION WITH IMPLANT PLACEMENT: HISTOMORPHOMETRIC COMPARISON OF TWO DIFFERENT GRAFTING MATERIALS. A MULTICENTER DOUBLE-BLIND PROSPECTIVE RANDOMIZED CONTROLLED CLINICAL TRIAL**  
 INT J ORAL MAXILLOFAC IMPLANTS, 2013 MAR-APR;28(2):543-9

WACHTEL H, FICKL S, HINZE M, BOLZ W, THALMAIR T  
**THE BONE LAMINA TECHNIQUE: A NOVEL APPROACH FOR LATERAL RIDGE AUGMENTATION - A CASE SERIES**  
 INT J PERIODONTICS RESTORATIVE DENT, 2013 JUL-AUG;33(4):491-7

RODRIGUEZ JG, ELIDIBANY RM  
**VERTICAL SPLITTING OF THE MANDIBULAR BODY AS AN ALTERNATIVE TO INFERIOR ALVEOLAR NERVE LATERALIZATION**  
 INT J ORAL MAXILLOFAC SURG, 2013 SEP;42(9):1060-6

FIGUEIREDO A, COIMBRA P, CABRITA A, GUERRA F, FIGUEIREDO M  
**COMPARISON OF A XENOGENEIC AND AN ALLOPLASTIC MATERIAL USED IN DENTAL IMPLANTS IN TERMS OF PHYSICO-CHEMICAL CHARACTERISTICS AND IN VIVO INFLAMMATORY RESPONSE**  
 MATER SCI ENG C, MATER BIOL APP, 2013 AUG 1;33(6):3506-13

FELICE P, PISTILLI R, PIATTELLI M, SOARDI E, CORVINO V, ESPOSITO M  
**POSTERIOR ATROPHIC JAWS REHABILITATED WITH PROSTHESES SUPPORTED BY 5 X 5 MM IMPLANTS WITH A NOVEL NANOSTRUCTURED CALCIUM-INCORPORATED TITANIUM SURFACE OR BY LONGER IMPLANTS IN AUGMENTED BONE. PRELIMINARY RESULTS FROM A RANDOMISED CONTROLLED TRIAL**  
 EUR J ORAL IMPLANTOL, SUMMER;5(2):149-61

TRAINI T, PIATTELLI A, CAPUTI S, DEGIDI M, MANGANO C, SCARANO A, PERROTTI V, IEZZI G  
**REGENERATION OF HUMAN BONE USING DIFFERENT BONE SUBSTITUTE BIOMATERIALS**  
 CLIN IMPLANT DENT RELAT RES, 2015 FEB;17(1):150-62

KOLMAS J, SZWAJA M, KOLODZIEJSKI W  
**SOLID-STATE NMR AND IR CHARACTERIZATION OF COMMERCIAL XENOGENEIC BIOMATERIALS USED AS BONE SUBSTITUTES**  
 J PHARM BIOMED ANAL, 2012 MAR 5;61:136-41

PISTILLI R, FELICE P, PIATTELLI M, GESSAROLI M, SOARDI E, BARAUSSE C, BUTI J, CORVINO V, ESPOSITO M  
**POSTERIOR ATROPHIC JAWS REHABILITATED WITH PROSTHESES SUPPORTED BY 5 X 5 MM IMPLANTS WITH A NOVEL NANOSTRUCTURED CALCIUM-INCORPORATED TITANIUM SURFACE OR BY LONGER IMPLANTS IN AUGMENTED BONE. ONE-YEAR RESULTS FROM A RANDOMISED CONTROLLED TRIAL**  
 EUR J ORAL IMPLANTOL, 2013 WINTER;6(4):343-357

FICKL S, NANNMARK U, SCHLAGENHAUF U, HÜRZELER M, KEBSCHULL M  
**PORCINE DERMAL MATRIX IN THE TREATMENT OF DEHISCENCE-TYPE DEFECTS - AN EXPERIMENTAL SPLIT-MOUTH ANIMAL TRIAL**  
 CLIN ORAL IMPLANTS RES, 2015 JUL;26(7):799-805

BARONE A, BORGIA V, COVANI U, RICCI M, PIATTELLI A, IEZZI G  
**FLAP VERSUS FLAPLESS PROCEDURE FOR RIDGE PRESERVATION IN ALVEOLAR EXTRACTION SOCKETS: A HISTOLOGICAL EVALUATION IN A RANDOMIZED CLINICAL TRIAL**  
 CLIN ORAL IMPLANTS RES, 2015 JUL;26(7):806-13

BARONE A, RICCI M, ROMANOS GE, TONELLI P, ALFONSI F, COVANI U  
**BUCCAL BONE DEFICIENCY IN FRESH EXTRACTION SOCKETS: A PROSPECTIVE SINGLE COHORT STUDY**  
 CLIN ORAL IMPLANTS RES, 2015 JUL;26(7):823-30

BARONE A, TOTI P, QUARANTA A, DERCHI G, COVANI U  
**THE CLINICAL OUTCOMES OF IMMEDIATE VERSUS DELAYED RESTORATION PROCEDURES ON IMMEDIATE IMPLANTS: A COMPARATIVE COHORT STUDY FOR SINGLE-TOOTH REPLACEMENT**  
 CLIN IMPLANT DENT RELAT RES, 2015 DEC;17(6):1114-26

CASSETTA M, PERROTTI V, CALASSO S, PIATTELLI A, SINJARI B, IEZZI G  
**BONE FORMATION IN SINUS AUGMENTATION PROCEDURES USING AUTOLOGOUS BONE, PORCINE BONE, AND A 50 : 50 MIXTURE: A HUMAN CLINICAL AND HISTOLOGICAL EVALUATION AT 2 MONTHS**  
 CLIN ORAL IMPLANTS RES, 2015 OCT;26(10):1180-4

GHEÑO E, PALERMO A, BUFFOLI B, RODELLA LF  
**THE EFFECTIVENESS OF THE USE OF XENOGENEIC BONE BLOCKS MIXED WITH AUTOLOGOUS CONCENTRATED GROWTH FACTORS (CGF) IN BONE REGENERATION TECHNIQUES**  
 J OSSEOINTEGRATION 2014;6(2):37-42

FALISI G, GALLI M, VITTORINI-VELASQUEZ P, GALLEGOS-RIVERA JC, MINASI R, DE BIAZE A, DI PAOLO C  
**USE OF 3D CARTILAGE SCAFFOLDS FOR THE STABILIZATION OF IMPLANTS AND BONE REGENERATION WITH THE FIT-LOCK TECHNIQUE**  
 ACTA ODONTOL LATINOAM 2013;26(3):167-172

FISCHER KR, FICKL S, MARDAS N, BOZEC L, DONOS N  
**STAGE-TWO SURGERY USING COLLAGEN SOFT TISSUE GRAFTS: CLINICAL CASES AND ULTRASTRUCTURAL ANALYSIS**  
 QUINTESSENCE INT, 2014 NOV-DEC; 45(10):853-60

FISCHER KR, STAVROPOULOS A, CALVO GUIRADO JL, SCHNEIDER D, FICKL S  
**INFLUENCE OF LOCAL ADMINISTRATION OF PAMIDRONATE ON EXTRACTION SOCKET HEALING - A HISTOMORPHOMETRIC PROOF-OF-PRINCIPLE PRE-CLINICAL IN VIVO EVALUATION**  
 CLIN ORAL IMPLANTS RES, 2015 OCT;26(10):1135-42


SCARANO A, MURMURA G, SINJARI B, ASSENZA B, SOLLAZZO V, SPINELLI G, CARINCI F  
**EXPANSION OF THE ALVEOLAR BONE CREST WITH ULTRASONIC SURGERY DEVICE: CLINICAL STUDY IN MANDIBLE**  
 INT J IMMUNOPATHOL PHARMACOL, 2011 APR-JUN; 24(2 SUPPL):71-5

SCARANO A, PIATTELLI A, MURMURA G, IEZZI G, ASSENZA B, MANCINO C  
**DELAYED EXPANSION OF THE ATROPHIC MANDIBLE BY ULTRASONIC SURGERY: A CLINICAL AND HISTOLOGICAL CASE SERIES**  
 INT J ORAL MAXILLOFAC IMPLANTS, 2015 JAN-FEB;30(1):144-9

LORENZON G, BUTTARELLO GM, CHESSA G  
**CASE REPORT: IMPLANT PLACEMENT AND IMMEDIATE LOADING WITH SIMULTANEOUS BONE REGENERATION FOLLOWING JAW ODONTOGENIC CYST ENUCLEATION**  
 DENTISTRY, 2015, 5:2

THALMAIR T, FICKL S, SCHNEIDER D, HINZE M, WACHTEL H  
**DIMENSIONAL ALTERATIONS OF EXTRACTION SITES AFTER DIFFERENT ALVEOLAR RIDGE PRESERVATION TECHNIQUES - A VOLUMETRIC STUDY**  
 J CLIN PERIODONTOL, 2013 JUL;40(7):721-7

MANESCU A, GIULIANI A, MOHAMMADI S, TROMBA G, MAZZONI S, DIOMEDE F, ZINI N, PIATTELLI A, TRUBIANI O  
**OSTEOGENIC POTENTIAL OF DUAL-BLOCKS CULTURED WITH HUMAN PERIODONTAL LIGAMENT STEM CELLS: IN VITRO AND SYNCHROTRON**  
 J PERIODONTAL RES, 2016 FEB;51(1):112-24



SEM image of OsteoBioL® cancellous block  
 Author: Prof Ulf Nannmark, University of Göteborg, Sweden



SCARANO A, PIATTELLI A, IEZZI G, VARVARA G  
**SPONTANEOUS BONE FORMATION ON THE MAXILLARY SINUS FLOOR IN ASSOCIATION WITH SURGERY TO REMOVE A MIGRATED DENTAL IMPLANT: A CASE REPORT**

MINERVA STOMATOL, 2014 OCT;63(10):351-9

BARONE A, TOTI P, QUARANTA A, ALFONSI F, CUCCHI A, CALVO GUIRADO JL, NEGRI B, DI FELICE R, COVANI U  
**VOLUMETRIC ANALYSIS OF REMODELLING PATTERN AFTER RIDGE PRESERVATION COMPARING USE OF TWO TYPES OF XENOGRAPHS. A MULTICENTRE RANDOMIZED CLINICAL TRIAL**

CLIN ORAL IMPLANTS RES, 2016 NOV;27(11):E105-E115

ESPOSITO M, GRUSOVIN MG, LAMBERT F, MATOS S, PIETRUSKA M, ROSSI R, SALHI L, BUTI J  
**THE EFFECTIVENESS OF A RESORBABLE BONE SUBSTITUTE WITH A RESORBABLE MEMBRANE IN THE TREATMENT OF PERIODONTAL INFRABONY DEFECT - A MULTICENTER RANDOMISED CONTROLLED TRIAL**

EUR J ORAL IMPLANTOL, 2015;8(3):233-244

OZEL B, FINDIKCIOGLU K, SEZGIN B, GUNAY K, BARUT I, OZMEN S

**A NEW OPTION FOR THE RECONSTRUCTION OF ORBITAL FLOOR DEFECTS WITH HETEROLOGOUS CORTICAL BONE**

J CRANIO-MAXILLOFAC SURG, 2015 OCT;43(8):1583-8

CORBELLA S, TASCHIERI S, WEINSTEIN R, DEL FABBRO M  
**HISTOMORPHOMETRIC OUTCOMES AFTER LATERAL SINUS FLOOR ELEVATION PROCEDURE: A SYSTEMATIC REVIEW OF THE LITERATURE AND META-ANALYSIS**

CLIN ORAL IMPLANTS RES, 2016 SEP;27(9):1106-22

BARONE A, MARCONCINI S, GIAMMARINARO E, MIJIRITSKY E, GELPI F, COVANI U

**CLINICAL OUTCOMES OF IMPLANTS PLACED IN EXTRACTION SOCKETS AND IMMEDIATELY RESTORED: A 7-YEAR SINGLE-COHORT PROSPECTIVE STUDY**

CLIN IMPLANT DENT RELAT RES, 2016 DEC;18(6):1103-1112

CAKIR M, KARACA IR, AYŞEGÜL F, KAYMAZ F, BOZKAYA S  
**EXPERIMENTAL EVALUATION OF THE EFFECTS OF ANKAFERD BLOOD STOPPER AND COLLAGENATED HETEROLOGOUS BONE GRAFT ON BONE HEALING IN SINUS FLOOR AUGMENTATION**

CLIN ORAL IMPLANTS RES, 2015 MAR-APR;30(2):279-85

LOPEZ MA, ANDREASI BASSI M, CONFALONE L, CARINCI F  
**REGENERATION OF ATROPHIC CRESTAL RIDGES WITH RESORBABLE LAMINA: TECHNICAL NOTE**

J BIOL REGUL HOMEOST AGENTS 2015 JUL-SEP;29(3 SUPPL 1):97-100

ETTORRE V, DE MARCO P, ZARA S, PERROTTI V, SCARANO A, DI CRESCENZO A, PETRINI M, HADAD C, BOSCO D, ZAVAN B, VALBONETTI L, SPOTO G, IEZZI G, PIATTELLI A, CATALDI A, FONTANA A

**IN VITRO AND IN VIVO CHARACTERIZATION OF GRAPHENE OXIDE COATED PORCINE BONE GRANULES**

CARBON, JULY 2016, VOLUME 103, PAGES 291-298

ROSSI R, RANCITELLI D, POLI PP, RASIA DAL POLO M, NANNMARK U, MAIORANA C

**THE USE OF A COLLAGENATED PORCINE CORTICAL LAMINA IN THE RECONSTRUCTION OF ALVEOLAR RIDGE DEFECTS. A CLINICAL AND HISTOLOGICAL STUDY**

MINERVA STOMATOL, 2016 OCT;65(5):257-68

SCARANO A, LORUSSO F, RAVERA L, MORTELLARO C, PIATTELLI A  
**BONE REGENERATION IN ILIAC CRESTAL DEFECTS: AN EXPERIMENTAL STUDY ON SHEEP**

BIOMED RES INT, 2016;2016:4086870

FELICE P, ZUCCHELLI G, CANNIZZARO G, BARAUSSE C, DIAZZI M, TRULLENQUE-ERIKSSON A, ESPOSITO M  
**IMMEDIATE, IMMEDIATE-DELAYED (6 WEEKS) AND DELAYED (4 MONTHS) POST-EXTRACTIVE SINGLE IMPLANTS: 4-MONTH POST-LOADING DATA FROM A RANDOMISED CONTROLLED TRIAL**

EUR J ORAL IMPLANTOL, 2016;9(3):233-247

ROMBOUS C, JEANNEAU C, CAMILLERI J, LAURENT P, ABOUT I  
**CHARACTERIZATION AND ANGIOGENIC POTENTIAL OF XENOGENEIC BONE GRAFTING MATERIALS: ROLE OF PERIODONTAL LIGAMENT CELLS**

DENT MATER J, 2016 DEC 1;35(6):900-907

BARONE A, TOTI P, MARCONCINI S, DERCHI G, MARCHIONNI S, COVANI U  
**ESTHETIC OUTCOME OF IMPLANTS PLACED IN FRESH EXTRACTION SOCKETS BY CLINICIANS WITH OR WITHOUT EXPERIENCE: A MEDIUM-TERM PROSPECTIVE EVALUATION**

INT J ORAL MAXILLOFAC IMPLANTS, 2016;31(6)

TALLARICO M, KHANARI E, PISANO M, DE RIU G, TULLIO A, MELONI SM

**SINGLE POST-EXTRACTIVE ULTRA-WIDE 7 MM-DIAMETER IMPLANTS VERSUS IMPLANTS PLACED IN MOLAR HEALED SITES AFTER SOCKET PRESERVATION FOR MOLAR REPLACEMENT: 6-MONTH POST-LOADING RESULTS FROM A RANDOMISED CONTROLLED TRIAL**

EUR J ORAL IMPLANTOL, 2016;9(3):263-275

LOPEZ MA, MANZULLI N, CASALE M, ORMIANER Z, CARINCI F

**THE USE OF RESORBABLE HETEROLOGOUS CORTICAL LAMINA AS A NEW SINUS LIFT FLOOR: A TECHNICAL NOTE**

J BIOL REGUL HOMEOST AGENTS, 2016 APR-JUN;30(2 SUPPL 1):75-79

LOPEZ MA, ANDREASI BASSI M, CONFALONE L, CARINCI F, ORMIANER Z, LAURITANO D

**THE USE OF RESORBABLE CORTICAL LAMINA AND MICRONIZED COLLAGENATED BONE IN THE REGENERATION OF ATROPHIC CRESTAL RIDGES: A SURGICAL TECHNIQUE. CASE SERIES**

J BIOL REGUL HOMEOST AGENTS, 2016 APR-JUN;30(2 SUPPL 1):81-85

ESPOSITO M, ZUCCHELLI G, BARAUSSE C, PISTILLI R, TRULLENQUE-ERIKSSON A, FELICE P

**FOUR MM-LONG VERSUS LONGER IMPLANTS IN AUGMENTED BONE IN POSTERIOR ATROPHIC JAWS: 4-MONTH POST-LOADING RESULTS FROM A RANDOMISED CONTROLLED TRIAL**

EUR J ORAL IMPLANTOL, 2016;9(4):393-409

BARONE A, TOTI P, QUARANTA A, ALFONSI F, CUCCHI A, NEGRI B, DI FELICE R, MARCHIONNI S, CALVO GUIRADO JL, COVANI U, NANNMARK U

**CLINICAL AND HISTOLOGICAL CHANGES AFTER RIDGE PRESERVATION WITH TWO XENOGRAPHS: PRELIMINARY RESULTS FROM A MULTICENTER RANDOMIZED CONTROLLED CLINICAL TRIAL**

J CLIN PERIODONTOL, 2017 FEB;44(2):204-214

ALFONSI F, BORGIA V, IEZZI G, PIATTELLI A, COVANI U, TONELLI P, BARONE A  
**MOLECULAR, CELLULAR AND PHARMACEUTICAL ASPECTS OF FILLING BIOMATERIALS DURING THE MANAGEMENT OF EXTRACTION SOCKETS**

CURR PHARM BIOTECHNOL, 2017;18(1):64-75

BARONE A, TOTI P, MENCHINI FABRIS GB, MARCHIONNI S, COVANI U  
**EARLY VOLUMETRIC CHANGES AFTER VERTICAL AUGMENTATION OF THE ATROPHIC POSTERIOR MANDIBLE WITH INTERPOSITIONAL BLOCK GRAFT VERSUS ONLAY BONE GRAFT: A RETROSPECTIVE RADIOLOGICAL STUDY**

J CRANIO-MAXILLOFAC, 2017 SEP;45(9):1438-1447

BARONE A, TOTI P, FUNEL N, CAMPANI D, COVANI U  
**EXPRESSION OF SP7, RUNX1, DLX5, AND CTNNB1 IN HUMAN MESENCHYMAL STEM CELLS CULTURED ON XENOGENEIC BONE SUBSTITUTE AS COMPARED WITH MACHINED TITANIUM**

IMPLANT DENT, 2014 AUG;23(4):407-15

ESPOSITO M, ZUCCHELLI G, CANNIZZARO G, CHECCHI L, BARAUSSE C, TRULLENQUE-ERIKSSON, FELICE P  
**IMMEDIATE, IMMEDIATE-DELAYED (6 WEEKS) AND DELAYED (4 MONTHS) POST-EXTRACTIVE SINGLE IMPLANTS: 1-YEAR POST-LOADING DATA FROM A RANDOMISED CONTROLLED TRIAL**

EUR J ORAL IMPLANTOL, 2017;10(1):11-26

SCARANO A, CRINCOLI V, DI BENEDETTO A, COZZOLINO V, LORUSSO F, PODALIRI VULPIANI M, GRANO M, KALEMAJ Z, MORI G, GRASSI FR

**BONE REGENERATION INDUCED BY BONE PORCINE BLOCK WITH BONE MARROW STROMAL STEM CELLS IN A MINIPIG MODEL OF MANDIBULAR "CRITICAL SIZE" DEFECT**

STEM CELLS INT, 2017;2017:9082869

SCARANO A  
**TRADITIONAL POSTEXTRACTIVE IMPLANT SITE PREPARATION COMPARED WITH PRE-EXTRACTIVE INTERRADICULAR IMPLANT BED PREPARATION IN THE MANDIBULAR MOLAR REGION, USING AN ULTRASONIC DEVICE: A RANDOMIZED PILOT STUDY**

INT J ORAL MAXILLOFAC IMPLANTS, 2017 MAY/JUN;32(3):655-660

BARONE A, TOTI P, MENCHINI-FABRIS GB, DERCHI G, MARCONCINI S, COVANI U  
**EXTRA ORAL DIGITAL SCANNING AND IMAGING SUPERIMPOSITION FOR VOLUME ANALYSIS OF BONE REMODELING AFTER TOOTH EXTRACTION WITH AND WITHOUT 2 TYPES OF PARTICULATE PORCINE MINERAL INSERTION: A RANDOMIZED CONTROLLED TRIAL**

CLIN IMPLANT DENT RELAT RES, 2017 AUG;19(4):750-759

GIULIANI A, IEZZI G, MAZZONI S, PIATTELLI A, PERROTTI V, BARONE A  
**REGENERATIVE PROPERTIES OF COLLAGENATED PORCINE BONE GRAFTS IN HUMAN MAXILLA: DEMONSTRATIVE STUDY OF THE KINETICS BY SYNCHROTRON RADIATION MICROTOMOGRAPHY AND LIGHT MICROSCOPY**

CLINICAL ORAL INVESTIGATIONS, 2018 JAN;22(1):505-513

IEZZI G, PIATTELLI A, GIULIANI A, MANGANO C, BARONE A, MANZONI L, DEGIDI M, SCARANO A, FILIPPONE A, PERROTTI V  
**MOLECULAR, CELLULAR AND PHARMACEUTICAL ASPECTS OF FILLING BIOMATERIALS DURING MAXILLARY SINUS-LIFT PROCEDURES. PART 2: DETAILED CHARACTERISTICS OF THE MATERIALS**

CURR PHARM BIOTECHNOL, 2017, 18, 33-44

FELICE P, BARAUSSE C, BARONE A, ZUCCHELLI G, PIATTELLI M, PISTILLI R, IPPOLITO DR, SIMION M  
**INTERPOSITIONAL AUGMENTATION TECHNIQUE IN THE TREATMENT OF POSTERIOR MANDIBULAR ATROPHIES: A RETROSPECTIVE STUDY COMPARING 129 AUTOGENOUS AND HETEROLOGOUS BONE BLOCKS WITH 2 TO 7 YEARS FOLLOW-UP**

INT J PERIODONTICS RESTORATIVE DENT, 2017 JUL/AUG;37(4):469-480

IIDA T, CARNEIRO MARTINS NETO E, BOTTICELLI D, APAZA ALCCAYHUAMAN KA, LANG NP, XAVIER SP  
**INFLUENCE OF A COLLAGEN MEMBRANE POSITIONED SUBJACENT THE SINUS MUCOSA FOLLOWING THE ELEVATION OF THE MAXILLARY SINUS. A HISTOMORPHOMETRIC STUDY IN RABBITS**

CLIN IMPLANT DENT RELAT RES, 2017 JUN 7

DE MARCO P, ZARA S, DE COLLI M, RADUNOVIC M, LAZOVIC V, ETTORRE V, DI CRESCENZO A, PIATTELLI A, CATALDI A, FONTANA A  
**GRAPHENE OXIDE IMPROVES THE BIOCOMPATIBILITY OF COLLAGEN MEMBRANES IN AN IN VITRO MODEL OF HUMAN PRIMARY GINGIVAL FIBROBLASTS**

BIOMED MATER, 2017 SEP 13;12(5):055005

MIJIRITSKY E, FERRONI L, GARDIN C, BRESSAN E, ZANETTE G, PIATTELLI A, ZAVAN B

**PORCINE BONE SCAFFOLDS ADSORB GROWTH FACTORS SECRETED BY MSCS AND IMPROVE BONE TISSUE REPAIR**

MATERIALS, 2017 SEP 8;10(9)

ROSSI R, FOCE E, SCOLAVINO S  
**THE CORTICAL LAMINA TECHNIQUE: A NEW OPTION FOR ALVEOLAR RIDGE AUGMENTATION. PROCEDURE, PROTOCOL, AND CASE REPORT**

J LEBANESE DENTAL ASS, 2017 JAN-JUN; 52(1):35-41

CHECCHI V, FELICE P, ZUCCHELLI G, BARAUSSE C, PIATTELLI M, PISTILLI R, GRANDI G, ESPOSITO M  
**WIDE DIAMETER IMMEDIATE POST-EXTRACTIVE IMPLANTS VS DELAYED PLACEMENT OF NORMAL-DIAMETER IMPLANTS IN PRESERVED SOCKETS IN THE MOLAR REGION: 1-YEAR POST-LOADING OUTCOME OF A RANDOMISED CONTROLLED TRIAL**

EUR J ORAL IMPLANTOL, 2017;10(3):263-278

CRESPI R, CAPPARÈ P, GHERLONE E  
**COMPARISON OF MAGNESIUM-ENRICHED HYDROXYAPATITE AND PORCINE BONE IN HUMAN EXTRACTION SOCKET HEALING: A HISTOLOGIC AND HISTOMORPHOMETRIC EVALUATION**

INT J ORAL MAXILLOFAC IMPLANTS, 2011 SEP-OCT;26(5):1057-62

CORBELLA S, TASCHIERI S, FRANCETTI L, WEINSTEIN R, DEL FABBRO M

**HISTOMORPHOMETRIC RESULTS AFTER POSTEXTRACTION SOCKET HEALING WITH DIFFERENT BIOMATERIALS: A SYSTEMATIC REVIEW OF THE LITERATURE AND META-ANALYSIS**


INT J ORAL MAXILLOFAC IMPLANTS, 2017 SEP/OCT;32(5):1001-1017

RADUNOVIC M, DE COLLI M, DE MARCO P, DI NISIO C, FONTANA A, PIATTELLI A, CATALDI A, ZARA S  
**GRAPHENE OXIDE ENRICHMENT OF COLLAGEN MEMBRANES IMPROVES DPSCS DIFFERENTIATION AND CONTROLS INFLAMMATION OCCURRENCE**

J BIOMED MATER RES A, 2017 AUG;105(8):2312-2320



- KILINC A, ATAOL M  
**HOW EFFECTIVE IS COLLAGEN RESORBABLE MEMBRANE PLACEMENT AFTER PARTIALLY IMPACTED MANDIBULAR THIRD MOLAR SURGERY ON POSTOPERATIVE MORBIDITY? A PROSPECTIVE RANDOMIZED COMPARATIVE STUDY**  
BMC ORAL HEALTH, 2017 OCT 5;17(1):126
- TROIANO G, ZHURAKIVSKA K, LO MUZIO L, LAINO L, CICCÌU M, LO RUSSO L  
**COMBINATION OF BONE GRAFT AND RESORBABLE MEMBRANE FOR ALVEOLAR RIDGE PRESERVATION: A SYSTEMATIC REVIEW, META-ANALYSIS AND TRIAL SEQUENTIAL ANALYSIS**  
J PERIODONTOL, 2017 SEP 12;117:1-17. Epub ahead of print
- ROSSI R, LONGO E, MIJIRITSKY E  
**A NEW INTERPRETATION OF GUIDED IMPLANT SURGERY TO ACHIEVE AN OPTIMAL RESULT IN THE AESTHETIC ZONES**  
MEDICAL RESEARCH ARCHIVES, 2017 APRIL, VOL. 5, ISSUE 4
- BOLLE C, FELICE P, BARAUSSE C, PISTILLI V, TRULLENQUE-ERIKSSON A, ESPOSITO M  
**FOUR MM-LONG VERSUS LONGER IMPLANTS IN AUGMENTED BONE IN POSTERIOR ATROPHIC JAWS: 1-YEAR POST-LOADING RESULTS FROM A MULTICENTRE RANDOMISED CONTROLLED TRIAL**  
EUR J ORAL IMPLANTOL, 2018;11(1):31-47
- ESPOSITO M, DAVÒ R, MARTI PAGES C, FERRER FUENTES A, BARAUSSE C, PISTILLI R, IPPOLITO DR, FELICE P  
**IMMEDIATELY LOADED ZYGOMATIC IMPLANTS VS CONVENTIONAL DENTAL IMPLANTS IN AUGMENTED ATROPHIC MAXILLAE: 4 MONTHS POST-LOADING RESULTS FROM A MULTICENTRE RANDOMISED CONTROLLED TRIAL**  
EUR J ORAL IMPLANTOL, 2018;11(1):11-28
- GASTALDI G, FELICE P, PISTILLI V, BARAUSSE C, IPPOLITO DR, ESPOSITO M  
**POSTERIOR ATROPHIC JAWS REHABILITATED WITH PROSTHESES SUPPORTED BY 5 × 5 MM IMPLANTS WITH A NANOSTRUCTURED CALCIUM-INCORPORATED TITANIUM SURFACE OR BY LONGER IMPLANTS IN AUGMENTED BONE. 3-YEAR RESULTS FROM A RANDOMISED CONTROLLED TRIAL**  
EUR J ORAL IMPLANTOL, 2018;11(1):49-61
- DIOMEDE F, D'AURORA M, GUGLIANDOLO A, MERCIARO I, ORSINI T, GATTA V, PIATTELLI A, TRUBIANI O, MAZZON E  
**BIOFUNCTIONALIZED SCAFFOLD IN BONE TISSUE REPAIR**  
INT J OF MOLECULAR SCIENCES, 2018, 19, 1022
- SCARANO A, LORUSSO F, SANTOS DE OLIVEIRA P, MURMURA G, CARINCI F  
**DENTAL IMPLANTS DISPLACED INTO THE MANDIBULAR CORPUS: CLINICAL NOTE**  
JOURNAL OF DENTISTRY AND ORAL CARE, 2018, 4(1):8-12
- FORABOSCO A, GHENO E, SPINATO S, GARUTI G, FORABOSCO E, CONSOLO U  
**CONCENTRATED GROWTH FACTORS IN MAXILLARY SINUS FLOOR AUGMENTATION: A PRELIMINARY CLINICAL COMPARATIVE EVALUATION**  
INT J OF GROWTH FACTORS AND STEM CELLS IN DENTISTRY, 2018;1:2-7
- PISTILLI R, FELICE P, CANNIZZARO G, PIATTELLI M, CORVINO V, BARAUSSE C, BUTI J, SOARDI E, ESPOSITO M  
**POSTERIOR ATROPHIC JAWS REHABILITATED WITH PROSTHESES SUPPORTED BY 6 MM LONG, 4 MM WIDE IMPLANTS OR BY LONGER IMPLANTS IN AUGMENTED BONE. ONE-YEAR POST-LOADING RESULTS FROM A PILOT RANDOMISED CONTROLLED TRIAL**  
EUR J ORAL IMPLANTOL, 2013;6(4):359-372
- FELICE P, BARAUSSE C, PISTILLI V, PIATTELLI M, IPPOLITO DR, ESPOSITO M  
**POSTERIOR ATROPHIC JAWS REHABILITATED WITH PROSTHESES SUPPORTED BY 6 MM LONG X 4 MM WIDE IMPLANTS OR BY LONGER IMPLANTS IN AUGMENTED BONE. 3-YEAR POST-LOADING RESULTS FROM A PILOT RANDOMISED CONTROLLED TRIAL**  
EUR J ORAL IMPLANTOL, 2018;11(2):175-187
- DAVÒ R, FELICE P, PISTILLI R, BARAUSSE C, MARTI PAGES C, FERRER FUENTES A, IPPOLITO DR, ESPOSITO M  
**IMMEDIATELY LOADED ZYGOMATIC IMPLANTS VS CONVENTIONAL DENTAL IMPLANTS IN AUGMENTED ATROPHIC MAXILLAE: 1-YEAR POST-LOADING RESULTS FROM A MULTICENTRE RANDOMISED CONTROLLED TRIAL**  
EUR J ORAL IMPLANTOL, 2018;11(1):145-161
- OMORI Y, SILVA ER, BOTTICELLI D, APAZA ALCCAYHUAMAN KA, LANG NP, XAVIER SP  
**REPOSITION OF THE BONE PLATE OVER THE ANTROSTOMY IN MAXILLARY SINUS AUGMENTATION: A HISTOMORPHOMETRIC STUDY IN RABBITS**  
CLIN ORAL IMPLANTS RES, 2018 AUG;29(8):821-834
- CANULLO L, GENOVA T, NAENNIC N, NAKAJIMA Y, MASUDA K, MUSSANO F  
**PLASMA OF ARGON ENHANCES THE ADHESION OF MURINE OSTEOBLASTS ON DIFFERENT GRAFT MATERIALS**  
ANN ANAT, 2018 JUL;218:265-270
- RAMÍREZ FERNÁNDEZ MP, MAZÓN P, GEHRKE SA, CALVO GUIRADO JL, DE AZA PN  
**COMPARISON OF TWO XENOGRAFT MATERIALS USED IN SINUS LIFT PROCEDURES: MATERIAL CHARACTERIZATION AND IN VIVO BEHAVIOR**  
MATERIALS, 2017 JUN 7;10(6)
- DIOMEDE F, ZINI N, GATTA V, FULLE S, MERCIARO I, D'AURORA M, LA ROVERE RM, TRAINI T, PIZZICANNELLA J, BALLERINI P, PIATTELLI A, TRUBIANI O  
**HUMAN PERIODONTAL LIGAMENT STEM CELLS CULTURED ONTO CORTICO-CANCELLOUS SCAFFOLD DRIVE BONE REGENERATIVE PROCESS**  
EUR CELL MATER, 2016 SEP 16;32:181-201
- DEVELIOGLU H, OZCAN G, GULTEKIN SE, SENGUVEN B, YILDIRIM A  
**THE SHORT-TERM EFFECTS OF VARIOUS XENOGRAFTS ON BONE HEALING IN RATS CRANIAL DEFECTS**  
BIOMEDICAL RESEARCH, 2018;29(8): 1598-1602
- BRUNELLI G, CARINCI F, GIRARDI A, PALMIERI A, CACCIANIGA GL, SOLLAZZO V  
**OSTEOBIOL® EFFECT ON DENTAL PULP DERIVED STEM CELLS**  
EUR J INFLAMMATION, VOL. 10, NO. 1 (S), 27-30
- MAZZONI S, MOHAMMADI S, TROMBA G, DIOMEDE F, PIATTELLI A, TRUBIANI O, GIULIANI A  
**ROLE OF CORTICO-CANCELLOUS HETEROLOGOUS BONE IN HUMAN PERIODONTAL LIGAMENT STEM CELLS XENO-FREE CULTURE STUDIED BY SYNCHROTRON RADIATION PHASE-CONTRAST MICRO TOMOGRAPHY**  
INT J MOL SCI, 2017;18;364
- LAURITANO D, CARINCI F, ZOLLINO I, HASSANIPOUR A, SAGGESE V, PALMIERI A, GIRARDI A, CURA F, PIRAS A, ZAMBONI P, BRUNELLI G  
**OSTEOBIOL® ENHANCES OSTEOGENIC DIFFERENTIATION IN BONE MARROW DERIVED STEM CELLS**  
EUR J INFLAMMATION, VOL. 10, NO. 1 (S3), 83-88
- MATÉ SANCHEZ DE VAL J, MAZON P, PIATTELLI A, CALVO GUIRADO J, MAREQUE BUENO J, GRANERO MARIN J, DE AZA P  
**COMPARISON AMONG THE PHYSICAL PROPERTIES OF CALCIUM PHOSPHATE-BASED BONE SUBSTITUTES OF NATURAL OR SYNTHETIC ORIGIN**  
INT J APPLIED CERAMIC TECHNOL, 2018;15(4):930-937
- NAKAJIMA Y, PIATTELLI A, IEZZI G, FORTICH MESA N, FERRI M, BOTTICELLI D  
**INFLUENCE OF THE PRESENCE OF ALVEOLAR MUCOSA AT IMPLANTS: A HISTOLOGICAL STUDY IN HUMANS**  
IMPLANT DENT, 2018;27(2):193-201
- CASCONI P, VELLONE V, RAMIERI V, BASILE E, TARSITANO A, MARCHETTI C  
**RECONSTRUCTION OF THE ADULT HEMIFACIAL MICROSMIA PATIENT WITH TEMPOROMANDIBULAR JOINT TOTAL JOINT PROSTHESIS AND ORTHOGNATIC SURGERY**  
CASE REP SURG, 2018;1-6
- AMR AEH, ABDEL GHAFAR KA, ABUEL-ELA HA, ABD ELHAMID ES  
**XENOGENIC FLEXIBLE BONE LAMINA GRAFT: A SUCCESSFUL ALTERNATIVE TO THE AUTOGENOUS ONLY BONE BLOCK GRAFT IN ALVEOLAR RIDGE AUGMENTATION: A CLINICAL, RADIOGRAPHIC AND HISTOLOGICAL EVALUATION**  
JOURNAL OF DENTAL TREATMENT AND ORAL CARE, 2017;1(1):104
- ATTIA AM  
**CLINICAL AND RADIOGRAPHIC EVALUATION OF CORTICO-CANCELLOUS BONE MIX XENOGRAFT (OSTEOBIOL GEN-OS) IN THE TREATMENT OF HUMAN PERIODONTAL INTRABONY DEFECTS**  
EGYPTIAN DENTAL JOURNAL, 2017;63;1-10
- BECHARA S, KUBILIUS R, VERONESI G, PIRES JT, SHIBLI JA, MANGANO FG  
**SHORT (6-MM) DENTAL IMPLANTS VERSUS SINUS FLOOR ELEVATION AND PLACEMENT OF LONGER (>10-MM) DENTAL IMPLANTS: A RANDOMIZED CONTROLLED TRIAL WITH A 3-YEAR FOLLOW-UP**  
CLIN ORAL IMPLANTS RES, 2017;28;1097-1107
- CHANDRASEKARAN B, SURESH N, MUTHUSAMY S  
**PLATELET-RICH FIBRIN WITH BONE GRAFTS FOR REGENERATION OF BONY DEFECT FOLLOWING EXTRACTION OF SUPERNUMERARY TEETH: A CASE REPORT**  
CHIN J DENT RES, 2017;20(4):231-234
- CHIRILA L, ROTARU C, FILIPOV I, SANDULESCU M  
**MANAGEMENT OF ACUTE MAXILLARY SINUSITIS AFTER SINUS BONE GRAFTING PROCEDURES WITH SIMULTANEOUS DENTAL IMPLANTS PLACEMENT - A RETROSPECTIVE STUDY**  
BMC INFECT DIS, 2016 MAR 8;16 SUPPL 1:94
- EKSTEIN J, TANDELICH M, NART J, CALVO GUIRADO JL, SHAPIRA L  
**MARGINAL BONE LEVEL AROUND CONICAL CONNECTION TAPERED IMPLANTS WITH PLATFORM SWITCHING: A MULTICENTRE RETROSPECTIVE STUDY AT 14 MONTHS FOLLOW-UP**  
J OSSEOINTEGRATION, 2016 JAN-APR;8(1)
- DEL CORSO M, DOHAN EHRENFEST DM  
**IMMEDIATE IMPLANTATION AND PERI-IMPLANT NATURAL BONE REGENERATION (NBR) IN THE SEVERELY RESORBED POSTERIOR MANDIBLE USING LEUKOCYTE- AND PLATELET-RICH FIBRIN (L-PRF): A 4-YEAR FOLLOW-UP**  
POSEIDO, 2013;1(2):85-92
- NOAMI SA, ELMOSY K, ASKAR N  
**EVALUATION OF PRE-HYDRATED COLLAGENATED CORTICO-CANCELLOUS GRANULES (MP3®) IN AUGMENTATION OF THE MAXILLARY SINUS (PRELIMINARY STUDY)**  
J OF DENT, ORAL DISORD AND THERAPY, 2014;2(3): 1-8
- BARONE A, TOTI P, PIATTELLI A, IEZZI G, DERCHI G, COVANI U  
**EXTRACTION SOCKET HEALING IN HUMANS AFTER RIDGE PRESERVATION TECHNIQUES: COMPARISON BETWEEN FLAPLESS AND FLAPPED PROCEDURES IN A RANDOMIZED CLINICAL TRIAL**  
J PERIODONTOL, 2014 JAN;85(1):14-23
- LOPEZ MA, LICO S, CASALE M, ORMANIER Z, CARINCI F  
**THE USE OF VARIOUS BIOMATERIALS IN COMPUTER-GUIDED CRESTAL SINUS LIFT PROCEDURES. A REPORT ON TWO CASE STUDIES WITH VOLUME COMPARISON**  
ORAL IMPLANTOL, 2016;APR-JUN 9(2)
- MEHL C, BOSCH T  
**IMPLANTOLOGY AND PROTHODONTICS AT CROSSROADS - CLASSIC VS MODERN TREATMENT CONCEPTS**  
JSM DENTAL SURGERY, 2016;1(1):1003
- NEMTOI A, DANILA V, DRAGAN E, PASCA S, NEMTOI A, CONSTANTIN M, SAVA A, HABA D  
**THE EFFECTS OF INSULIN AND STRONTIUM RANELATE ON GUIDED BONE REGENERATION IN DIABETIC RATS**  
REVISTA DE CHIMIE, 2017;68(4):693-697
- RAMIREZ FERNANDEZ MP, GEHRKE SA, ALBACETE MARTINEZ CP, CALVO GUIRADO JL, DE AZA PN  
**SEM-EDX STUDY OF THE DEGRADATION PROCESS OF TWO XENOGRAFT MATERIALS USED IN SINUS LIFT PROCEDURES**  
MATERIALS, 2017;10;542
- ZITA GOMES R, DE VASCONCELOS MR, LOPES GUERRA IM, DE ALMEIDA RAB, DE CAMPOS FELINO AC  
**IMPLANT STABILITY IN THE POSTERIOR MAXILLA: A CONTROLLED CLINICAL TRIAL**  
BIOMED RES INT, 2017;6825213
- COVANI U, CANULLO L, TOTI P, ALFONSI F, BARONE A  
**TISSUE STABILITY OF IMPLANTS PLACED IN FRESH EXTRACTION SOCKETS: A 5-YEAR PROSPECTIVE SINGLE-COHORT STUDY**  
J PERIODONTOL, 2014 SEP;85(9):E323-32
- DIOMEDE F, D'AURORA M, GUGLIANDOLO A, MERCIARO I, ETORRE V, BRAMANTI A, PIATTELLI A, GATTA V, MAZZON E, FONTANA A, TRUBIANI O  
**A NOVEL ROLE IN SKELETAL SEGMENT REGENERATION OF EXTRACELLULAR VESICLES RELEASED FROM PERIODONTAL-LIGAMENT STEM CELLS**  
INT J NANOMEDICINE, 2018 JUN 29;13:3805-3825



**Histology showing a granule of OsteoBiol® mp3® partially resorbed and in contact with new bone tissue. Htx-eosine**  
Author: Prof Ulf Nanmark, University of Göteborg, Sweden

KAWAKAMI S, LANG NP, IIDA T, FERRI M, APAZA ALCCAYHUAMAN KA, BOTTICELLI D  
**INFLUENCE OF THE POSITION OF THE ANTROSTOMY IN SINUS FLOOR ELEVATION ASSESSED WITH CONE-BEAM COMPUTED TOMOGRAPHY: A RANDOMIZED CLINICAL TRIAL**  
JOURNAL INVESTIG CLIN DENT, 2018 NOV;9(4):E12362

SCARANO A, DE OLIVEIRA PS, TRAINI T, LORUSSO F  
**SINUS MEMBRANE ELEVATION WITH HETEROLOGOUS CORTICAL LAMINA: A RANDOMIZED STUDY OF A NEW SURGICAL TECHNIQUE FOR MAXILLARY SINUS FLOOR AUGMENTATION WITHOUT BONE GRAFT**  
MATERIALS, 2018 AUG 17;11(8)

FIGLIUZZI MM, GIUDICE A, CRISTOFARO MG, PACIFICO D, BIAMONTE P, FORTUNATO L  
**POSTEXTRACTIVE IMPLANTS IN AESTHETIC AREAS: EVALUATION OF PERIMPLANT BONE REMODELING OVER TIME**  
ANN STOMATOL, 2015 MAY 18;6(1):29-34

KIVOVICS M, SZABÓ BT, NÉMETH O, TARI N, DORI F, NAGY P, DOBÓ-NAGY C, SZABÓ G  
**MICROARCHITECTURAL STUDY OF THE AUGMENTED BONE FOLLOWING RIDGE PRESERVATION WITH A PORCINE XENOGRAFT AND A COLLAGEN MEMBRANE: PRELIMINARY REPORT OF A PROSPECTIVE CLINICAL, HISTOLOGICAL, AND MICRO-COMPUTED TOMOGRAPHY ANALYSIS**  
INT J ORAL MAXILLOFAC SURG, 2017 FEB;46(2):250-260

MARCONCINI S, GIAMMARINARO E, DERCHI G, ALFONSI F, COVANI U, BARONE A  
**CLINICAL OUTCOMES OF IMPLANTS PLACED IN RIDGE-PRESERVED VERSUS NONPRESERVED SITES: A 4-YEAR RANDOMIZED CLINICAL TRIAL**  
CLIN IMPL DENT RELAT RES, 2018 DEC;20(6):906-914

MARCONCINI S, COVANI U, GIAMMARINARO E, VELASCO-ORTEGA E, DE SANTIS D, ALFONSI F, BARONE A  
**CLINICAL SUCCESS OF DENTAL IMPLANTS PLACED IN POSTERIOR MANDIBLE AUGMENTED WITH INTERPOSITIONAL BLOCK GRAFT: 3-YEAR RESULTS FROM A PROSPECTIVE COHORT CLINICAL STUDY**  
J ORAL MAXILLOFAC SURG. 2019 FEB;77(2):289-298

IIDA T, RICARDO SILVA E, LANG NP, APAZA ALCCAYHUAMAN KA, DANIELE BOTTICELLI D, XAVIER SP  
**HISTOLOGICAL AND MICRO-COMPUTED TOMOGRAPHY EVALUATIONS OF NEWLY FORMED BONE AFTER MAXILLARY SINUS AUGMENTATION USING A XENOGRAFT WITH SIMILAR DENSITY AND MINERAL CONTENT OF BONE: AN EXPERIMENTAL STUDY IN RABBITS**  
CLIN EXP DENT RES, 2018;1-7

BERNARDI S, GATTO R, SEVERINO M, BOTTICELLI G, CARUSO S, RASTELLI C, LUPI E, ROIAS AQ, IACOMINO E, FALISI G  
**SHORT VERSUS LONGER IMPLANTS IN MANDIBULAR ALVEOLAR RIDGE AUGMENTED USING OSTEOGENIC DISTRACTION: ONE-YEAR FOLLOW-UP OF A RANDOMIZED SPLIT-MOUTH TRIAL**  
J ORAL IMPLANTOL, 2018 JUN;44(3):184-191

DIKER N, SARICA H, CUMBUL A, KILIC E  
**EFFECTS OF SYSTEMIC ERYTHROPOIETIN TREATMENT AND HETEROGENEOUS XENOGRAFT IN COMBINATION ON BONE REGENERATION OF A CRITICAL-SIZE DEFECT IN AN EXPERIMENTAL MODEL**  
J CRANIOMAXILLOFAC SURG, 2018 NOV;46(11):1919-1923

SCARANO A, MURMURA G, MASTRANGELO F, LORUSSO F, GRECO LUCCHINA A, CARINCI F  
**A NOVEL TECHNIQUE TO PREVENT SINUS MEMBRANE COLLAPSE DURING MAXILLARY SINUS FLOOR AUGMENTATION WITHOUT BONE GRAFT: TECHNICAL NOTE**  
J BIOL REGUL HOMEOST AGENTS, 2018 NOV-DEC;32(6):1589-1592

SENESE O, BOUTREMAN E, GOSSIAUX C, LOEB I, DEQUANTER D  
**RETROSPECTIVE ANALYSIS OF 79 PATIENTS WITH ORBITAL FLOOR FRACTURE: OUTCOMES AND PATIENT-REPORTED SATISFACTION**  
ARCH CRANIOFAC SURG, 2018 JUN;19(2):108-113

GENOVA T, PESCE P, MUSSANO F, TANAKA K, CANULLO L  
**THE INFLUENCE OF BONE-GRAFT BIO-FUNCTIONALIZATION WITH PLASMA OF ARGON ON BACTERIAL CONTAMINATION**  
J BIOMED MATER RES A, 2019 JAN;107(1):67-70

ESPOSITO M, BARAUSSE C, PISTILLI R, PIATTELLI M, DI SIMONE S, IPPOLITO DR, FELICE P  
**POSTERIOR ATROPHIC JAWS REHABILITATED WITH PROSTHESES SUPPORTED BY 5 X 5 MM IMPLANTS WITH A NOVEL NANOSTRUCTURED CALCIUM-INCORPORATED TITANIUM SURFACE OR BY LONGER IMPLANTS IN AUGMENTED BONE. FIVE-YEAR RESULTS FROM A RANDOMISED CONTROLLED TRIAL**  
INT J OF ORAL IMPLANTOL, 2019;12(1):39-54

FELICE P, PISTILLI R, BARAUSSE C, PIATTELLI M, BUTI J, ESPOSITO M  
**POSTERIOR ATROPHIC JAWS REHABILITATED WITH PROSTHESES SUPPORTED BY 6 MM LONG X 4 MM WIDE IMPLANTS OR BY LONGER IMPLANTS IN AUGMENTED BONE. FIVE-YEAR POST-LOADING RESULTS FROM A WITHIN-PERSON RANDOMISED CONTROLLED TRIAL**  
INT J OF ORAL IMPLANTOL, 2019;12(1):57-72

DI CARLO R, ZARA S, VENTRELLA A, SIANI G, DA ROS T, IEZZI G, CATALDI A, FONTANA A  
**COVALENT DECORATION OF CORTICAL MEMBRANES WITH GRAPHENE OXIDE AS A SUBSTRATE FOR DENTAL PULP STEM CELLS**  
NANOMATERIALS, 2019;9:604

CABALLÉ-SERRANO J, MUNAR-FRAUA A, DELGADO L, PÉREZ R, HERNÁNDEZ-ALFARO F  
**PHYSICO-CHEMICAL CHARACTERIZATION OF BARRIER MEMBRANES FOR BONE REGENERATION**  
J MECH BEHAV BIOMED, 2019;97:13-20

MATOH U, PETELIN M, GASPERIC R  
**SPLIT-MOUTH COMPARISON OF CORONALLY ADVANCED FLAP WITH CONNECTIVE TISSUE GRAFT OR COLLAGEN MATRIX FOR TREATMENT OF ISOLATED GINGIVAL RECESSIONS**  
INT J PERIODONTICS RESTORATIVE DENT, 2019;39(3):439-446

AMBROZEWICZ E, MUSZYNSKA M, TOKAJUK G, GRYNKIEWICZ G, ŻARKOVIC N, SKRZYDLEWSKA  
**BENEFICIAL EFFECTS OF VITAMINS K AND D3 ON REDOX BALANCE OF HUMAN OSTEOBLASTS CULTURED WITH HYDROXYAPATITE-BASED BIOMATERIALS CELLS**  
2019;8:325

POLIS-YANES C, CADENAS-SEBASTIÁN C, GUAL-VAQUÉS P, AYOUSO-MONTERO R, MARÍ-ROIG A, LÓPEZ-LÓPEZ J  
**GUIDED BONE REGENERATION OF AN ATROPHIC MAXILLA USING HETEROLOGOUS CORTICAL LAMINA**  
CASE REP DENT., 2019; 5216362

KAWAKAMI S, LANG NP, FERRI M, APAZA ALCCAYHUAMAN KA, BOTTICELLI D  
**INFLUENCE OF THE HEIGHT OF THE ANTROSTOMY IN SINUS FLOOR ELEVATION ASSESSED BY CONE BEAM COMPUTED TOMOGRAPHY: A RANDOMIZED CLINICAL TRIAL**  
INT J ORAL MAXILLOFAC IMPLANTS, 2019;34(1):223-232

CHECCHI V, GASPARRO R, PISTILLI R, CANULLO L, FELICE P  
**CLINICAL CLASSIFICATION OF BONE AUGMENTATION PROCEDURE FAILURES IN THE ATROPHIC ANTERIOR MAXILLAE: ESTHETIC CONSEQUENCES AND TREATMENT OPTIONS**  
BIOMED RES INT., 2019:4386709

FISCHER K R, TESTORI T, WACHTEL H, MÜHLEMANN S, HAPPE A, DEL FABBRO M  
**SOFT TISSUE AUGMENTATION APPLYING A COLLAGENATED PORCINE DERMAL MATRIX DURING SECOND STAGE SURGERY: A PROSPECTIVE MULTICENTER CASE SERIES**  
CLIN IMPLANT DENT RELAT RES.,2019;1-8

HIROTA A, LANG NP, FERRI M, FORTICH MESA N, APAZA ALCCAYHUAMAN KA, BOTTICELLI D  
**TOMOGRAPHIC EVALUATION OF THE INFLUENCE OF THE PLACEMENT OF A COLLAGEN MEMBRANE SUBJACENT TO THE SINUS MUCOSA DURING MAXILLARY SINUS FLOOR AUGMENTATION: A RANDOMIZED CLINICAL TRIAL**  
INT J IMPLANT DENT.,2019; 5(1):31

TANAKA K, IEZZI G, PIATTELLI A, FERRI M, MESA NF, APAZA ALCCAYHUAMAN KA, BOTTICELLI D  
**SINUS FLOOR ELEVATION AND ANTROSTOMY HEALING: A HISTOMORPHOMETRIC CLINICAL STUDY IN HUMANS**  
IMPLANT DENT., 2019; AUG 14

JEANNEAU C, LE FOURNIS C, ABOUT I  
**XENOGENIC BONE FILLING MATERIALS MODULATE MESENCHYMAL STEM CELL RECRUITMENT: ROLE OF THE COMPLEMENT CSA**  
CLIN. ORAL INVESTIG.; 2019 OCT 23

RAMANAUSKAITE A, BORGES T, ALMEIDA BL, CORREIA A  
**DENTAL IMPLANT OUTCOMES IN GRAFTED SOCKETS: A SYSTEMATIC REVIEW AND META-ANALYSIS**  
J ORAL MAXILLOFAC RES., 2019 SEP 5;10(3):E8

FARIA-ALMEIDA R, ASTRAMSKAITE-JANUSEVICIENE I, PUISYS A, CORREIA F  
**EXTRACTION SOCKET PRESERVATION WITH OR WITHOUT MEMBRANES, SOFT TISSUE INFLUENCE ON POST EXTRACTION ALVEOLAR RIDGE PRESERVATION: A SYSTEMATIC REVIEW**  
J ORAL MAXILLOFAC RES., 2019 SEP 5;10(3):E5

ADIOGLU S, GIRAY CB, KULAC I, USUBUTUN A, AKTAS A  
**CLINICAL AND HISTOPATHOLOGICAL COMPARATIVE STUDY OF TWO EQUINE-DERIVED BONE GRAFT: A HUMAN STUDY**  
J PAK MED ASSOC., 2019; 69: 1617

ROSSI R, FOCE E  
**RECONSTRUCTION OF A HORIZONTAL AND VERTICAL BONE DEFECT USING THE CORTICAL LAMINA TECHNIQUE**  
MED. RES. ARCH, 2019; VOL 7 (11)

CANULLO L, GENOVA T, RAKIC M, SCULEAN A, MIRON R, MUZZI M, CAROSSA S, MUSSANO F  
**EFFECTS OF ARGON PLASMA TREATMENT ON THE OSTEOCONDUCTIVITY OF BONE GRAFTING MATERIALS**  
CLIN ORAL INVESTIG, 2019 NOV 20




# OsteoBiol® product codes



PRODUCT	PACKAGING	TYPE	SIZE	PORCINE CODE	EQUINE CODE
<b>BONE SUBSTITUTES</b>					
Gen-Os®	1 Vial	DRIED GRANULES	0.25 g	M1052FS	M1052FE
Gen-Os®	1 Vial	DRIED GRANULES	0.5 g	M1005FS	M1005FE
Gen-Os®	1 Vial	DRIED GRANULES	1.0 g	M1010FS	M1010FE
Gen-Os®	1 Vial	DRIED GRANULES	2.0 g	M1020FS	M1020FE
Gen-Os® 1000-2000	1 Vial	DRIED GRANULES	1.0 g	M0210FS	
Gen-Os® 1000-2000	1 Vial	DRIED GRANULES	2.0 g	M0220FS	
TSV Gel	1 Syringe	TSV GEL	0.5 g	TSV005S <small>in kit with M1005FS or A1005FS</small>	TSV005E <small>in kit with M1005FE or A1005FE</small>
TSV Gel	1 Syringe	TSV GEL	1.0 g	TSV010S <small>in kit with M1010FS or A1010FS</small>	TSV010E <small>in kit with M1010FE or A1010FE</small>
mp3®	1 Syringe	BONE MIX	0.5 cc	A3095FS	A3095FE
mp3®	1 Syringe	BONE MIX	1.0 cc	A3005FS	A3005FE
mp3®	3 Syringes	BONE MIX	3x0.25 cc (0.75 cc)	A3075FS	
mp3®	3 Syringes	BONE MIX	3x0.5 cc (1.5 cc)	A3015FS	A3015FE
mp3®	3 Syringes	BONE MIX	3x1.0 cc (3.0 cc)	A3030FS	A3030FE
mp3®	1 Syringe (wide tip)	BONE MIX	2.0 cc	A3010FS	A3010FE
mp3® 1000-2000	1 Syringe (wide tip)	BONE MIX	2.0 cc	A3210FS	A3210FE
GTO®	1 Syringe	BONE MIX + TSV Gel	0.5 cc	MU0005S	MU0005E
GTO®	1 Syringe	BONE MIX + TSV Gel	2.0 cc	MU0020S	MU0020E
Putty	1 Syringe	BONE PASTE	0.25 cc	HPT52S	
Putty	1 Syringe	BONE PASTE	0.5 cc	HPT09S	HPT09E
Putty	3 Syringes	BONE PASTE	3x0.25 cc (0.75 cc)	HPT32S	HPT32E
Putty	3 Syringes	BONE PASTE	3x0.5 cc (1.5 cc)	HPT35S	HPT35E
Putty	1 Syringe (wide tip)	BONE PASTE	1.0 cc	HPT61S	HPT61E
Gel 40	1 Syringe	BONE GEL	0.5 cc	05GEL40S	05GEL40E
Gel 40	3 Syringes	BONE GEL	3x0.5 cc (1.5 cc)	15GEL40S	15GEL40E
Apatos Mix	1 Vial	DRIED GRANULES	0.5 g	A1005FS	A1005FE
Apatos Mix	1 Vial	DRIED GRANULES	1.0 g	A1010FS	A1010FE
Apatos Mix	1 Vial	DRIED GRANULES	2.0 g	A1020FS	A1020FE
Apatos Cortical	1 Vial	DRIED GRANULES	0.5 g	AC1005FS	
Apatos Cortical	1 Vial	DRIED GRANULES	1.0 g	AC1010FS	
Apatos Mix 1000-2000	1 Vial	DRIED GRANULES	1.0 g	A0210FS	A0210FE




PRODUCT	PACKAGING	TYPE	SIZE	PORCINE CODE	EQUINE CODE
<b>BLOCKS</b>					
Sp-Block	1 Blister	DRIED BLOCK / NORM	10x10x10 mm		BN0E
Sp-Block	1 Blister	DRIED BLOCK / NORM	10x10x20 mm		BN1E
Sp-Block	1 Blister	DRIED BLOCK / NORM	10x20x20 mm		BN2E
Sp-Block	1 Blister	DRIED BLOCK / NORM	35x10x5 mm		BN8E
Dual-Block CURVED	1 Blister	DRIED BLOCK / SOFT	20x15x5 mm	STS7S	
Dual-Block CURVED	1 Blister	DRIED BLOCK / NORM	20x10x5 mm	STN5S	
<b>MEMBRANES AND BARRIERS</b>					
Evolution	3 Blister	DRIED / X-FINE	30x30x (0.2) mm	EM33XS	
Evolution	1 Blister	DRIED / FINE	20x20x (0.3) mm		EV02LLE
Evolution	1 Blister	DRIED / FINE	30x30x (0.3) mm		EV03LLE
Evolution	1 Blister	DRIED / FINE	Oval 25x35x (0.3) mm		EVOLLE
Evolution	1 Blister	DRIED / FINE	40x40x (0.3) mm		EV04LLE
Evolution	1 Blister	DRIED / FINE	80x60x (0.3) mm		EV06LLE
Evolution	1 Blister	DRIED / STANDARD	20x20x (0.4) mm	EM02HS	EV02HHE
Evolution	1 Blister	DRIED / STANDARD	30x30x (0.4) mm	EM03HS	EV03HHE
Evolution	1 Blister	DRIED / STANDARD	Oval 25x35x (0.4) mm	EM00HS	
Derma	1 Blister	DRIED	20x20x (0.5) mm	ED02LS	
Derma	1 Blister	DRIED	Oval 12x8x (0.9) mm	ED21FS	
Derma	1 Blister	DRIED	25x25x (0.9) mm	ED25FS	
Derma	1 Blister	DRIED	50x50x (0.9) mm	ED05FS	
Derma	1 Blister	DRIED	7x50x (2.0) mm	ED75SS	
Derma	1 Blister	DRIED	15x5x (2.0) mm	ED15SS	
Derma	1 Blister	DRIED	30x30x (2.0) mm	ED03SS	
Derma	1 Blister	DRIED	50x50x (2.0) mm	ED05SS	
Soft Cortical Lamina	1 Blister	DRIED	25x25x (0.5) mm	LS25FS	LS25FE
Soft Cortical Lamina	1 Blister	DRIED	Oval 25x35x (0.5) mm	LS23FS	LS23FE
Soft Cortical Lamina	1 Blister	DRIED	20x40x (1.0) mm	LS24LS	
Curved Lamina	1 Blister	DRIED	35x35x (1.0) mm	LS10HS	LS10HE
Soft Cortical Lamina	1 Blister	DRIED	30x30x (3.0) mm	LS03SS	LS03SE
Cortical Lamina	1 Blister	DRIED	35x15x (0.7) mm	LS15LS	
Cortical Lamina	1 Blister	DRIED	35x35x (1.0) mm	LS35LS	
Special	1 Blister	DRIED	20x20x (0.2) mm		EM02LE
Special	1 Blister	DRIED	30x30x (0.2) mm		EM03LE







# OsteoBiol<sup>®</sup> by Tecross




Tecross s.r.l. is an innovative, globally active company that develops, produces and documents premium-quality xenogenic biomaterials by the brands Tecross<sup>®</sup> and OsteoBiol<sup>®</sup>.

Its 20 years of research led to its patent-protected production process that ensures neutralization of antigenic components in order to achieve biocompatibility, while preserving the natural collagen matrix inside the biomaterial.

Tecross<sup>®</sup> products comply with highest quality standards such as ISO 13485 and European laws.

[osteobiol.com](http://osteobiol.com)

DISTRIBUTED BY



MK-CAT0120EN

**Tecross<sup>®</sup> s.r.l.**  
info@tecross.com  
www.tecross.com

BIOMATERIALS ENGINEERING

**Tecross<sup>®</sup> Dental s.r.l.**  
info@tecross-dental.com  
www.osteobiol.com

INTERNATIONAL SALES & MARKETING