



OsteoBiol[®]
by TecnoSS

Sp-Block | Dual-Block

HIGHLY OSTEOCONDUCTIVE SCAFFOLDS

REGENERATION SCIENCE

INSPIRED BY NATURE



A unique biotechnology

TECNOSS®: A UNIQUE PROCESS THAT ACCELERATES AND GUIDES NATURAL BONE REGENERATION

Tecnoss® developed and patented a unique biotechnology that prevents the ceramization phase of natural bone and preserves the tissue collagen, allowing an osteoclastic-type remodelling of the biomaterial similar to physiological bone turnover and delivering a product endowed with characteristics very similar to human mineral bone⁽¹⁾.

The combination of these factors allows a consistent new bone formation and a close contact between neo-formed bone and biomaterial (Fig. A).

COLLAGEN: A KEY FACTOR FOR BONE REGENERATION

Collagen has a key role in bone regeneration process in that:

- it acts as a valid substrate for platelet activation and aggregation
- it serves to attract and differentiate the mesenchymal stem cells present in the bone marrow⁽²⁾
- it increases the proliferation rate of the osteoblasts up to 2/3 times⁽³⁾
- it stimulates the activation of the platelets, osteoblasts and osteoclasts in the tissue healing process

OSTEOBIOL®: UNIQUE COLLAGENATED BIOMATERIALS

Thanks to the innovative Tecness® technology, the OsteoBiol® line has the following important characteristics:

- absence of a foreign body response⁽⁴⁾
- gradual resorption over time^(5,6)
- stimulation/acceleration of physiological tissue healing process⁽²⁾
- protection of the grafting site from infection (membranes)⁽⁷⁾
- capability of carrying medication to the surgical site⁽⁸⁾

The Tecness® new generation of biomaterials, thanks to a revolutionary technology, goes beyond the simple role of aiding natural bone regrowth by stimulating and accelerating this vital physiological process.

OsteoBiol®
by Tecness



Fig. A




Fig. B




Fig. C




Fig. A – Histology at 3 months. Human mandible grafted with OsteoBiol® Sp-Block. Biopsy courtesy of Dr P Felice, Bologna, Italy. Histology courtesy of Prof U Nannmark, University of Göteborg, Sweden

Fig. B – SEM image of an OsteoBiol® bone matrix colonized by osteoblasts from a cell-line (MG63). Courtesy of Prof U Nannmark, University of Göteborg, Sweden

Fig. C – SEM image of an OsteoBiol® cancellous block trabeculae

(1) Figueiredo M et al. J Biomed Mater Res B Appl Biomater, 2010 Feb; 92(2):409-419

(2) Brunelli G et al. Eur J Inflamm, 2011, Vol. 9, no. 3 (S), 103-107

(3) Hsu FY et al. Biomaterials, 1999, 20:1931-1936

(4) Crespi R et al. Int J Oral Maxillofac Implants, 2011 Jul - Aug; 26(4):866-72

(5) Nannmark U, Sennerby L. Clin Implant Dent Relat Res, 2008 Dec; 10(4):264-70


(6) Barone A et al. Clin Implant Dent Relat Res, 2012 Jun; 14(3):373-9

(7) Barone A et al. Clin Oral Implants Res, 2013 Nov; 24(11):1231-7

(8) Fischer K et al. Clin Oral Implants Res, 2014 Sep 15 Epub ahead of print



A highly osteoconductive scaffold



CHARACTERISTICS


Sp-Block and *Dual-Block* support new bone formation^(9,10) because of their extremely osteoconductive surface: thanks to the rigid consistency these blocks are able to maintain in time the original graft volume, which is particularly important in case of large regenerations. Moreover the collagen content facilitates blood clotting and the subsequent invasion of regenerative and repairing cells, favoring *restitutio ad integrum* of missing bone.

HANDLING

Sp-Block and *Dual-Block* must be hydrated before use for with sterile lukewarm physiological solution or with antibiotics (5/10 minutes for *Sp-Block* and *Dual-Block Soft*, up to 40 minutes for *Dual-Block Norm*). Afterwards, the block can be adapted to the receiving site which must be accurately decorticated in order to guarantee maximum contact; the blocks should be always fixed with osteosynthesis microscrews. In case of vertical augmentation with inlay technique *Sp-Block* should be fixed also with miniplates. Protection with *OsteoBiol® Evolution* membrane is recommended.

CLINICAL INDICATIONS OVERVIEW

Sp-Block is indicated in cases where a vertical gain in posterior mandible is required^(11,12,13), to achieve an augmentation of maximum 5 mm, by means of the inlay technique. *Dual-Block* can be grafted with the onlay technique only to augment horizontally heavily resorbed maxilla. Whatever is the applied technique, it is recommended to fill the gaps around the block with a biomaterial in granules to achieve the desired volume and contour of the augmented recipient site.



(9) Scarano A et al. Clin Oral Implants Res, 2011 Oct; 22(10):1125-30
(10) Felice P et al. Int J Periodontics Restorative Dent, 2013 Mar-Apr;33(2):159-66
(11) Esposito M et al. Eur J Oral Implantol, 2012 Spring;5(1):19-33
(12) Felice P et al. Eur J Oral Implantol, 2012 Summer;5(2):149-61
(13) Pistilli R et al. Eur J Oral Implantol, 2013 Winter;6(4):343-357

Tissue of origin

Sp-Block: Cancellous bone
Dual-Block: Cortico-cancellous bone

Tissue collagen

Preserved

Physical form

Rigid dried block

Re-entry time

About 8 months, variable depending on characteristics and irrigation grade of grafting site and on clinical conditions of patient

Packaging

Sterile blister


Product codes

Sp-Block
BNOE | 10x10x10 mm | Equine
BN1E | 10x10x20 mm | Equine
BN2E | 10x20x20 mm | Equine
BN8E | 35x10x5 mm | Equine

Dual-Block
STS7S | 20x15x5 mm | Soft | Porcine curved
STN5S | 20x10x5 mm | Norm | Porcine curved

GMDN code

38746





Excellent clinical performances




Fig. 1



Fig. 2




Fig. 3




Fig. 4




Fig. 5




Fig. 6



Fig. 7



Fig. 8



Fig. 9




Fig. 10




Fig. 12




Fig. 11

CASE REPORT

VERTICAL REGENERATION

Vertical regeneration with inlay technique

Sex: **Female** | Age: **60**

Fig. 1 Pre-operative x-ray

Fig. 2 After one horizontal and two vertical osteotomies, the bone fragment is moved towards the coronal direction

Fig. 3 Space obtained after moving the bone fragment

Fig. 4 Positioning of OsteoBiol® Sp-Block

Fig. 5 Post-surgery x-ray

Fig. 6 Clinical appearance of the graft during re-opening, after 3 months

Fig. 7 Preparation of implant sites

Fig. 8 Positioning of the implants

Fig. 9 Positioning of the implants

Fig. 10 Histology after 3 months*

Fig. 11 Histology detail*


Fig. 12 Histology detail*

Documentation provided by
Dr **Pietro Felice**
University of Bologna, Italy
E-mail: pietro.felice@unibo.it

Prof **Ulf Nannmark***
Göteborg University, Sweden

Block: **OsteoBiol® Sp-Block**

- Barone A, Crespi R, Aldini NN, Fini M, Giardino R, Covani U
MAXILLARY SINUS AUGMENTATION: HISTOLOGIC AND HISTOMORPHOMETRIC ANALYSIS
 INT J ORAL MAXILLOFAC IMPLANTS, 2005 JUL-AUG; 20(4):519-25
- Barone A, Ameri S, Covani U
IMMEDIATE POSTEXTRACTION IMPLANTS: TREATMENT OF RESIDUAL PERI-IMPLANT DEFECTS. A RETROSPECTIVE ANALYSIS
 EUR J IMPLANT PROSTHODONTICS, 2006,2: 99-106
- Barone A, Santini S, Sbordone L, Crespi R, Covani U
A CLINICAL STUDY OF THE OUTCOMES AND COMPLICATIONS ASSOCIATED WITH MAXILLARY SINUS AUGMENTATION
 INT J ORAL MAXILLOFAC IMPLANTS, 2006 JAN-FEB; 21(1):81-5
- Covani U, Barone A, Cornelini R, Crespi R
CLINICAL OUTCOME OF IMPLANTS PLACED IMMEDIATELY AFTER IMPLANT REMOVAL
 J PERIODONTOL, 2006 APR;77(4):722-7
- Orsini G, Scarano A, Piattelli M, Piccirilli M, Caputi S, Piattelli A
HISTOLOGIC AND ULTRASTRUCTURAL ANALYSIS OF REGENERATED BONE IN MAXILLARY SINUS AUGMENTATION USING A PORCINE BONE-DERIVED BIOMATERIAL
 J PERIODONTOL, 2006 DEC;77(12):1984-90
- Trubiani O, Scarano A, Orsini G, Di Iorio D, D'Arcangelo C, Piccirilli M, Sigismondo M, Caputi S
THE PERFORMANCE OF HUMAN PERIODONTAL LIGAMENT MESENCHYMAL STEM CELLS ON XENOGENIC BIOMATERIALS
 INT J IMMUNOPATHOL PHARMACOL, 2007 JAN-MAR; 20(1 SUPPL 1):87-91
- Barone A, Aldini NN, Fini M, Giardino R, Calvo Guirado JL, Covani U
XENOGRAFT VERSUS EXTRACTION ALONE FOR RIDGE PRESERVATION AFTER TOOTH REMOVAL: A CLINICAL AND HISTOMORPHOMETRIC STUDY
 J PERIODONTOL, 2008 AUG;79(8):1370-7
- Covani U, Cornelini R, Barone A
BUCCAL BONE AUGMENTATION AROUND IMMEDIATE IMPLANTS WITH AND WITHOUT FLAP ELEVATION: A MODIFIED APPROACH
 INT J ORAL MAXILLOFAC IMPLANTS, 2008 SEP-OCT; 23(5):841-6
- Cardaropoli D, Cardaropoli G
PRESERVATION OF THE POSTEXTRACTION ALVEOLAR RIDGE: A CLINICAL AND HISTOLOGIC STUDY
 INT J PERIODONTICS RESTORATIVE DENT, 2008 OCT; 28(5):469-77
- Nannmark U, Sennerby L
THE BONE TISSUE RESPONSES TO PREHYDRATED AND COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS: A STUDY IN RABBIT MAXILLARY DEFECTS
 CLIN IMPLANT DENT RELAT RES, 2008 DEC;10(4):264-70. EPUB 2008 JAN 30
- Covani U, Marconcini S, Crespi R, Barone A
IMMEDIATE IMPLANT PLACEMENT AFTER REMOVAL OF A FAILED IMPLANT: A CLINICAL AND HISTOLOGICAL CASE REPORT
 J ORAL IMPLANTOL, 2009; 35(4):189-95
- Calvo Guirado JL, Gomez Moreno G, Lopez Mari L, Ortiz Ruiz AJ, Guardia J
ATRAUMATIC MAXILLARY SINUS ELEVATION USING THREADED BONE DILATORS FOR IMMEDIATE IMPLANTS. A THREE-YEAR CLINICAL STUDY
 MED ORAL PATOL ORAL CIR BUCAL, 2010 MAR 1; 15(2):E366-70
- Figueiredo M, Henriques J, Martins G, Guerra F, Judas F, Figueiredo H
PHYSICO-CHEMICAL CHARACTERIZATION OF BIOMATERIALS COMMONLY USED IN DENTISTRY AS BONE SUBSTITUTES - COMPARISON WITH HUMAN BONE
 J BIOMED MATER RES B APPL BIOMATER, 2010 FEB; 92(2):409-19
- Crespi R, Capparè P, Gherlone E
DENTAL IMPLANTS PLACED IN EXTRACTION SITES GRAFTED WITH DIFFERENT BONE SUBSTITUTES: RADIOGRAPHIC EVALUATION AT 24 MONTHS
 J PERIODONTOL, 2009 OCT; 80(10):1616-1621
- Nannmark U, Azarmehr I
SHORT COMMUNICATION: COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS. A STUDY IN RABBIT MAXILLARY DEFECTS
 CLIN IMPLANT DENT RELAT RESE, 2010 JUN 1; 12(2):161-3. EPUB 2010 APR 9
- Scarano A, Piattelli A, Assenza B, Quaranta A, Perrotti V, Piattelli M, Iezzi G
PORCINE BONE USED IN SINUS AUGMENTATION PROCEDURES: A 5-YEAR RETROSPECTIVE CLINICAL EVALUATION
 J ORAL MAXILLOFAC SURG, 2010 AUG; 68(8):1869-73. EPUB 2010 MAY 10
- Pagliani L, Andersson P, Lanza M, Nappo A, Verocchi D, Volpe S, Sennerby L
A COLLAGENATED PORCINE BONE SUBSTITUTE FOR AUGMENTATION AT NEOS IMPLANT SITES: A PROSPECTIVE 1-YEAR MULTICENTER CASE SERIES STUDY WITH HISTOLOGY
 CLIN IMPLANT DENT RELAT RES, 2012 OCT;14(5):746-58. EPUB 2010 OCT 26
- Scarano A, Piattelli A, Perrotti V, Manzoni L, Iezzi G
MAXILLARY SINUS AUGMENTATION IN HUMANS USING CORTICAL PORCINE BONE: A HISTOLOGICAL AND HISTOMORPHOMETRIC EVALUATION AFTER 4 AND 6 MONTHS
 CLIN IMPLANT DENT RELAT RES, 2011 MAR; 13(1):13-18
- Barone A, Ricci M, Calvo Guirado JL, Covani U
BONE REMODELLING AFTER REGENERATIVE PROCEDURES AROUND IMPLANTS PLACED IN FRESH EXTRACTION SOCKETS: AN EXPERIMENTAL STUDY IN BEAGLE DOGS
 CLIN ORAL IMPLANTS RES, 2011 OCT;22(10):1131-7
- Scarano A, Carinci F, Assenza B, Piattelli M, Murmura G, Piattelli A
VERTICAL RIDGE AUGMENTATION OF ATROPHIC POSTERIOR MANDIBLE USING AN INLAY TECHNIQUE WITH A XENOGRAFT WITHOUT MINISCREWS AND MINIPLATES: CASE SERIES
 CLIN ORAL IMPLANTS RES, 2011 OCT;22(10):1125-30. EPUB 2011 JAN 20
- Ramirez Fernandez Mp, Calvo Guirado JL, Maté Sanchez De Val Je, Delgado Ruiz Ra, Negri B, Barona Dorado C
ULTRASTRUCTURAL STUDY BY BACKSCATTERED ELECTRON IMAGING AND ELEMENTAL MICROANALYSIS OF BONE-TO-BIOMATERIAL INTERFACE AND MINERAL DEGRADATION OF PORCINE XENOGRAFTS USED IN MAXILLARY SINUS FLOOR ELEVATION
 CLIN ORAL IMPLANTS RES, 2013 MAY;24(5):523-30. EPUB 2012 JAN 26
- Esposito M, Cannizaro G, Soardi E, Pistilli R, Piattelli M, Corvino V, Felice P
POSTERIOR ATROPHIC JAWS REHABILITATED WITH PROSTHESES SUPPORTED BY 6 MM-LONG, 4 MM-WIDE IMPLANTS OR BY LONGER IMPLANTS IN AUGMENTED BONE. PRELIMINARY RESULTS FROM A PILOT RANDOMISED CONTROLLED TRIAL
 EUR J ORAL IMPLANTOL, 2012 SPRING;5(1):19-33
- Felice P, Piana L, Checchi L, Pistilli R, Pellegrino G
VERTICAL RIDGE AUGMENTATION OF THE ATROPHIC POSTERIOR MANDIBLE WITH A 2-STAGE INLAY TECHNIQUE: A CASE REPORT
 IMPLANT DENT, 2012 JUN;21(3):190-5
- Bottini Lp, Ricci L, Piattelli A, Perrotti V, Iezzi G
BUCCO-LINGUAL CRESTAL BONE CHANGES AROUND IMPLANTS IMMEDIATELY PLACED IN FRESH SOCKETS IN ASSOCIATION OR NOT WITH PORCINE BONE: A NON-BLINDED RANDOMIZED CONTROLLED TRIAL IN HUMANS
 J PERIODONTOL, 2012 OCT 29, EPUB AHEAD OF PRINT
- Brunelli G, Sollazzo V, Carinci F, Palmieri A, Girardi R, Monguzzi R
OSTEOBIOL® INFLUENCES OSTEOGENIC DIFFERENTIATION OF ADIPOSE DERIVED STEM CELLS
 EUR J INFLAMMAT, 2011, VOL. 9, NO. 3 (S), 103-107
- Felice P, Piana L, Checchi L, Corvino V, Nannmark U, Piattelli M
VERTICAL RIDGE AUGMENTATION OF ATROPHIC POSTERIOR MANDIBLE WITH AN INLAY TECHNIQUE AND CANCELLOUS EQUINE BONE BLOCK: A CASE REPORT
 INT J PERIODONTICS RESTORATIVE DENT, 2013 MAR;33(2):159-66
- Wachtel H, Fickl S, Hinze M, Bolz W, Thalmair T
THE BONE LAMINA TECHNIQUE: A NOVEL APPROACH FOR LATERAL RIDGE AUGMENTATION - A CASE SERIES
 INT J PERIODONTICS RESTORATIVE DENT, 2013 JUL-AUG;33(4):491-7
- Felice P, Pistilli R, Piattelli M, Soardi E, Corvino V, Esposito M
POSTERIOR ATROPHIC JAWS REHABILITATED WITH PROSTHESES SUPPORTED BY 5 X 5 MM IMPLANTS WITH A NOVEL NANOSTRUCTURED CALCIUM-INCORPORATED TITANIUM SURFACE OR BY LONGER IMPLANTS IN AUGMENTED BONE. PRELIMINARY RESULTS FROM A RANDOMISED CONTROLLED TRIAL
 EUR J ORAL IMPLANTOL, SUMMER;5(2):149-61
- Calvo Guirado JL, Gomez Moreno G, Guardia J, Ortiz Ruiz A, Piattelli A, Barone A, Martinez González Jm, Mesenguer Olmo L, López Mari L, Dorado Cb
BIOLOGICAL RESPONSE TO PORCINE XENOGRAFT IMPLANTS: AN EXPERIMENTAL STUDY IN RABBITS
 IMPLANT DENT, 2012 APR;21(2):112-7
- Figueiredo A, Coimbra P, Cabrita A, Guerra F, Figueiredo M
COMPARISON OF A XENOGENIC AND AN ALLOPLASTIC MATERIAL USED IN DENTAL IMPLANTS IN TERMS OF PHYSICO-CHEMICAL CHARACTERISTICS AND IN VIVO INFLAMMATORY RESPONSE
 MATER SCI ENG C, MATER BIOL APP, 2013 AUG 1;33(6):3506-13
- Pistilli R, Felice P, Piattelli M, Gessaroli M, Soardi E, Barausse C, Buti J, Corvino V, Esposito M
POSTERIOR ATROPHIC JAWS REHABILITATED WITH PROSTHESES SUPPORTED BY 5 X 5 MM IMPLANTS WITH A NOVEL NANOSTRUCTURED CALCIUM-INCORPORATED TITANIUM SURFACE OR BY LONGER IMPLANTS IN AUGMENTED BONE. ONE-YEAR RESULTS FROM A RANDOMISED CONTROLLED TRIAL
 EUR J ORAL IMPLANTOL, 2013 WINTER;6(4):343-357
- Cassetta M, Perrotti V, Calasso S, Piattelli A, Sinjari B, Iezzi G
BONE FORMATION IN SINUS AUGMENTATION PROCEDURES USING AUTOLOGOUS BONE, PORCINE BONE, AND A 50 : 50 MIXTURE: A HUMAN CLINICAL AND HISTOLOGICAL EVALUATION AT 2 MONTHS
 CLIN ORAL IMPLANTS RES, 2014 MAY 26 EPUB AHEAD OF PRINT
- Gheno E, Palermo A, Buffoli B, Rodella LF
THE EFFECTIVENESS OF THE USE OF XENOGENIC BONE BLOCKS MIXED WITH AUTOLOGOUS CONCENTRATED GROWTH FACTORS (CGF) IN BONE REGENERATION TECHNIQUES
 J OSSEOINTEGRATION 2014;6(2):37-42
- Fischer Kr, Stavropoulos A, Calvo Guirado JL, Schneider D, Fickl S
INFLUENCE OF LOCAL ADMINISTRATION OF PAMIDRONATE ON EXTRACTION SOCKET HEALING - A HISTOMORPHOMETRIC PROOF-OF-PRINCIPLE PRE-CLINICAL IN VIVO EVALUATION
 CLIN ORAL IMPLANTS RES, 2014 SEP 15 EPUB AHEAD OF PRINT
- Scarano A, Piattelli A, Murmura G, Iezzi G, Assenza B, Mancino C
DELAYED EXPANSION OF THE ATROPHIC MANDIBLE BY ULTRASONIC SURGERY: A CLINICAL AND HISTOLOGIC CASE SERIES
 INT J ORAL MAXILLOFAC IMPLANTS, 2015 JAN-FEB;30(1):144-9
- Thalmair T, Fickl S, Schneider D, Hinze M, Wachtel H
DIMENSIONAL ALTERATIONS OF EXTRACTION SITES AFTER DIFFERENT ALVEOLAR RIDGE PRESERVATION TECHNIQUES - A VOLUMETRIC STUDY
 J CLIN PERIODONTOL, 2013 JUL;40(7):721-7



Histology at 3 months. Human mandible grafted with OsteoBiol® Sp-Block. Biopsy courtesy of Dr P Felice, Bologna, Italy. Histology courtesy of Prof U Nannmark, University of Göteborg, Sweden

Sp-Block | Dual-Block

HIGHLY OSTEOCONDUCTIVE SCAFFOLDS



Tecnoss s.r.l. is an innovative, globally active company that develops, produces and documents premium-quality xenogenic biomaterials by the brands Tecnoss® and OsteoBiol®.

Its 20 years of research led to its patent-protected production process that ensures neutralization of antigenic components in order to achieve biocompatibility, while preserving the natural collagen matrix inside the biomaterial.

Tecnoss® products comply with highest quality standards such as ISO 10993, ISO 13485 (notified body Kiwa Cermet) and 93/42/EC (notified body CE 0373).

osteobiol.com


Authorized Distributor

Tecnoss® Dental

Via Torino, 23
10044 Pianezza (TO) | Italy
Tel +39 011 9682823
Fax +39 011 9787577
info@tecnoss-dental.com

www.osteobiol.com

International Sales & Marketing



MKT-BLS0415EN

