

Suturer för alla typer av ingrepp

Suturer är en grundläggande del av all kirurgi - från enklare ingrepp till avancerade operationer. Med rätt suturmaterial säkerställs stabil vävnadsanpassning, optimal läkning och hög patientkomfort. Ett brett sortiment ger flexibilitet och precision i varje steg av behandlingen.



dab

Needle Curvings



1/2 Circle




3/8 circle



straight

INFORMATION FÖRPACKNING




ideal
SOLUTION

- | | |
|---|----------------------------|
| 1. Brand Name 1. Artikel | 10. EP metric/USP size |
| 2. Product code | 9. Fadenlänge |
| 2. Produktcode (Nadelcode) | 10. EP metric/USP size |
| 3. Needle curving | 11. Produktbezeichnung |
| 3. Nadelkurve | 12. Sterilization method |
| 4. Needle length | 12. Sterilisierungsmethode |
| 4. Nadellänge | 13. Manufacture date |
| 5. LOT number | 13. Herstellungsdatum |
| 5. Nadellänge | 14. Number of sutures |
| 6. Notified body/ authorization number | 14. Anzahl der Fäden |
| 6. bekannte Stelle/ Genehmigungsnummer | 15. Expiration date |
| 7. Do not reuse | 15. Verfallsdatum |
| 7. Nicht wieder verwendbar (Einmalgebrauch) | 16. Needle point type |
| 8. Read instructions for use | 16. Spitzentyp der Nadel |
| 8. Siehe Gebrauchsanweisungen | 17. Needle drawing |
| 9. Thread length | 18.. Zugfestigkeit |
| 9. Fadenlänge | |

SUTURES


ABSORBABLE

Observera att i alla suturer finns fler 0-typersamt längder och storlekar.




P.G.A.

Braided, violet, synthetic absorbable polyglycolic acid surgical suture, PGA-PCl-copolymer coated.



Neosorb (PGLA)

Braided, violet, synthetic absorbable polyglactin 910 surgical suture, coated.



Neosorb Rapid (RPGLA)

Braided, undyed synthetic rapid absorbable polyglactin 910 surgical suture, coated.

SUTURES


NON-ABSORBABLE

Observera att i alla suturer finns fler 0-typersamt längder och storlekar.



Propylen

Blue monofilament non absorbable polypropylen surgical suture.




Polyamid

Blue monofilament non absorbable polyamid surgical suture (nylon).




Silk

Black, braided, non absorbable, wax coated, silk surgical suture.




Monosorb

Monofilament, violet, synthetic absorbable, poly (p-dioxanone) surgical suture-PDO.



Monofast

Monofilament, violet, synthetic absorbable, polyglycolic acid-co-polycaprolactone surgical suture-PGCL.



Polyester

Green, braided, non absorbable, silicone coated, polyester surgical suture.




Profimed

Undyed monofilament non absorbable, PTFE suture.




Steel

Stainless steel surgical suture 316LVM.



Supramid


Black pseudo-monofilament or monofilament non absorbable polyamide surgical suture.



P.G.A

Absorbable polyglycolic acid, braided violet, coated


Art nr	USP	EP	längd	nål	beskrivning
9609	7/0	0.5	45cm	10mm	3/8 circle, reverse cutting
9610	6/0	0.7	45cm	12mm	3/8 circle, reverse cutting
9611	5/0	1.0	45cm	12mm	3/8 circle, reverse cutting
9620	6/0	0.7	45cm	16mm	3/8 circle, reverse cutting
9621	5/0	1.0	45cm	16mm	3/8 circle, reverse cutting
9622	4/0	1.5	45cm	16mm	3/8 circle, reverse cutting
9623	3/0	2.0	45cm	16mm	3/8 circle, reverse cutting
9631	5/0	1.0	45cm	19mm	3/8 circle, reverse cutting
9632	4/0	1.5	45cm	19mm	3/8 circle, reverse cutting
9633	3/0	2.0	75cm	19mm	3/8 circle, reverse cutting



NEOSORB

Absorbable polyglactin 910, braided, violet, coated


Art nr	USP	EP	längd	nål	beskrivning
N610	6/0	0.7	45cm	12mm	3/8 circle, reverse cutting
N611	5/0	1.0	45cm	12mm	3/8 circle, reverse cutting
N620	6/0	0.7	45cm	16mm	3/8 circle, reverse cutting
N621	5/0	1.0	45cm	16mm	3/8 circle, reverse cutting
N622	4/0	1.5	45cm	16mm	3/8 circle, reverse cutting
N631	5/0	1.0	45cm	19mm	3/8 circle, reverse cutting
N632	4/0	1.5	45cm	19mm	3/8 circle, reverse cutting



NEOSORB

Absorbable polyglactin 910, braided, violet, coated


N633	3/0	2.0	75cm	19mm	1/2 circle, round bodied
N431D	5/0	1.0	45cm	19mm	1/2 circle, round bodied
N922	4/0	1.5	75cm	20mm	1/2 circle, round bodied
N923	3/0	2.0	75cm	20mm	1/2 circle, round bodied
N642	4/0	1.5	75cm	24mm	1/2 circle, round bodied
N643	3/0	2.0	75cm	24mm	1/2 circle, round bodied
N644	2/0	3.0	75cm	24mm	1/2 circle, round bodied
N142	4/0	1.5	75cm	26mm	1/2 circle, round bodied
N143	3/0	2.0	75cm	26mm	1/2 circle, round bodied
N144	2/0	3.0	75cm	26mm	1/2 circle, round bodied
N145	0	3.5	75cm	26mm	1/2 circle, round bodied



NEOSORB RAPID RPGLA

polyglactin 910, braided undyed, coated


Art nr	USP	EP	längd	nål	beskrivning
NR610	6/0	0.7	45cm	12mm	3/8 circle, reverse cutting
NR611	5/0	1.0	45cm	12mm	3/8 circle, reverse cutting
NR620	6/0	0.7	45cm	16mm	3/8 circle, reverse cutting
NR621	5/0	1.0	45cm	16mm	3/8 circle, reverse cutting
NR622	4/0	1.5	45cm	16mm	3/8 circle, reverse cutting
NR631	5/0	1.0	45cm	19mm	3/8 circle, reverse cutting
NR632	4/0	1.5	45cm	19mm	3/8 circle, reverse cutting



MONOFAST

Absorbable PGCL monofilament, violet


Art nr	USP	EP	längd	nål	beskrivning
18610	6/0	0.7	45cm	12mm	3/8 circle, reverse cutting
18611	5/0	1.0	45cm	12mm	3/8 circle, reverse cutting
18612	4/0	1.5	45cm	12mm	3/8 circle, reverse cutting
18621	5/0	1.0	45cm	16mm	3/8 circle, reverse cutting
18622	4/0	1.5	45cm	16mm	3/8 circle, reverse cutting
18631	5/0	1.0	45cm	19mm	3/8 circle, reverse cutting
18632	4/0	1.5	45cm	19mm	3/8 circle, reverse cutting



Profimed

Non absorbable PTFE,
monofilament, undyed

Art nr	USP	EP	längd	nål	beskrivning
PR609S	7/0	0.5	45cm	10mm	3/8 circle, reverse EXTRA cutting
PR610	6/0	0.7	45cm	12mm	3/8 circle, reverse cutting
PR611	5/0	1.0	45cm	12mm	3/8 circle, reverse cutting
PR620	6/0	0.7	45cm	16mm	3/8 circle, reverse cutting
PR621	5/0	1.0	45cm	16mm	3/8 circle, reverse cutting
PR622	4/0	1.5	45cm	16mm	3/8 circle, reverse cutting
PR631	5/0	1.0	45cm	19mm	3/8 circle, reverse cutting
PR632	4/0	1.5	45cm	19mm	3/8 circle, reverse cutting
PR631T	5/0	1.0	45cm	18mm	3/8 circle, round- diamond point
PR632T	4/0	1.5	45cm	18mm	3/8 circle, round- diamond point



Polyamid

Non absorbable
polyamid (nylon)
monofilament, blue

Art nr	USP	EP	längd	nål	beskrivning
4609	7/0	0.5	45cm	10mm	3/8 circle, reverse cutting
4600	6/0	0.7	45cm	10mm	3/8 circle, reverse cutting
4600S	6/0	0.7	45cm	10mm	3/8 circle, reverse EXTRA cutting
4610	6/0	0.7	45cm	12mm	3/8 circle, reverse cutting
4611	5/0	1.0	45cm	12mm	3/8 circle, reverse cutting
4620	6/0	0.7	45cm	16mm	3/8 circle, reverse cutting
4620S	6/0	0.7	45cm	16mm	3/8 circle, reverse EXTRA cutting
4621	5/0	1.0	45cm	16mm	3/8 circle, reverse cutting
4622	4/0	1.5	45cm	16mm	3/8 circle, reverse cutting
4631	5/0	1.0	45cm	19mm	3/8 circle, reverse cutting
4632	4/0	1.5	45cm	19mm	3/8 circle, reverse cutting



Propylen

Non absorbable
polypropylen
monofilament, blue

Art nr	USP	EP	längd	nål	beskrivning
6610	6/0	0.7	45cm	12mm	3/8 circle, reverse cutting
6611	5/0	1.0	45cm	12mm	3/8 circle, reverse cutting
6612	4/0	1.5	45cm	12mm	3/8 circle, reverse cutting
6620	6/0	0.7	45cm	16mm	3/8 circle, reverse cutting
6621	5/0	1.0	45cm	16mm	3/8 circle, reverse cutting
6622	4/0	1.5	45cm	16mm	3/8 circle, reverse cutting
6631	5/0	1.0	45cm	19mm	3/8 circle, reverse cutting
6632	4/0	1.5	45cm	19mm	3/8 circle, reverse cutting

The Dab logo is a green square with the word "dab" in white, lowercase, rounded letters.The background of the entire page is a photograph of a dental procedure. A dentist, wearing a blue surgical cap, glasses, and a blue face mask, is focused on a patient's mouth. The patient is lying back, and the dentist is using a dental instrument. The scene is brightly lit, typical of a dental clinic.

dab

För mer
information
kontakta oss eller
beställ vår
produktkatalog
för Kirurgi
på Dab!

KATALOG KIRURGI
DAB DENTAL

Det här är är
stegen för att
lägga en order
på DAB Dental
Kirurgi:

Här är vi som
jobbar på DAB
Dental Kirurgi



Kundtjänst

Vi hjälper till med eventuella frågor och beställningar. Ring eller maila gärna oss!

Telefon: 08-506 505 05

Mail: kontakt@dabdental.se

WebShop

Logga in eller skapa konto

Klicka på Logga in högst upp på sidan. Om du är ny kund, välj Skapa konto och fyll i dina uppgifter.

Navigera till webbutiken

Använd menyn för att välja kategori: Kirurgi

Välj produkter

Klicka på en produkt för att se detaljer. Välj antal och klicka på Lägg i varukorg.

Gå till varukorgen

Klicka på varukorgsikon. Kontrollera att allt stämmer: antal, pris, produktnamn.

Till kassan

Klicka på Till kassan. Fyll i leveransinformation och välj fraktsätt.

Välj betalningsmetod

Välj mellan tillgängliga alternativ (t.ex. faktura eller kort).

Bekräfta ordern

Granska all information och klicka på Bekräfta köp.

Orderbekräftelse

Du får en bekräftelse via e-post med ordernummer och detaljer.

Leveranstid

Du beräknas ha din order inom 24 h.



Ninni Teljstedt

Försäljningschef
Kirurgi
08-506 505 30
0708-59 34 40
ninni.teljstedt
@dabdental.se



Anna Bäckmark Olsson

Kundtjänst
Förbrukning
Tandsköterska
08-506 505 14
anna.backmark.olsson
@dabdental.se



Jeanette Nilsson

Distriktsansvarig
Säljare Kirurgi
södra & västra Sverige
0708-59 34 54
jeanette.nilsson
@dabdental.se



Tommy Holmén

Senior Advisor
tommy.holmen
@dabdental.se

Följ oss!

Följ vårt nyhetsflöde i sociala medier och håll dig uppdaterad om allt spännande som händer hos oss på DAB.

På Facebook, Instagram och LinkedIn delar vi bland annat produktnyheter, lediga tjänster, exklusiva erbjudanden och förhandsinformation om våra resor, kurser och event. Du får även en inblick i våra tidigare evenemang genom bilder och uppdateringar. Välkommen att hänga med oss!



Användarnamn:
dab_dental



Användarnamn:
dab_dental



linkedin.com/
company/
dab-dental-ab

Kundtjänst förbrukning

08-506 505 05

Kundtjänst utrustning

08-506 505 37

Växel

08-506 505 00

Mail

kontakt@dabdental.se

Hemsida

dabdental.se

e-handel

dabdental.se

Kvalitetsstyrning

Vår samling på kvalitetsstyrning innebär att du får den goda service vi föresätter oss att leverera till dig. Du kan också vara säker på att vi har tänkt igenom hur vi förhåller oss till den miljöpåverkan som vår verksamhet innebär.

DAB Dental AB
Box 423
194 04 Upplands Väsby
Besöksadress: Finvids Väg 8